

product information


| | | | | |
|--|---|---|--|--|
|  Electric or Hydraulic height adjustment 45-96cm |  Maximum user weight (35st.) 225kg |  Gas strut assisted multi-position head and foot sections |  Back rest and leg rest adjustment +70° / -35° |  Easy to use handset (Electric models only) |
|  Retractable positioning castors |  Adjustable height levelling foot |  Fully upholstered head & foot sections |  Double layer density foam for support & comfort |  Luxury 65mm deep contoured upholstery |
|  Anti-bacterial/anti-microbial vinyl |  Paper roll holder (supplied as standard) |  Easy clean design and materials |  3 year warranty on all moving components |  Lifetime warranty on frame |

accessories

* these accessories can only be fitted when ordered with the couch

| | | | | |
|---|---|--|---|---|
|  Sun-FVHC/LEG1 Leg Supports |  Sun-FVHC/PHB1 Phlebotomy Arm |  Sun-FVHC/FSW1* Foot Switch |  Sun-FVHC/SUP1 Gynae Stirrups |  Sun-CHP1 Head Pillow |
|  Sun-FVHC/CBH1* Breathing Hole |  Sun-FVHC/CSR Folding Side Rails (pair) |  Sun-FVHC/IV2 IV Pole Attachment |  Sun-FVHC/RAIL Medi-Rail Attachment |  Sun-FVHC/BB* Battery Backup |
|  Sun-FVHC/HEO* Electrically Operated Head Section Upgrade |  Sun-FVHC/FEO* Electrically Operated Foot Section Upgrade |  Sun-INST1* Installation Service & Packaging Recycling |  Sun-DISPOSAL1* Removal & Disposal of Old Couch | |


dimensions



185cm(L) 65cm(W)




anti-bacterial vinyl upholstery colours - please specify when ordering please see inside back cover for full specification of our vinyl range



The new generation of Fusion Classic has been designed to assist in making treatment procedures even easier and more efficient than ever. Fusion Classic is simply, strong, stable, aesthetically pleasing and easy to clean.



Fully upholstered head & foot sections with gas assisted adjustment



Electric or Hydraulic height adjustment



Increased maximum user weight to 225kg



Code: Sun-FVHC3E
Fusion Classic 3 Section Electric Couch (Shown with optional accessories)


Optional
Electric Head & Foot Sections (electric models only)



Normally delivered within 7-10 working days (UK Only)



Hydraulic height adjustment





Code: Sun-FVHC3H
Fusion Classic 3 Section Hydraulic Couch (Shown with optional pillow)

product information


| | | | | |
|--|---|---|--|--|
|  Electric or Hydraulic height adjustment 45-96cm |  Maximum user weight (35st.) 225kg |  Gas strut assisted multi-position head section |  Back rest adjustment +70° / -35° |  Easy to use handset (Electric models only) |
|  Retractable positioning castors |  Adjustable height levelling foot |  Fully upholstered head section |  Double layer density foam for support & comfort |  Luxury 65mm deep contoured upholstery |
|  Anti-bacterial/anti-microbial vinyl |  Paper roll holder (supplied as standard) |  Easy clean design and materials |  3 year warranty on all moving components |  Lifetime warranty on frame |

accessories

* these accessories can only be fitted when ordered with the couch

| | | | | |
|---|--|---|---|---|
|  Sun-FVHC/LEG1 Leg Supports |  Sun-FVHC/PHB1 Phlebotomy Arm |  Sun-FVHC/FSW1* Foot Switch |  Sun-FVHC/SUP1 Gynae Stirrups |  Sun-CHP1 Head Pillow |
|  Sun-FVHC/CBH1* Breathing Hole |  Sun-FVHC/CSR Folding Side Rails (pair) |  Sun-FVHC/IV2 IV Pole Attachment |  Sun-FVHC/RAIL Medi-Rail Attachment |  Sun-FVHC/BB* Battery Backup |
|  Sun-FVHC/HEO* Electrically Operated Head Section Upgrade |  Sun-INST1* Installation Service & Packaging Recycling |  Sun-DISPOSAL1* Removal & Disposal of Old Couch | | |

dimensions



185cm(L) 65cm(W)




anti-bacterial vinyl upholstery colours - please specify when ordering please see inside back cover for full specification of our vinyl range



The new generation of Fusion Classic has been designed to assist in making treatment procedures even easier and more efficient than ever. Fusion Classic is simply, strong, stable, aesthetically pleasing and easy to clean.



Fully upholstered head section with gas assisted adjustment




Electric or Hydraulic height adjustment



Increased maximum user weight to 225kg



Code: Sun-FVHC2E
Fusion Classic 2 Section Electric Couch


Optional
Electric Head Section
(electric models only)



Normally delivered within 7-10 working days (UK Only)



Code: Sun-FVHC2H
Fusion Classic 2 Section Hydraulic Couch



Hydraulic height adjustment