



Medical Equipment Product Range





UK's leading **manufacturer**
of healthcare furniture



Based in Braford, West Yorkshire we are proud of our UK manufacturing facility which undertakes its design, manufacturing, assembly and distribution all from one central site. Having the ability to design new innovative products, manufacture and supply them from start to finish allows us to provide bespoke solutions and outstanding levels of service.



our promise

To create smart, contemporary, functional furniture, we start with intuitive design and solid engineering. Then, we use the best materials and subject our products to rigorous testing. That's what makes a Sunflower product and why our warranty and service is a **real promise**.



We **manufacture**
on site and in the UK



Our CNC, machining, drilling, spraying and pressing facilities enable us to manufacture accurately and quickly in response to customer requirements, whether it be for a single item or a large quantity in a specific colour or finish, while still being flexible on delivery times.



In-house design with
3D CAD technology



Our specialised product designers have many years of experience in the healthcare market. Designing from 'concept to reality' using the latest solid modelling CAD software our designers can create anything from a new single product, something bespoke or a whole new room or environment to your exact requirements.

our products & service

We design medical furniture to help healthcare become more efficient, combined with high service levels, extensive warranties and consideration for the environment.

We believe passionately in innovative design and continuous development, utilise modern materials and technology to produce stylish and hygienic new products and solutions.

Ease of cleaning and hygienic design is carefully incorporated into everything we do.

We are committed to reducing our impact on the environment and utilise sustainable materials and sources wherever possible, achieving ISO and FSC standards. We are committed to recognised industry standards, giving you confidence in our company, our products and our people.



Certified Number 5288
ISO 9001



Certified Number 5288
ISO 14001



FSC The mark of
responsible forestry



Membership No:
00334186



FIRA Standard 6250
2005 Specification:
Adhesion Performance
General/Severe



BS EN 16139: 2013
Furniture Strength,
Durability & Safety



Anti-bacterial Vinyl Strength Test
BS EN ISO 13934-1: 2013
BS EN ISO 13973-3: 2000
BS EN ISO 12947-2: 1998



Fire (Flammability) Testing
Flammability BS 7176 Medium Hazard
BS 5852 : 2006 Clause 11
(composites)



Suppliers to the NHS

NHS Supply Chain Contract:
Medical Furniture & Sack Holders (FAG15887–OJEU ref: 2014/S 198-349358)

NHS Supply Chain Contract:
Ward & Residential Furniture (FAG000016544–OJEU ref: 2017/S 056-104552)

NHS Shared Business Services Contract:
Medical Furniture (SBS/16/SB/TBP/9050/12 – OJEU 2016/S 174-312172)

For more details visit my.supplychain.nhs.uk/catalogue



Warranty

We provide a minimum 1 year warranty on all products, with the majority having between 2 to 5 years and some components for life, giving you complete peace of mind.



Express Delivery

Need an item quickly? - We offer a range of products through our express system with delivery normally within 3-5 working days, just look for the Express symbol.



Delivery & Installation

We can deliver and install your new furniture anywhere in the UK using our own vehicles and trained installers, please enquire for further details.



New Build, Equipping and Refurbishment

We offer a bespoke design, manufacture, fitting and project management service for new or existing healthcare premises. Visit our website for our full service solutions and case studies.



Export

We can export our products worldwide and have established distributors around the globe. Contact our export department at export@sunflowermedical.co.uk for further details, or if you would like to become part of our global distributor network.

Contact Us & How To Order

T: +44 (0)1274 684004

F: +44 (0)1274 684006

E: sales@sunflowermedical.co.uk

Contact us with your requirements or for details of our approved national distributors.

Alternatively buy online at: www.sunflowermedical.co.uk

new products



We believe designing new products is a process of evolution, and is one element that makes Sunflower different. We start with detailed customer research, so we understand your needs and requirements, and then we use the latest CAD technology to model initial concepts, build prototypes, and then test them in the environment prior to production. That's why and how we design **new products.**



contents

Couches, Plinths & Accessories	Page 1 - 43
Phlebotomy Chairs, Operator Chairs & Stools	Page 44 - 54
Vista Trolleys & Storage	Page 55 - 91
Dressing & Surgical Trolleys, Theatre Equipment	Page 92 - 106
Medicine & Controlled Drug Cabinets	Page 107 - 125
Bins & Sack Holders	Page 126 - 130
Mobile Screens	Page 131 - 135
Mirrors & Mobile Drip Stands	Page 136 - 140
Furniture Packages	Page 141 - 143
Cabinet Ranges	Page 144 - 148
Integrated Plumbing Systems, Sinks & Taps	Page 149 - 161
Visitor, Patient & Waiting Room Seating	Page 162 - 190
Desks, Chairs, Filing, & Administration Storage	Page 191 - 207
Drug, Medicine, Records Trolleys, & Pharmacy Fridges	Page 208 - 224
Overbed Tables, Bedside Furniture & Patient Seating	Page 225 - 251
Recommended Cleaning Guidelines	Page 252 - 253
Upholstery and KYDEX® Colours	Page 254 - 256

clean
simple
examination

Couches & Accessories

We design couches for most medical applications, from treatment, gynecological and examination to resting and portable use, all with optional accessories to make procedures easier.

	Fusion Classic Plus+ Variable Height Couches	PAGE 2-3		Treatment Couch Systems	PAGE 28-35
	Fusion Classic Variable Height Couches	4-7		Deluxe Examination Couches	36-37
	Fusion Eco Variable Height Couches	8-11		Examination Couches	38-39
	Fusion Gynae Variable Height Couches	12-15		Rest Couches	40-41
	Fusion Drop End Multi-Discipline Couches	16-17		Portable Couches	42-43
	Fusion Podiatry Chairs	18-19		Couch Steps	42-43
	Fusion Treatment Chairs	20-21			
	Treatment Chairs	22-23			
	Fusion Phlebotomy Chairs	24-25			
	Phlebotomy Chairs	26-27			

product information

Fully powered height adjustment 45-96cm

Fully powered multi-position head section +70° / -35°

Fully powered multi-position foot section +70° / -35°

Maximum user weight (35st.)

Easy to use handset

Retractable positioning castors

Adjustable height levelling foot

Fully upholstered head & foot sections

Double layer density foam for support & comfort

Luxury 65mm deep contoured upholstery

Anti-bacterial/anti-microbial vinyl

Paper roll holder (supplied as standard)

Easy clean design and materials

FIRA Independently tested by FIRA

West Yorkshire Materials Testing Service

Anti-bacterial Vinyl

Anti-bacterial Vinyl Strength Test BS EN ISO 13934-1: 2013 BS EN ISO 13973-3: 2000 BS EN ISO 12947-2: 1998

Fire (Flammability) Testing BS 5852 : 2006 Clause 11 (composites)

accessories * these accessories can only be fitted when ordered with the couch

Sun-FVHC/LEG1 Leg Supports

Sun-FVHC/PHB1 Phlebotomy Arm (each)

Sun-FVHC/FSW1* Foot Switch

Sun-FVHC/FSW2* Foot Switch 2 Section

Sun-FVHC/FSW3* Foot Switch 3 Section

Sun-FVHC/SUP1 Gynae Stirrups

Sun-CHP1 Head Pillow

Sun-FVHC/CBH1* Breathing Hole

Sun-FVHC/CSR* Folding Side Rails (pair)

Sun-FVHC/IV2* IV Pole Attachment

Sun-FVHC/RAIL* Medi-Rail Attachment


Sun-FVHC/BB* Battery Backup

Sun-INST1* Installation Service & Packaging Recycling


Sun-DISPOSAL1** Removal & Disposal of Old Couch

** Removal from ground floor only - weight & removal restrictions apply, please enquire

dimensions



anti-bacterial vinyl upholstery colours - please specify when ordering
please see inside back cover for full specification of our vinyl range



Fusion Classic Plus+
2 & 3 Section Fully Powered Electric Couches

The new Fusion Classic Plus+ extends our already popular Fusion range, providing enhanced user functionality and easier patient positioning through its fully powered height, head and foot sections. Fusion Classic Plus+ is simply, strong, stable, aesthetically pleasing and easy to clean.



Fully powered height adjustment



Fully powered head section



Fully powered foot section



Code: **Sun-FVHC+3E**
Fusion Classic Plus+ 3 Section Electric Couch with powered height, head & foot section (Shown with optional accessories)



Normally delivered within 7-10 working days (*UK Only)



Code: **Sun-FVHC+2E**
Fusion Classic Plus+ 2 Section Electric Couch with powered height & head section



3 year warranty on all moving components



Lifetime warranty on frame



product information

Electric or Hydraulic height adjustment 45-96cm

Maximum user weight (35st.)

Gas strut assisted multi-position head and foot sections

Back rest and leg rest adjustment +70° / -35°

Easy to use handset (Electric models only)

Retractable positioning castors

Adjustable height levelling foot

Fully upholstered head & foot sections

Double layer density foam for support & comfort

Luxury 65mm deep contoured upholstery

Anti-bacterial/anti-microbial vinyl

Paper roll holder (supplied as standard)

Easy clean design and materials

Anti-bacterial Vinyl Strength Test
BS EN ISO 13934-1: 2013
BS EN ISO 13973-3: 2000
BS EN ISO 12947-2: 1998

Fire (Flammability) Testing
Flammability BS 7176 Medium Hazard
BS 5852 : 2006 Clause 11 (composites)

Anti-bacterial Vinyl

accessories * these accessories can only be fitted when ordered with the couch

Sun-FVHC/LEG1
Leg Supports

Sun-FVHC/PHB1
Phlebotomy Arm (each)

Sun-FVHC/FSW1*
Foot Switch

Sun-FVHC/SUP1
Gynae Stirrups

Sun-CHP1
Head Pillow

Sun-FVHC/CBH1*
Breathing Hole

Sun-FVHC/CSR*
Folding Side Rails (pair)

Sun-FVHC/IV2*
IV Pole Attachment

Sun-FVHC/RAIL*
Medi-Rail Attachment

Sun-FVHC/BB*
Battery Backup

Sun-FVHC/HEO*
Electrically Operated Head Section Upgrade

Sun-FVHC/FEO*
Electrically Operated Foot Section Upgrade

Sun-INST1*
Installation Service & Packaging Recycling

Sun-DISPOSAL1**
Removal & Disposal of Old Couch

** Removal from ground floor only - weight & removal restrictions apply, please enquire

dimensions

185cm(L) 65cm(W)

HEAVY ITEM

anti-bacterial vinyl upholstery colours - please specify when ordering
please see inside back cover for full specification of our vinyl range

black grey beige primrose walnut red wine salmon lilac navy mid blue sky blue cool blue green mint white

Fusion Classic - 3 Section Variable Height Couches

The new generation of Fusion Classic has been designed to assist in making treatment procedures even easier and more efficient than ever. Fusion Classic is simply, strong, stable, aesthetically pleasing and easy to clean.

Fully upholstered head & foot sections with gas assisted adjustment

Electric or Hydraulic height adjustment

Increased maximum user weight to 225kg

Optional
Electric Head &
Foot Sections
(electric models only)

3 year warranty on all moving components

Lifetime warranty on frame

FSC
www.fsc.org
FSC® C122861
The mark of responsible forestry



Code: **Sun-FVHC3E**
Fusion Classic 3 Section Electric Couch (Shown with optional accessories)



Code: **Sun-FVHC3H**
Fusion Classic 3 Section Hydraulic Couch (Shown with optional pillow)

product information



Electric or Hydraulic height adjustment 45-96cm



225kg
Maximum user weight (35st.)



Gas strut assisted multi-position head and foot sections



Back rest and leg rest adjustment +70° / -35°



Easy to use handset (Electric models only)



Retractable positioning castors



Adjustable height levelling foot



Fully upholstered head & foot sections



Double layer density foam for support & comfort



Luxury 65mm deep contoured upholstery



Anti-bacterial/anti-microbial vinyl



Paper roll holder (supplied as standard)



Easy clean design and materials




Anti-bacterial Vinyl Strength Test
BS EN ISO 13934-1: 2013
BS EN ISO 13973-3: 2000
BS EN ISO 12947-2: 1998




Fire (Flammability) Testing
Flammability BS 7176 Medium Hazard
BS 5852 : 2006 Clause 11 (composites)


accessories * these accessories can only be fitted when ordered with the couch




Sun-FVHC/LEG1
Leg Supports




Sun-FVHC/PHB1
Phlebotomy Arm (each)




Sun-FVHC/FSW1*
Foot Switch




Sun-FVHC/SUP1
Gynae Stirrups




Sun-CHP1
Head Pillow




Sun-FVHC/CBH1*
Breathing Hole




Sun-FVHC/CSR*
Folding Side Rails (pair)




Sun-FVHC/IV2*
IV Pole Attachment




Sun-FVHC/RAIL*
Medi-Rail Attachment




Sun-FVHC/BB*
Battery Backup



Sun-FVHC/HEO*
Electrically Operated Head Section Upgrade



installation
Sun-INST1*
Installation Service & Packaging Recycling




Sun-DISPOSAL1*/**
Removal & Disposal of Old Couch

** Removal from ground floor only - weight & removal restrictions apply, please enquire

dimensions




185cm(L) 65cm(W)



anti-bacterial vinyl upholstery colours - please specify when ordering

please see inside back cover for full specification of our vinyl range




black grey beige primrose walnut red wine salmon lilac navy mid blue sky blue cool blue green mint white

The new generation of Fusion Classic has been designed to assist in making treatment procedures even easier and more efficient than ever. Fusion Classic is simply, strong, stable, aesthetically pleasing and easy to clean.



Fully upholstered head section with gas assisted adjustment



Electric or Hydraulic height adjustment



Increased maximum user weight to 225kg



Optional Electric Head Section (electric models only)



3 year warranty on all moving components



Lifetime warranty on frame



Code: Sun-FVHC2E
Fusion Classic 2 Section Electric Couch



Code: Sun-FVHC2H
Fusion Classic 2 Section Hydraulic Couch



Hydraulic height adjustment

product information



Electric or Hydraulic height adjustment 45-96cm



Maximum user weight (40st.) 250kg



Gas strut assisted multi-position head and foot sections



Back rest and leg rest adjustment +70° / -35°



Easy to use handset (Electric models only)



Retractable positioning castors



Adjustable height levelling foot



Fully upholstered head & foot sections



Double layer density foam for support & comfort



Luxury 65mm deep contoured upholstery



Anti-bacterial/anti-microbial vinyl



Paper roll holder (supplied as standard)



Easy clean design and materials



FIRA Independently tested by FIRA



Anti-bacterial Vinyl


Anti-bacterial Vinyl Strength Test
BS EN ISO 13934-1: 2013
BS EN ISO 13973-3: 2000
BS EN ISO 12947-2: 1998

Fire (Flammability) Testing
Flammability BS 7176 Medium Hazard
BS 5852 : 2006 Clause 11 (composites)


accessories * these accessories can only be fitted when ordered with the couch




Sun-FECO/LEG1 Leg Supports




Sun-FECO/PHB1 Phlebotomy Arm (each)




Sun-FECO/FSW1* Foot Switch




Sun-FECO/SUP1 Gynae Stirrups



Sun-CHP1 Head Pillow




Sun-FECO/CBH1* Breathing Hole



Sun-FECO/CSR* Folding Side Rails (pair)



Sun-FECO/IV2* IV Pole Attachment




Sun-FECO/RAIL* Medi-Rail Attachment




Sun-FECO/BB* Battery Backup




Sun-FECO/HEO* Electrically Operated Head Section Upgrade



Sun-FECO/FEO* Electrically Operated Foot Section Upgrade




installation



Sun-DISPOSAL1** Removal & Disposal of Old Couch

** Removal from ground floor only - weight & removal restrictions apply, please enquire

dimensions




185cm(L) 65cm(W)

HEAVY ITEM

anti-bacterial vinyl upholstery colours - please specify when ordering

please see inside back cover for full specification of our vinyl range




black grey beige primrose walnut red wine salmon lilac navy mid blue sky blue cool blue green mint white

Fusion Eco - 3 Section Variable Height Couches

The new Fusion Eco couch range provides most of the functionality you would expect for our standard Fusion range at a more affordable price to suit individual budgets. Fusion Eco is simply, strong, stable, aesthetically pleasing and easy to clean.



Fully upholstered head & foot sections with gas assisted adjustment



Electric or Hydraulic height adjustment



Increased maximum user weight to 250kg



Optional Electric Head & Foot Sections (electric models only)



3 year warranty on all moving components



Lifetime warranty on frame



Code: Sun-FECO3E
Fusion Eco 3 Section Electric Couch (Shown with optional accessories)



Code: Sun-FECO3H
Fusion Eco 3 Section Hydraulic Couch



Hydraulic height adjustment

product information



Electric or Hydraulic height adjustment 45-96cm



250kg
Maximum user weight (40st.)



Gas strut assisted multi-position head and foot sections



Back rest and leg rest adjustment +70° / -35°



Easy to use handset (Electric models only)



Retractable positioning castors



Adjustable height levelling foot



Fully upholstered head & foot sections



Double layer density foam for support & comfort



Luxury 65mm deep contoured upholstery



Anti-bacterial/anti-microbial vinyl



Paper roll holder (supplied as standard)



Easy clean design and materials



FIRA
Independently tested by FIRA




West Yorkshire Materials Testing Service
Anti-bacterial Vinyl


Anti-bacterial Vinyl Strength Test
BS EN ISO 13934-1: 2013
BS EN ISO 13973-3: 2000
BS EN ISO 12947-2: 1998

Fire (Flammability) Testing
Flammability BS 7176 Medium Hazard
BS 5852 : 2006 Clause 11 (composites)


accessories * these accessories can only be fitted when ordered with the couch




Sun-FECO/LEG1
Leg Supports




Sun-FECO/PHB1
Phlebotomy Arm (each)




Sun-FECO/FSW1*
Foot Switch




Sun-FECO/SUP1
Gynae Stirrups




Sun-CHP1
Head Pillow




Sun-FECO/CBH1*
Breathing Hole




Sun-FECO/CSR*
Folding Side Rails (pair)




Sun-FECO/IV2*
IV Pole Attachment




Sun-FECO/RAIL*
Medi-Rail Attachment




Sun-FECO/BB*
Battery Backup



Sun-FECO/HEO*
Electrically Operated Head Section Upgrade



Sun-INST1*
Installation Service & Packaging Recycling



Sun-DISPOSAL1**
Removal & Disposal of Old Couch

** Removal from ground floor only - weight & removal restrictions apply, please enquire

dimensions



185cm(L) 65cm(W)

45-96cm(H)

+70°


-35°

0°

HEAVY ITEM

anti-bacterial vinyl upholstery colours - please specify when ordering

please see inside back cover for full specification of our vinyl range




black grey beige primrose walnut red wine salmon lilac navy mid blue sky blue cool blue green mint white

Fusion Eco - 2 Section Variable Height Couches

The new Fusion Eco couch range provides most of the functionality you would expect for our standard Fusion range at a more affordable price to suit individual budgets. Fusion Eco is simply, strong, stable, aesthetically pleasing and easy to clean.



Fully upholstered head section with gas assisted adjustment



Electric or Hydraulic height adjustment



Increased maximum user weight to 250kg



Optional Electric Head Section (electric models only)



3 year warranty on all moving components



Lifetime warranty on frame



Code: Sun-FECO2E
Fusion Eco 2 Section Electric Couch (Shown with optional accessories)



Code: Sun-FECO2H
Fusion Eco 2 Section Hydraulic Couch



Hydraulic height adjustment

product information



Electric or Hydraulic height adjustment 45-98cm



Maximum user weight (35st.)



Gas assisted back (-20 - +80°)



Easy to use handset (Electric models only)



Rear retractable positioning castors



100mm Front braked castors



Cut away lithotomy section for easier access



100mm Deep stainless steel operators tray



Fully adjustable leg supports



Fold down arms for patient comfort & security



Double layer density foam for support & comfort



Anti-bacterial/anti-microbial vinyl



Easy clean design and materials



Anti-bacterial Vinyl Strength Test
BS EN ISO 13934-1: 2013
BS EN ISO 13973-3: 2000
BS EN ISO 12947-2: 1998




Fire (Flammability) Testing
Flammability BS 7176 Medium Hazard
BS 5852: 2006 Clause 11 (composites)


accessories * these accessories can only be fitted when ordered with the couch



Sun-FGYN/FSW1*
Foot Switch (For Gynae 1)




Sun-FGYN/CHP1
Head Pillow



Sun-FGYN/CBH1*
Breathing Hole




Sun-FGYN/IV2*
IV Pole Attachment




Sun-FGYN/SUP1
Stirrups (pair)




Sun-FGYN/FS
Detachable Foot Section Extension




Sun-FGYN/WB
Wall Bracket for Foot Section



Sun-FGYN/BB*
Battery Backup




Sun-INST1*
Installation Service & Packaging Recycling




Sun-DISPOSAL1*/**
Removal & Disposal of Old Couch

** Removal from ground floor only - weight and removal restrictions apply, please enquire


dimensions




-20 - +80°




45-98cm(H)




Detachable Foot Section Extension
50cm(L) 65cm(W)




anti-bacterial vinyl upholstery colours - please specify when ordering
please see inside back cover for full specification of our vinyl range




black




grey



beige




primrose




walnut




red wine




salmon




lavender




navy




mid blue




sky blue




cool blue



green



mint



white

The new Gynae couch range features an ergonomic design providing comfort for the patient and ease of access and patient positioning for the clinician. Fusion Gynae is simply, strong, stable, aesthetically pleasing and easy to clean.



Code: **Sun-FGYNE1**
Fusion Gynae1 - 2 section couch
Electric powered height, gas assisted back,
no seat tilt with arm supports & leg supports



3 year warranty on all moving components



Lifetime warranty on frame



Normally delivered within 7-10 working days (*UK Only)

Code: **Sun-FGYNH1**
Fusion Gynae1 - 2 section couch
Hydraulic powered height, gas assisted back,
no seat tilt with arm supports & leg supports



Hydraulic height adjustment

product information



Electric height adjustment
45-98cm



225kg
Maximum user weight (35st.)



Electric Powered Tilting Seat 0 to -15°
(Trendelenburg)



Fully powered head section
(Sun-FGYNE3)



Easy to use handset



Rear retractable positioning castors



100mm Front braked castors



Cut away lithotomy section for easier access



100mm Deep stainless steel operators tray



Fully adjustable leg supports



Fold down arms for patient comfort & security



Double layer density foam for support & comfort



Easy clean design and materials



FIRA
Independently tested by FIRA



West Yorkshire Materials Testing Service
Anti-bacterial Vinyl

Anti-bacterial Vinyl Strength Test
BS EN ISO 13934-1: 2013
BS EN ISO 13973-3: 2000
BS EN ISO 12947-2: 1998

Fire (Flammability) Testing
Flammability BS 7176 Medium Hazard
BS 5852 : 2006 Clause 11 (composites)

accessories * these accessories can only be fitted when ordered with the couch



Sun-FGYN/FSW2*
Foot Switch (Gynae 2)



Sun-FGYN/FSW3*
Foot Switch (Gynae 3)



Sun-FGYN/CHP1
Head Pillow



Sun-FGYN/CBH1*
Breathing Hole




Sun-FGYN/IV2*
IV Pole Attachment



Sun-FGYN/SUP1
Stirrups (pair)



Sun-FGYN/FS
Detachable Foot Section Extension



Sun-FGYN/WB
Wall Bracket for Foot Section



Sun-FGYN/BB*
Battery Backup



Sun-INST1*
Installation Service & Packaging Recycling



Sun-DISPOSAL1/**
Removal & Disposal of Old Couch

** Removal from ground floor only - weight & removal restrictions apply, please enquire

dimensions



137cm(L) 65cm(W)
187cm with foot extension

45-98cm(H)
-20 - +80°
Seat Tilt 0 to -15°

Detachable Foot Section Extension
50cm(L) 65cm(W)



HEAVY ITEM

anti-bacterial vinyl upholstery colours - please specify when ordering

please see inside back cover for full specification of our vinyl range



black grey beige primrose walnut red wine salmon lilac navy mid blue sky blue cool blue green mint white

Fusion Gynae 2 & 3 - Electric Couches

The new Gynae couch range features an ergonomic design providing comfort for the patient and ease of access and patient positioning for the clinician. Fusion Gynae is simply, strong, stable, aesthetically pleasing and easy to clean.



Code: **Sun-FGYNE2**
Fusion Gynae2 - 2 section couch
Electric powered height, gas assisted back, powered tilting seat with arm supports & leg supports



3 Year Warranty
on all moving components



Lifetime Warranty
on frame



CE



British Made
Quality Assured



7-10 DAYS
DELIVERY*

Normally delivered within
7-10 working days (*UK Only)




FSC
www.fsc.org
FSC® C122861
The mark of responsible forestry

Code: **Sun-FGYNE3**
Fusion Gynae3 - 2 section couch
Electric powered height, back and tilting seat with arm supports & leg supports




Fusion Drop End Multi-Discipline Couch


product information




Electric or Hydraulic height adjustment




Maximum user weight (40st.)




Optional breathing hole & plug in foot section for physio use




50mm Retractable positioning castors




Adjustable height levelling foot




Gas assisted head & foot sections as standard




Electric head/foot section adjustment (Sun-FDECE2)




Tilting Seat for Trendelenburg (Sun-FDECE2)




Hoist clearance 13cm




Range of additional accessories for flexibility




Fully upholstered head & foot sections




Double layer density foam for support & comfort




Luxury 65mm deep contoured upholstery




Anti-bacterial/anti-microbial vinyl




Paper roll holder (supplied as standard)



Easy clean design and materials



Anti-bacterial Vinyl Strength Test
BS EN ISO 13934-1: 2013
BS EN ISO 13973-3: 2000
BS EN ISO 12947-2: 1998



Fire (Flammability) Testing
Flammability BS 7176 Medium Hazard
BS 5852 : 2006 Clause 11 (composites)

accessories * these accessories can only be fitted when ordered with the couch



Sun-CHP1
Head Pillow



Sun-FDEC/CBH2*
Breathing Hole - in foot section



Sun-FDEC/LEG1
Leg Supports (pair)



Sun-FDEC/PHB1
Phlebotomy Arm (pair)



Sun-FDEC/PAS1
Padded Arm Supports (pair)



Sun-FDEC/SUP1
Stirrups (pair)



Sun-FDEC/CSR*
Folding Side Rails (pair)



Sun-FDEC/IV2*
IV Pole Attachment




Sun-FDEC/FSW1*
Foot Switch



Sun-DISPOSAL1**
Removal & Disposal of Old Couch


** Removal from ground floor only - weight and removal restrictions apply, please enquire

dimensions



anti-bacterial vinyl upholstery colours - please specify when ordering

please see inside back cover for full specification of our vinyl range



Fusion Drop End Multi-Discipline Couch

The new Fusion Drop End Couch range has been designed using our extensive knowledge to provide a trusted, ergonomic design providing comfort for the patient, ease of access via the drop end design and comfortable patient positioning for the clinician. The couch offers amazing flexibility and can be configured into a wide range of specialist treatment uses including physio, gynae, phlebotomy or treatment chair via the addition of relevant accessories.



Code: **Sun-FDECE2**
Fusion Drop End Multi-Discipline Couch
Electric Height Adjustment
Electric back and foot sections



Code: **Sun-FDECE1**
Fusion Drop End Multi-Discipline Couch
Electric Height Adjustment
Gas assisted head and foot sections



3 year warranty on all moving components



Lifetime warranty on frame




Normally delivered within 7-10 working days (*UK Only)

Code: **Sun-FDECH1**
Fusion Drop End Multi-Discipline Couch
Hydraulic Height Adjustment
Gas assisted head and foot sections




Hydraulic height adjustment


product information




Electric or Hydraulic height adjustment




Maximum user weight (40st.)




50mm Retractable positioning castors




Gas assisted head and split foot sections




Electric head section adjustment (Sun-FPODE2)




Adjustable leg length extension upto 15cm




Tilting Seat for Trendelenburg (Sun-FPODE1 & Sun-FPODE2)




Stay level swing back padded arm rests for easy patient access




Adjustable height levelling foot




Hoist clearance 13cm




Fully upholstered head & foot sections




Double layer density foam for support & comfort




Luxury 65mm deep contoured upholstery




Anti-bacterial/anti-microbial vinyl




Paper roll holder (supplied as standard)



Easy clean design and materials



Anti-bacterial Vinyl Strength Test
BS EN ISO 13934-1: 2013
BS EN ISO 13973-3: 2000
BS EN ISO 12947-2: 1998



Fire (Flammability) Testing
Flammability BS 7176 Medium Hazard
BS 5852 : 2006 Clause 11 (composites)

accessories * these accessories can only be fitted when ordered with the couch



Sun-CHP1 Head Pillow



Sun-FPOD/IV2* IV Pole Attachment




Sun-FPOD/FSW1* Foot Switch (FPODE1)



Sun-FPOD/FSW2* Foot Switch (FPODE2)



Sun-FPOD/FSW2* Foot Switch (FPODE2)




installation



Sun-DISPOSAL1*/** Removal & Disposal of Old Couch

** Removal from ground floor only - weight & removal restrictions apply, please enquire

dimensions



anti-bacterial vinyl upholstery colours - please specify when ordering

please see inside back cover for full specification of our vinyl range



The new Fusion Podiatry Chair range has been designed using our extensive knowledge to provide a trusted, ergonomic design providing comfort for the patient, ease of access via the split leg design, Trendelenburg movement as standard on electric models and comfortable patient positioning for the clinician.



Code: **Sun-FPODE1**
Fusion Podiatry Chair
Electric Height Adjustment
Gas assisted head and split adjustable length foot sections
(Powered Trendelenburg -20° Movement)

Code: **Sun-FPODE2** (shown above)
Fusion Podiatry Chair
Electric Height Adjustment
Electric head and split adjustable length foot sections
(Powered Trendelenburg -20° Movement)



3 year warranty on all moving components



Lifetime warranty on frame



Normally delivered within 7-10 working days (*UK Only)

Code: **Sun-FPODH1**
Fusion Podiatry Chair
Hydraulic Height Adjustment
Gas assisted head and split adjustable length foot sections
(No Trendelenburg Movement)



Hydraulic height adjustment

product information



Powered tilting seat
Sun-FTRTE2/3/5/6 only)



-20°
TREDELENBURG
Tilting Seat for Trendelenburg
Sun-FTRTE2/3/5/6 only)



Head adjustment
+80° / -20°



Adjustable foot section
+20° / -90°



2 Adjustable/removable
padded arm supports



250kg
Maximum user
weight (40st.)



Easy to use handset



50mm Retractable
positioning castors



Adjustable height
levelling foot



Dual sided foot pedals
(Hydraulic models only)



Hoist clearance
13cm



Fully upholstered head
& foot sections



Double layer density
foam for support & comfort



Luxury 65mm deep
contoured upholstery



Anti-bacterial/
anti-microbial vinyl



Paper roll holder
(supplied as standard)



Easy clean design
and materials



FIRA
Independently
tested by FIRA




West Yorkshire
Materials
Testing Service
Anti-bacterial Vinyl

Anti-bacterial Vinyl Strength Test
BS EN ISO 13934-1: 2013
BS EN ISO 13973-3: 2000
BS EN ISO 12947-2: 1998

Fire (Flammability) Testing
Flammability BS 7176 Medium Hazard
BS 5852 : 2006 Clause 11
(composites)


accessories * these accessories can only be fitted when ordered with the couch




Sun-FTRT/PHB1
Phlebotomy Arm (pair)




Sun-CHP1
Head Pillow




Sun-FTRT/CBH1*
Breathing Hole - in foot section
(Sun-FTRTE1,2,3 & Sun-FTRTH1 only)




Sun-FTRT/FSW1*
Foot Switch
(Sun-FTRTE1 & 4)



Sun-FTRT/FSW2*
Foot Switch 2 Section
(Sun-FTRTE2 & 5)



Sun-FTRT/FSW3*
Foot Switch 3 Section
(Sun-FTRTE3 & 6)



Sun-FTRT/IV2*
IV Pole Attachment




Sun-INST1*
Installation Service &
Packaging Recycling



Sun-DISPOSAL1*/**
Removal & Disposal
of Old Couch

** Removal from
ground floor only -
weight & removal
restrictions apply,
please enquire

dimensions



anti-bacterial vinyl upholstery colours - please specify when ordering
please see inside back cover for full specification of our vinyl range



The new Fusion Treatment Chair range has been designed using our extensive knowledge to provide a trusted, ergonomic design providing comfort for the patient, ease of access via the split leg design and comfortable patient positioning for the clinician.



Powered tilting seat



Adjustable foot section



Recovery position



Code: **Sun-FTRTE3**
Fusion Treatment Chair - Electric Height Single Foot Section
Powered head section and powered tilting seat

Electric Height Single Foot Section

Code: **Sun-FTRTE1** - Gas assisted head section and fixed seat
Code: **Sun-FTRTE2** - Gas assisted head section and powered tilting seat
Code: **Sun-FTRTE3** - Powered head section and powered tilting seat

Electric Height Split Foot Section

Code: **Sun-FTRTE4** - Gas assisted head section and fixed seat
Code: **Sun-FTRTE5** - Gas assisted head section and powered tilting seat
Code: **Sun-FTRTE6** - Powered head section and powered tilting seat

Hydraulic Height Adjustment with Fixed Seat

Code: **Sun-FTRTH1** - Gas assisted head section and single foot section
Code: **Sun-FTRTH4** - Gas assisted head section and split foot section



3 year warranty on all
moving components



Lifetime warranty
on frame







Normally delivered within
7-10 working days (*UK Only)





Code: **Sun-FTRTH4**
Fusion Treatment Chair - Hydraulic Height Adjustment Fixed Seat
Gas assisted head section and split foot section

Fixed Height Reclining Treatment Chairs

product information



50cm Fixed seat height



Maximum user weight (35st.)



50mm Locking castors for stability & easy movement (Sun-TREA2)



4 Adjustable feet (Sun-TREA1)



Split leg design with locking bars for leg positioning



2 Adjustable/removable padded arm supports



Luxury 60mm thick padded upholstery for comfort




Self locking gas assisted back rest section



Hoist clearance 13cm



Head and foot, up or down -10°/+10° for patient emergencies



Double layer density foam for support & comfort



Anti-bacterial/anti-microbial vinyl



Paper roll holder (supplied as standard)



Easy clean design and materials



Anti-bacterial Vinyl Strength Test
BS EN ISO 13934-1: 2013
BS EN ISO 13973-3: 2000
BS EN ISO 12947-2: 1998



Fire (Flammability) Testing
Flammability BS 7176 Medium Hazard
BS 5852 : 2006 Clause 11 (composites)

accessories


* these accessories can only be fitted when ordered with the couch



Sun-TREA/PHB1
Phlebotomy Arm (each)



Sun-CHP1
Head Pillow



Sun-TREA/IV2*
IV Pole Attachment




Sun-INST1*
Installation Service & Packaging Recycling



Sun-DISPOSAL1**
Removal & Disposal of Old Couch

** Removal from ground floor only - weight and removal restrictions apply, please enquire

dimensions



170cm(L) 65cm(W) 50cm(H)



anti-bacterial vinyl upholstery colours - please specify when ordering

please see inside back cover for full specification of our vinyl range



black grey beige primrose walnut red wine salmon lilac navy mid blue sky blue cool blue green mint white

Fixed Height Reclining Treatment Chairs

The new Fixed Height Reclining Treatment Chair range has been designed using our extensive knowledge to provide a trusted, ergonomic design providing comfort for the patient and comfortable patient positioning for the clinician. The chair features a robust locking bar and gas strut backrest arrangement allowing the patient to rapidly go from a seated position to head down, feet up position in only 3 movements. The split leg design with locking bars is easy to adjust and gives a wide range of positions.



Split leg design with locking bars for leg positioning



In emergency situations for patient recovery chair goes head down -10° and foot up +10°



Code: **Sun-TREA2**
Fixed Height Treatment Chair - with 4 locking castors



2 year warranty on all moving components



Normally delivered within 7-10 working days (*UK Only)



Code: **Sun-TREA1**
Fixed Height Treatment Chair - with adjustable feet



product information



Electric or Hydraulic height adjustment



Maximum user weight (40st.)



Gas assisted head & foot sections



Electric head section & foot sections adjustment option (Sun-FPHBE2 & 3)



Tilting Seat for Trendelenburg (Sun-FPHBE3)



Two fully adjustable phlebotomy arms as standard



50mm Retractable positioning castors



Adjustable height levelling foot




Dual sided foot pedals (Hydraulic models only)



Hoist clearance 13cm



Fully upholstered head & foot sections



Double layer density foam for support & comfort



Luxury 65mm deep contoured upholstery



Anti-bacterial/anti-microbial vinyl



Paper roll holder (supplied as standard)



Easy clean design and materials



Anti-bacterial Vinyl Strength Test
BS EN ISO 13934-1: 2013
BS EN ISO 13973-3: 2000
BS EN ISO 12947-2: 1998



Fire (Flammability) Testing
BS 5852 : 2006 Clause 11 (composites)

accessories

* these accessories can only be fitted when ordered with the couch



Sun-CHP1 Head Pillow




Sun-FPHB/FSW1* Foot Switch (Sun-FPHBE1)



Sun-FPHB/FSW3* Foot Switch 3 Section (Sun-FPHBE2)



Sun-FPHB/FSW4* Foot Switch 4 Section (Sun-FPHBE3)



Sun-FPHB/IV2* IV Pole Attachment




Sun-INST1* Installation Service & Packaging Recycling



Sun-DISPOSAL1** Removal & Disposal of Old Couch

** Removal from ground floor only - weight & removal restrictions apply, please enquire

dimensions



170cm(L) 65cm(W)

HEAVY ITEM

anti-bacterial vinyl upholstery colours - please specify when ordering

please see inside back cover for full specification of our vinyl range



black grey beige primrose walnut red wine salmon lilac navy mid blue sky blue cool blue green mint white

The new Fusion Phlebotomy Chair range has been designed using our extensive knowledge to provide a trusted, ergonomic design providing comfort for the patient and comfortable patient positioning for the clinician.



Easy removal of the phlebotomy arms allows the chair to be used as a standard drop end couch



In emergency situations the chair goes head down and feet up for patient recovery



Code: **Sun-FPHBE3**
Fusion Phlebotomy Chair - Electric Height Adjustment
Electric back and foot sections and tilting seat

Electric Height Adjustment

Code: **Sun-FPHBE1** - Electric Height Adjustment gas assisted head and foot sections
Code: **Sun-FPHBE2** - Electric Height Adjustment electric back and foot sections
Code: **Sun-FPHBE3** - Electric back and foot sections and tilting seat



3 year warranty on all moving components



Lifetime warranty on frame




Normally delivered within 7-10 working days (*UK Only)




Code: **Sun-FPHBH1**
Fusion Phlebotomy Chair - Hydraulic Height Adjustment
Gas assisted head and foot sections




product information




50cm Fixed height seat




225kg
Maximum user weight (35st.)




50mm Locking castors for stability & easy movement (PHLEB2)




4 Adjustable feet (PHLEB1)




Optional padded arm supports




Luxury 60mm thick padded upholstery for comfort




Self locking gas assisted back section adjustment




Hoist clearance 13cm




Head and foot, up or down -10°/+10° for patient emergencies




One fully adjustable phlebotomy arm as standard (can be used either side)




Double layer density foam for support & comfort




Anti-bacterial/anti-microbial vinyl




Paper roll holder (supplied as standard)



Easy clean design and materials



FIRA
Independently tested by FIRA



West Yorkshire Materials Testing Service
Anti-bacterial Vinyl

Fire (Flammability) Testing
Flammability BS 7176 Medium Hazard
BS 5852 : 2006 Clause 11 (composites)

accessories * these accessories can only be fitted when ordered with the couch



Sun-PHLEB/PHB1
Additional Phlebotomy Arm (each)



Sun-CHP1
Head Pillow



Sun-PHLEB/IV2*
IV Pole Attachment



Sun-PHLEB/PAS1
Padded Arm Supports (each)




Sun-INST1*
Installation Service & Packaging Recycling



** Removal from ground floor only - weight & removal restrictions apply, please enquire

Sun-DISPOSAL1**
Removal & Disposal of Old Couch

dimensions




170cm(L) 65cm(W) 50cm(H)




anti-bacterial vinyl upholstery colours - please specify when ordering

please see inside back cover for full specification of our vinyl range




black



grey



beige



primrose




walnut



red wine




salmon




lavender




navy




mid blue




sky blue



cool blue



green



mint



white

The new Fixed Height Reclining Treatment Chair range has been designed using our extensive knowledge to provide a trusted, ergonomic design providing comfort for the patient and comfortable patient positioning for the clinician. The chair features a robust locking bar and gas strut backrest arrangement allowing the patient to rapidly go from a seated position to head down, feet up position in only 2 movements. The leg section with locking bars is easy to adjust and gives a wide range of positions.



In emergency situations for patient recovery chair goes head down -10° and foot up +10°



Code: **Sun-PHLEB1**
Fixed Height Phlebotomy Chair - with adjustable feet



Code: **Sun-PHLEB2**
Fixed Height Phlebotomy Chair - with 4 locking castors




2 year warranty on all moving components




Normally delivered within 7-10 working days (*UK Only)




product information




4 Cabinet finishes




Maximum user weight (26st.)




Multi-position self locking backrest



Luxury 65mm deep contoured upholstery




Anti-bacterial/anti-microbial vinyl




Paper roll holder (supplied as standard)




2 Pull-out work surfaces




Cupboards include pull-out work surface & basket, shelf and lock




Locking doors supplied with 2 keys




Durable silver epoxy coated frame




Cabinets supplied fully assembled




Easy clean design and materials



Anti-bacterial Vinyl Strength Test
BS EN ISO 13934-1: 2013
BS EN ISO 13973-3: 2000
BS EN ISO 12947-2: 1998



Fire (Flammability) Testing
Flammability BS 7176 Medium Hazard
BS 5852 : 2006 Clause 11 (composites)



FIRA Standard 6250 2005
Specification:
Adhesion Performance
General/Severe

accessories * these accessories can only be fitted when ordered with the couch



Sun-CBH1* Breathing Hole



Sun-CHP1 Head Pillow



Sun-CDT1 Drawer Insert (each)



Sun-LEG1* Leg Supports



Sun-UPH1* Special Upholstery Colour



Sun-INST1* Installation Service & Packaging Recycling



Sun-DISPOSAL1*/** Removal & Disposal of Old Couch



Sun-STP1 Single Tier Step (oval)




Sun-STP2 Double Tier Step (oval)




Sun-STP3 Single Tier Step (rectangular)

** Removal from ground floor only - weight and removal restrictions apply, please enquire

dimensions



83cm(H) 193cm(L) 60cm(W)
with paper roll holder



anti-bacterial vinyl upholstery colours - please specify when ordering

please see inside back cover for full specification of our vinyl range



black grey beige primrose walnut red wine salmon lilac navy mid blue sky blue cool blue green mint white

Specialist 2 Section Treatment Couch Systems

Combining functionality with storage and comfort, Specialist Couches assist in making examination and treatment procedures easier and more efficient.



Code: **Sun-CS1**
Specialist treatment couch system with 2 drawerline units, 1 drawer pack, 2 section upholstery with optional head pillow and leg supports



White (high gloss) cabinet finish



Titanium (high gloss) cabinet finish















Beech effect cabinet finish



Maple effect cabinet finish



product information

 4 Cabinet finishes	 Maximum user weight (26st.)	 Multi-position self locking backrest	 Luxury 65mm deep contoured upholstery	 Anti-bacterial/anti-microbial vinyl
 Paper roll holder (supplied as standard)	 2 Pull-out work surfaces	 Cupboards include pull-out work surface & basket, shelf and lock	 Locking doors supplied with 2 keys	 Durable silver epoxy coated frame
 Cabinets supplied fully assembled	 Easy clean design and materials	 Anti-bacterial Vinyl Strength Test BS EN ISO 13934-1: 2013 BS EN ISO 13973-3: 2000 BS EN ISO 12947-2: 1998	 Fire (Flammability) Testing Flammability BS 7176 Medium Hazard BS 5852 : 2006 Clause 11 (composites)	 FIRA Standard 6250 2005 Specification: Adhesion Performance General/Severe

accessories * these accessories can only be fitted when ordered with the couch

 Sun-CBH1* Breathing Hole	 Sun-CHP1 Head Pillow	 Sun-CDT1 Drawer Insert (each)	 Sun-LEG1* Leg Supports	 Sun-UPH1* Special Upholstery Colour
 Sun-INST1* Installation Service & Packaging Recycling	 Sun-DISPOSAL1*/** Removal & Disposal of Old Couch	 Sun-STP1 Single Tier Step (oval)	 Sun-STP2 Double Tier Step (oval)	 Sun-STP3 Single Tier Step (rectangular)

** Removal from ground floor only - weight and removal restrictions apply, please enquire

dimensions


















83cm(H) 193cm(L) 60cm(W)
with paper roll holder

HEAVY ITEM

anti-bacterial vinyl upholstery colours - please specify when ordering

please see inside back cover for full specification of our vinyl range

														
black	grey	beige	primrose	walnut	red wine	salmon	lavender	navy	mid blue	sky blue	cool blue	green	mint	white

Specialist 2 Section Treatment Couch Systems

Combining functionality with storage and comfort, Specialist Couches assist in making examination and treatment procedures easier and more efficient.



Code: **Sun-CS3**
Specialist treatment couch system with 1 drawerline unit, 1 drawer pack, 2 section upholstery with optional head pillow and leg supports



White (high gloss) cabinet finish



Titanium (high gloss) cabinet finish




Beech effect cabinet finish




Maple effect cabinet finish


product information




4 Cabinet finishes




Maximum user weight (26st.)




Multi-position self locking backrest




Luxury 65mm deep contoured upholstery




Anti-bacterial/anti-microbial vinyl




Paper roll holder (supplied as standard)




2 Pull-out work surfaces




Cupboards include pull-out work surface & basket, shelf and lock




Locking doors supplied with 2 keys




Durable silver epoxy coated frame




Cabinets supplied fully assembled




Easy clean design and materials



Anti-bacterial Vinyl Strength Test
BS EN ISO 13934-1: 2013
BS EN ISO 13973-3: 2000
BS EN ISO 12947-2: 1998



Fire (Flammability) Testing
Flammability BS 7176 Medium Hazard
BS 5852 : 2006 Clause 11 (composites)



FIRA Standard 6250 2005
Specification:
Adhesion Performance
General/Severe

accessories * these accessories can only be fitted when ordered with the couch



Sun-CBH1* Breathing Hole



Sun-CHP1 Head Pillow



Sun-CDT1 Drawer Insert (each)



Sun-LEG1* Leg Supports



Sun-UPH1* Special Upholstery Colour



Sun-INST1* Installation Service & Packaging Recycling



Sun-DISPOSAL1** Removal & Disposal of Old Couch



Sun-STP1 Single Tier Step (oval)




Sun-STP2 Double Tier Step (oval)



Sun-STP3 Single Tier Step (rectangular)

** Removal from ground floor only - weight and removal restrictions apply, please enquire

dimensions



anti-bacterial vinyl upholstery colours - please specify when ordering

please see inside back cover for full specification of our vinyl range

black	grey	beige	primrose	walnut	red wine	salmon	lavender	navy	mid blue	sky blue	cool blue	green	mint	white
-------	------	-------	----------	--------	----------	--------	----------	------	----------	----------	-----------	-------	------	-------

Specialist 3 Section Treatment Couch Systems

Combining functionality with storage and comfort, Specialist Couches assist in making examination and treatment procedures easier and more efficient.



Code: **Sun-CS1/3S**
Specialist treatment couch system with 2 drawerline units, 1 drawer pack, 3 section upholstery with optional head pillow



White (high gloss) cabinet finish



Titanium (high gloss) cabinet finish















Beech effect cabinet finish



Maple effect cabinet finish

product information


 4 Cabinet finishes	 Maximum user weight (26st.)	 Multi-position self locking backrest	 Luxury 65mm deep contoured upholstery	 Anti-bacterial/anti-microbial vinyl
 Paper roll holder (supplied as standard)	 2 Pull-out work surfaces	 Cupboards include pull-out work surface & basket, shelf and lock	 Locking doors supplied with 2 keys	 Durable silver epoxy coated frame
 Cabinets supplied fully assembled	 Easy clean design and materials	 Anti-bacterial Vinyl Strength Test BS EN ISO 13934-1: 2013 BS EN ISO 13973-3: 2000 BS EN ISO 12947-2: 1998	 Fire (Flammability) Testing Flammability BS 7176 Medium Hazard BS 5852 : 2006 Clause 11 (composites)	 FIRA Standard 6250 2005 Specification: Adhesion Performance General/Severe

accessories * these accessories can only be fitted when ordered with the couch

 Sun-CBH1* Breathing Hole	 Sun-CHP1 Head Pillow	 Sun-CDT1 Drawer Insert (each)	 Sun-LEG1* Leg Supports	 Sun-UPH1* Special Upholstery Colour
 Sun-INST1* Installation Service & Packaging Recycling	 Sun-DISPOSAL1*/** Removal & Disposal of Old Couch	 Sun-STP1 Single Tier Step (oval)	 Sun-STP2 Double Tier Step (oval)	 Sun-STP3 Single Tier Step (rectangular)
















** Removal from ground floor only - weight and removal restrictions apply, please enquire

dimensions



anti-bacterial vinyl upholstery colours - please specify when ordering

please see inside back cover for full specification of our vinyl range

														
black	grey	beige	primrose	walnut	red wine	salmon	lavender	navy	mid blue	sky blue	cool blue	green	mint	white

Specialist 3 Section Treatment Couch Systems

Combining functionality with storage and comfort, Specialist Couches assist in making examination and treatment procedures easier and more efficient.



Code: **Sun-CS3/3S**
Specialist treatment couch system with 1 drawerline unit, 1 drawer pack, 3 section upholstery with optional head pillow



White (high gloss) cabinet finish



Titanium (high gloss) cabinet finish



Beech effect cabinet finish



Maple effect cabinet finish

product information



4 Cabinet finishes



Maximum user weight (26st.)



Multi-position self locking backrest



Luxury 65mm deep contoured upholstery



Anti-bacterial/anti-microbial vinyl



Paper roll holder (supplied as standard)



Durable silver epoxy coated frame



Cupboards include pull-out work surface & basket, shelf and lock



60cm diameter legs for stability with adjustable feet for uneven floors



Choice of integrated storage solutions



Optional couch steps (see pages 42-43)



Easy clean design and materials



Anti-bacterial Vinyl Strength Test
BS EN ISO 13934-1: 2013
BS EN ISO 13973-3: 2000
BS EN ISO 12947-2: 1998



Fire (Flammability) Testing
Flammability BS 7176 Medium Hazard
BS 5852 : 2006 Clause 11 (composites)



FIRA Standard 6250 2005
Specification:
Adhesion Performance
General/Severe

accessories * these accessories can only be fitted when ordered with the couch



Sun-CBH1* Breathing Hole



Sun-CHP1 Head Pillow



Sun-MCD1 Single Drawer Unit



Sun-CDT2 Drawer Insert (each)



Sun-CP52 Pull-out Work Surface



Sun-CMU1 Mobile Storage Unit (with locking door)



Sun-PB1* Pull-out Basket for Storage Unit



Sun-LEG1* Leg Supports (pair)




Sun-INST1* Installation Service & Packaging Recycling



Sun-DISPOSAL1** Removal & Disposal of Old Couch

** Removal from ground floor only - weight and removal restrictions apply, please enquire

dimensions



0° +65°



HEAVY ITEM

79cm(H) 193cm(L) 60cm(W)
with paper roll holder




Mobile Storage Unit

58cm(H) 50cm(W) 60cm(D)

anti-bacterial vinyl upholstery colours - please specify when ordering

please see inside back cover for full specification of our vinyl range



black grey beige primrose walnut red wine salmon lilac navy mid blue sky blue cool blue green mint white

Practitioner Deluxe Examination Couches

Stylish, practical and comfortable, Practitioner Couches make examination and treatment procedures easier, combined with contemporary aesthetics.



Code: **Sun-MEC1**
Practitioner deluxe examination couch
with optional pull-out work surface, drawer units,
mobile storage unit, leg supports, step and head pillow



White (high gloss) cabinet finish Titanium (high gloss) cabinet finish



Beech effect cabinet finish



Maple effect cabinet finish

product information



Maximum user weight (26st.)



Multi-position self locking backrest



Single layer density foam for support & comfort



Luxury 65mm deep contoured upholstery



Anti-bacterial/anti-microbial vinyl



Paper roll holder (supplied as standard)



Durable silver epoxy coated frame



Stylish 60cm diameter legs for stability



4 non marking leg protectors



Special order upholstery colours available



Easy clean design and materials



Anti-bacterial Vinyl Strength Test
BS EN ISO 13934-1: 2013
BS EN ISO 13973-3: 2000
BS EN ISO 12947-2: 1998



Fire (Flammability) Testing
Flammability BS 7176 Medium Hazard
BS 5852 : 2006 Clause 11 (composites)

accessories * these accessories can only be fitted when ordered with the couch



Sun-CBH1*
Breathing Hole



Sun-CHP1
Head Pillow



Sun-LEG1*
Leg Supports (pair)



Sun-INST1*
Installation Service & Packaging Recycling



Sun-DISPOSAL1*/**
Removal & Disposal of Old Couch



Sun-MPT51
Vista 15 Under Couch Trolley



Sun-MPT52
Vista 15 Under Couch Trolley



Sun-MPT53
Vista 15 Under Couch Trolley



Vista 15 Under Couch Trolleys (see pages 56-57)



Couch Steps (see pages 42-43)

** Removal from ground floor only - weight and removal restrictions apply, please enquire

dimensions



HEAVY ITEM

79cm(H) 193cm(L) 60cm(W)
with paper roll holder

anti-bacterial vinyl upholstery colours - please specify when ordering

please see inside back cover for full specification of our vinyl range

black	grey	beige	primrose	walnut	red wine	salmon	lavender	navy	mid blue	sky blue	cool blue	green	mint	white

Colenso Examination Couches

Simple, strong and durable, Colenso Couches are designed for everyday general healthcare, providing comfort with contemporary design.



Code: **Sun-MEC2**
Colenso examination couch with optional head pillow and couch step



Code: **Sun-MPT51**
Vista 15 under couch trolley
4 single depth trays



Code: **Sun-MPT52**
Vista 15 under couch trolley
2 single/1 double depth trays



Code: **Sun-MPT53**
Vista 15 under couch trolley
2 double depth trays

product information



Low height
only 50cm high



Maximum user
weight (26st.)



Multi-position
self locking backrest



Single layer density foam
for support & comfort



Luxury 65mm deep
contoured upholstery



Anti-bacterial/
anti-microbial vinyl



Paper roll holder
(supplied as standard)



Durable silver epoxy
coated frame



60cm diameter
legs for stability



4 adjustable feet
for uneven floors



Special order upholstery
colours available



Easy clean design
and materials



Anti-bacterial Vinyl Strength Test
BS EN ISO 13934-1: 2013
BS EN ISO 13973-3: 2000
BS EN ISO 12947-2: 1998



Fire (Flammability) Testing
Flammability BS 7176 Medium Hazard
BS 5852 : 2006 Clause 11
(composites)

accessories * these accessories can only be fitted when ordered with the couch



Sun-CBH1*
Breathing Hole



Sun-CHP1
Head Pillow



Sun-INST1*
Installation Service &
Packaging Recycling



Sun-DISPOSAL1**
Removal & Disposal
of Old Couch



Couch Steps
(see pages 42-43)

** Removal from ground floor only - weight and removal restrictions apply, please enquire

dimensions


















HEAVY ITEM

50cm(H) 193cm(L) 60cm(W)
with paper roll holder

anti-bacterial vinyl upholstery colours - please specify when ordering

please see inside back cover for full specification of our vinyl range

														
black	grey	beige	primrose	walnut	red wine	salmon	lavender	navy	mid blue	sky blue	cool blue	green	mint	white

Alberti Rest Couches

Designed for resting, recovery or where a reduced height application is required, Alberti Couches are strong, durable, easy to clean while being exceptionally comfortable.



Code: **Sun-RC1**
Alberti rest couch with optional head pillow



product information



Maximum user weight (24st.)



Lightweight design only 19kg (excludes carry case)



Multi-position self locking backrest



Breathing hole and plug



Adjustable height range 64cm - 80cm



Head cradle



Carrying case



Quick and easy to erect



2 year warranty




Delivery normally within 3-5 working days *

* Please contact us regarding availability for quantities of 10 or more.

dimensions



vinyl upholstery colours - please specify when ordering



product information



Maximum user weight (26st.)



Anti-slip surfaces



Anti-tipping designs



Rubber feet prevent movement



Unique oval design (STP1 & STP2 only)



Easy clean design and materials



2 year warranty




British made



Delivery normally within 3-5 working days *

* Please contact us regarding availability for quantities of 10 or more.

dimensions




Portable Couches

Lightweight but strong, durable and easy to clean, Portable Couches are ideal for occasional or temporary use as well as for mobile treatment and procedures.



Code: **Sun-PC1**
Portable couch with breathing hole and plug, head cradle and carrying case



Couch Steps

Designed for using with fixed height couches, couch steps are the perfect solution for assisting patients to raise or lower themselves from or to the floor.



Code: **Sun-STP1**
Single tier couch step (oval)






Code: **Sun-STP2**
Double tier couch step (oval)

Code: **Sun-STP3**
Single tier couch step (rectangular)

stable
adjustable
support

Phlebotomy Chairs, Operator Chairs & Examination Stools

Hygienically designed and with patient comfort in mind, phlebotomy chairs are a simple solution for taking blood, whilst operator chairs and examination stools provide support for clinicians.

	PAGE
	Reclining Phlebotomy Chair 45-46
	Phlebotomy/Treatment Chairs 47-48
	Operator Chairs 49-50
	Examination Stools 51-52
	Examination Stools with Back Rests 53-54

product information



50cm Fixed height seat



225kg
Maximum user weight (35st.)



50mm Locking castors for stability & easy movement (PHLEB2)



4 Adjustable feet (PHLEB1)



Optional padded arm supports



Luxury 60mm thick padded upholstery for comfort



Self locking gas assisted back section adjustment



Hoist clearance 13cm



Head and foot, up or down -10°/+10° for patient emergencies



One fully adjustable phlebotomy arm as standard (can be used either side)



Double layer density foam for support & comfort



Anti-bacterial/anti-microbial vinyl



Paper roll holder (supplied as standard)



Easy clean design and materials



FIRA
Independently tested by FIRA



West Yorkshire Materials Testing Service
Anti-bacterial Vinyl

Fire (Flammability) Testing
Flammability BS 7176 Medium Hazard
BS 5852 : 2006 Clause 11 (composites)

Anti-bacterial Vinyl Strength Test
BS EN ISO 13934-1: 2013
BS EN ISO 13973-3: 2000
BS EN ISO 12947-2: 1998

accessories * these accessories can only be fitted when ordered with the couch



Sun-PHLEB/PHB1
Additional Phlebotomy Arm (each)



Sun-CHP1
Head Pillow



Sun-PHLEB/IV2*
IV Pole Attachment



Sun-PHLEB/PAS1
Padded Arm Supports (each)




Sun-INST1*
Installation Service & Packaging Recycling




Sun-DISPOSAL1**
Removal & Disposal of Old Couch

** Removal from ground floor only - weight & removal restrictions apply, please enquire

dimensions



170cm(L) 65cm(W) 50cm(H)



anti-bacterial vinyl upholstery colours - please specify when ordering

please see inside back cover for full specification of our vinyl range



black grey beige primrose walnut red wine salmon lilac navy mid blue sky blue cool blue green mint white

The new Fixed Height Reclining Treatment Chair range has been designed using our extensive knowledge to provide a trusted, ergonomic design providing comfort for the patient and comfortable patient positioning for the clinician. The chair features a robust locking bar and gas strut backrest arrangement allowing the patient to rapidly go from a seated position to head down, feet up position in only 2 movements. The leg section with locking bars is easy to adjust and gives a wide range of positions.



Code: **Sun-PHLEB1**
Fixed Height Phlebotomy Chair - with adjustable feet



2 Year Warranty
2 Year Warranty
2 year warranty on all moving components



Normally delivered within 7-10 working days (*UK Only)



Code: **Sun-PHLEB2**
Fixed Height Phlebotomy Chair - with 4 locking castors

product information

110kg

Maximum user weight (17st.)

Lightweight design only 8kgs

Multi-position arm rest, made from a hygienic PU moulding

Arm rest can be positioned either side (2 brackets supplied)

Ergonomically designed moulded seat & back (Sun-PCHA only)

Anti-bacterial/anti-microbial vinyl (Sun-PCHA/VYL only)

40mm Deep foam for support and comfort (Sun-PCHA/VYL only)

Easy clean design and materials

1 year warranty

Supplied fully assembled

Independently tested by FIRA

Independently tested

Independently tested by FIRA

Anti-bacterial Vinyl Strength Test
BS EN ISO 13934-1: 2013
BS EN ISO 13973-3: 2000
BS EN ISO 12947-2: 1998

West Yorkshire Materials Testing Service
Anti-bacterial Vinyl

Fire (Flammability) Testing
Flammability BS 7176 Medium Hazard
BS 5852 : 2006 Clause 11 (composites)

British made

CE marked

accessories

Sun-PCHA/Arm
Additional Arm Rest

dimensions

83cm(H) 70cm(W) 56cm(D)
Including arm brackets (excluding arm)
Seat Height 43.5cm(H)

Arm Rest

35cm(D)

49cm(H)

Maximum user weight 110kg (17st.)

12cm(D)

35cm(H)

moulded seat & back colours - please specify when ordering

Black

Blue

Grey

Beige

anti-bacterial vinyl upholstery colours - please specify when ordering

please see inside back cover for full specification of our vinyl range

black

grey

beige

primrose

walnut

red wine

salmon

lavender

navy

mid blue

sky blue

cool blue

green

mint

white

Durable, stable and ergonomically designed, phlebotomy chairs are easy to use and adjust and can be used left, right or dual sided with an additional arm, whilst also being easy to clean.



Code: Sun-PCHA
Phlebotomy chair
Moulded seat and back, one arm rest



CODE: Sun-PCHA/VYL
Phlebotomy chair
Vinyl upholstery, one arm rest

product information

 110kg Maximum user weight (17st.)	 Anti-bacterial/anti-microbial vinyl	 Comfortable padded seat and back	 40mm Deep foam for support and comfort	 Hygienic rounded corners
 Strong and durable ergonomic arms (CHA2 & CHA4)	 Adjustable height chrome foot ring provides leg support (CHA3, CHA4 & CHA5)	 Gas-assisted height adjustment	 Five castor base for stability and manoeuvrability	 Glides reduce lateral movement (CHA5 only)
 HEIGHT RANGE 42-57cm Sun-CHA1	 HEIGHT RANGE 42-57cm Sun-CHA2	 HEIGHT RANGE 48-68cm Sun-CHA3	 HEIGHT RANGE 48-68cm Sun-CHA4	 HEIGHT RANGE 55-81cm Sun-CHA5
 Easy clean design and materials	 Anti-bacterial Vinyl Strength Test BS EN ISO 13934-1: 2013 BS EN ISO 13973-3: 2000 BS EN ISO 12947-2: 1998	 Fire (Flammability) Testing Flammability BS 7176 Medium Hazard BS 5852 : 2006 Clause 11 (composites)	 2 year warranty	 British made
 Delivery normally within 3-5 working days *				

* Please contact us regarding availability for quantities of 10 or more.

dimensions

 Sun-CHA1 Seat 46cm(W) 43cm(D) Base 60cm Ø	 Sun-CHA2 Seat 46cm(W) 43cm(D) Width with arms 60cm Base 60cm Ø	 Sun-CHA3 Seat 46cm(W) 43cm(D) Foot ring 46cm Ø Base 60cm Ø	 Sun-CHA4 Seat 46cm(W) 43cm(D) Width with arms 60cm Foot ring 46cm Ø Base 60cm Ø	 Sun-CHA5 Seat 46cm(W) 43cm(D) Foot ring 46cm Ø Base 60cm Ø
---	---	---	---	---

anti-bacterial vinyl upholstery colours - please specify when ordering

please see inside back cover for full specification of our vinyl range

black	grey	beige	primrose	walnut	red wine	salmon	lavender	navy	mid blue	sky blue	cool blue	green	mint	white

Operator Chairs

Adjustable for support and comfort, Operator Chairs are designed for everyday general healthcare use, therefore are comfortable, strong, durable and easy to clean.



Code: Sun-CHA1 - Gas-lift chair



Code: Sun-CHA2 - Gas-lift chair with arms



Code: Sun-CHA3 - Gas-lift chair & foot ring



Code: Sun-CHA4 - Gas-lift chair with arms & foot ring



Code: Sun-CHA5 - High level gas-lift chair & foot ring

product information

Maximum user weight (17st.)

Anti-bacterial/ anti-microbial vinyl

Luxury 65mm deep upholstery

2 layer foam for support and comfort

Adjustable height chrome foot ring (STO3 & STO3/GLIDES)

Gas-assisted height adjustment

HEIGHT RANGE
42-58cm
Sun-STO2

HEIGHT RANGE
39-55cm
Sun-STO2/GLIDES

HEIGHT RANGE
51-69cm
Sun-STO3

HEIGHT RANGE
48-66cm
Sun-STO3/GLIDES

Five castor base for stability and manoeuvrability

Glides reduce lateral movement (STO2/GLIDES & STO3/GLIDES only)

Easy clean design and materials

Anti-bacterial Vinyl Strength Test
BS EN ISO 13934-1: 2013
BS EN ISO 13973-3: 2000
BS EN ISO 12947-2: 1998

Fire (Flammability) Testing
BS 5852 : 2006 Clause 11 (composites)

Anti-bacterial Vinyl


2 year warranty

British made

Delivery normally within 3-5 working days *


* Please contact us regarding availability for quantities of 10 or more.

dimensions



anti-bacterial vinyl upholstery colours - please specify when ordering

please see inside back cover for full specification of our vinyl range



Examination Stools

Designed for everyday general healthcare use, Examination Stools can be static or mobile, are simple to adjust and comfortable, while being strong, durable and easy to clean.



Code: **Sun-STO2**
Gas-lift stool with 5 castor base



Code: **Sun-STO2/GLIDES**
Gas-lift stool with 5 glides



Code: **Sun-STO3**
Gas-lift stool & foot ring with 5 castor base




Code: **Sun-STO3/GLIDES**
Gas-lift stool & foot ring with 5 glides


product information




Maximum user weight (17st.)




Anti-bacterial/anti-microbial vinyl



Comfortable padded back rest



Luxury 65mm deep seat with upholstered back rest




2 layer foam for support and comfort



Adjustable height chrome foot ring (Sun-STO5 & Sun-STO6)



Gas-assisted height adjustment




Sun-STO4



Sun-STO5



Sun-STO6



Five castor base for stability & manoeuvrability (Sun-STO4 & Sun-STO6)




Glides reduce lateral movement (Sun-STO5)



Easy clean design and materials




Anti-bacterial Vinyl Strength Test
BS EN ISO 13934-1: 2013
BS EN ISO 13973-3: 2000
BS EN ISO 12947-2: 1998




Fire (Flammability) Testing
Flammability BS 7176 Medium Hazard
BS 5852 : 2006 Clause 11 (composites)



2 year warranty




British made



Delivery normally within 3-5 working days *

* Please contact us regarding availability for quantities of 10 or more.

dimensions



anti-bacterial vinyl upholstery colours - please specify when ordering

please see inside back cover for full specification of our vinyl range



Examination Stools with Back Rests

Designed for everyday general healthcare use, Examination Stools with back rests can be static or mobile, are simple to adjust and comfortable, while being strong, durable and easy to clean.



Code: **Sun-STO4**
Gas -lift stool with back rest



Code: **Sun-STO6**
Gas -lift stool with back rest and foot ring









Code: **Sun-STO5**
High level gas -lift stool with back rest, foot ring and glides

visual
accessible
storage

Vista® Trolleys & Storage

Flexible, modular and colourful, Vista incorporates a mobile or fixed storage solution for every healthcare environment, using premium materials which are all easy to clean. A lockable option is available on most of the range.

	PAGE
 Vista Trolleys	56-71
 Vista Clinical Trolleys	72-79
 Vista Colour Concept Clinical Trolleys	80-81
 Vista Action Station Trolleys	82-85
 Vista Storage Modules	86-89
 Vista Base & Larder Cabinets	90-91

product information



Tilt & locking tray runner system



Trays integrate with all other Vista products



100mm castors with front brakes (Vista 15 only)



Optional tamper resistant locking door (Vista 15 only)



Extra low level (68cm), ideal for under couches (Vista 15 only)



Stainless steel handle (Vista 15 only)



Removable top shelf with integrated upstand (Vista 15 only)



Easy clean design and materials



2 year warranty



British Made

accessories

* These accessories can only be fitted when ordered with the trolley



Sun-LC3A *
4 Antistatic Castors (front 2 braked)



Sun-BUF2 *
4 Protective Buffers (Vista 10 Trolleys)



Sun-BUF3 *
4 Protective Buffers (Vista 15 Trolleys)



Sun-MPT51 *
Folding Shelf



Sun-SBOX1 */**
Sharps Box Bracket (Daniels - All Sizes)



Sun-SBOX2 */**
Sharps Box Bracket (Frontier) Upto 1 Litre



Sun-SBOX3 */**
Sharps Box Bracket (Frontier) 2 & 3 Litre



Sun-MPT-AGH1 */**
Alcohol Gel Holder



Sun-GBH1 */**
Glove Box Holder



Sun-MPT/Rail1 *
Medi-rail



Sun-MPT/BIN *
2.5 Litre Bin (Requires Medi-rail)



Sun-MPT51/LD */***
Tamper Resistant Locking Door (Vista 15 Trolleys)



Sun-DIV1
Single Depth Tray Dividers (per set) Vista 10 & 15



Sun-DIV4
Double Depth Tray Dividers (per set) Vista 15



Sun-MPT1/CV
Trolley Cover (Vista 10 Trolleys)





Sun-MPT51/CV
Trolley Cover (Vista 15 Trolleys)

** Includes Medi-rail

*** Locking door option includes a secure top and back panel


dimensions



Single Depth Trays
8cm(H) 47cm(W) 42cm(D)
Double Depth Trays
16cm(H) 47cm(W) 42cm(D)

tray colours

please specify when ordering
(Tray labels and covers included)



Vista® 10 & 15 Trolleys

Vista 10 is an economy priced trolley without any compromise on quality, while Vista 15 is designed to fit under a couch, both trolleys are functional, durable and easy to clean.



Code: **Sun-TRO1**
Vista 10 Trolley - 1 tray & 1 shelf



Code: **Sun-MPT51**
Vista 15 Trolley - 4 single depth trays



Code: **Sun-MPT51/LD**
Locking door for Vista 15 trolley



Code: **Sun-MPT52**
Vista 15 Trolley - 2 single/1 double depth trays



Code: **Sun-MPT53**
Vista 15 Trolley - 2 double depth trays

product information



Tilt & locking tray runner system



Trays integrate with all other Vista products



Narrow width (39cm) compact design



100mm castors with front brakes



Optional tamper resistant locking door



Removable top shelf with integrated upstand



Easy clean design and materials



2 year warranty



Supplied fully assembled



British made

accessories * These accessories can only be fitted when ordered with the trolley



Sun-LC3A *
4 Antistatic Castors (front 2 braked)



Sun-BUF2 *
4 Protective Buffers



Sun-MPTS1 *
Folding Shelf



Sun-SBOX1 */**
Sharps Box Bracket (Daniels - All Sizes)



Sun-SBOX2 */**
Sharps Box Bracket (Frontier) Upto 1 Litre



Sun-SBOX3 */**
Sharps Box Bracket (Frontier) 2 & 3 Litre



Sun-MPT/AGH1 */**
Alcohol Gel Holder



Sun-GBH1 */**
Glove Box Holder



Sun-MPT/Rail1 *
Medi-rail



Sun-MPT/BIN *
2.5 Litre Bin (Requires Medi-rail)



Sun-MPT5/LD */***
Tamper Resistant Locking Door



Sun-DIV2
Double Depth Tray Dividers (per set)



Sun-DIV3
Single Depth Tray Dividers (per set)



Sun-MPT5/CV
Trolley Cover

** Includes Medi-rail *** Locking door option includes a secure top and back panel

dimensions



101cm(H) 39cm(W) 52cm(D)

Single Depth Trays
8cm(H) 31cm(W) 42cm(D)
Double Depth Trays
16cm(H) 31cm(W) 42cm(D)

tray colours

please specify when ordering
(Tray labels and covers included)



clear blue green

Vista® 30 Trolleys

Designed for environments with limited space, Vista 30 provides standard working height with compact design, optimum storage, and is functional, durable and easy to clean.



Code: Sun-MPT10
Vista 30 Trolley - 6 single depth trays



Code: Sun-MPT5/LD
Locking door for Vista 30 trolley



Code: Sun-MPT5
Vista 30 Trolley - 2 single/2 double depth trays



Code: Sun-MPT5/LD
Locking door for Vista 30 trolley



Code: Sun-MPT6
Vista 30 Trolley - 3 double depth trays



Code: Sun-MPT5/LD
Locking door for Vista 30 trolley

product information



Tilt & locking tray runner system



Trays integrate with all other Vista products



Low level (80cm), ideal for under worktops



100mm castors with front brakes



Optional tamper resistant locking door



Removable top shelf with integrated upstand



Stainless steel handle



Easy clean design and materials



2 year warranty



British made

accessories * These accessories can only be fitted when ordered with the trolley



Sun-LC3A *
4 Antistatic Castors (front 2 braked)



Sun-BUF3 *
4 Protective Buffers



Sun-MPTS1 *
Folding Shelf



Sun-SBOX1 */**
Sharps Box Bracket (Daniels - All Sizes)



Sun-SBOX2 */**
Sharps Box Bracket (Frontier) Upto 1 Litre



Sun-SBOX3 */**
Sharps Box Bracket (Frontier) 2 & 3 Litre



Sun-MPT/AGH1 */**
Alcohol Gel Holder



Sun-GBH1 */**
Glove Box Holder



Sun-MPT/Rail1 *
Medi-rail



Sun-MPT/BIN *
2.5 Litre Bin (Requires Medi-rail)



Sun-MPT40/LD */***
Tamper Resistant Locking Door



Sun-DIV1
Single Depth Tray Dividers (per set)




Sun-DIV4
Double Depth Tray Dividers (per set)



Sun-MPT40/CV
Trolley Cover

** Includes Medi-rail *** Locking door option includes a secure top and back panel


dimensions



Single Depth Trays
8cm(H) 47cm(W) 42cm(D)
Double Depth Trays
16cm(H) 47cm(W) 42cm(D)

tray colours

please specify when ordering
(Tray labels and covers included)



clear blue green

Vista® 40 Trolleys

Designed to fit under a worktop, Vista 40 can be stored neatly away, while providing optimum number of trays for storage, and is functional, durable and easy to clean.



Code: **Sun-MPT40** Vista 40 Trolley
5 single depth trays



Code: **Sun-MPT40/LD**
Locking door for Vista 40 trolley



Code: **Sun-MPT41** Vista 40 Trolley
3 single/1 double depth trays



Code: **Sun-MPT40/LD**
Locking door for Vista 40 trolley



Code: **Sun-MPT42** Vista 40 Trolley
1 single/2 double depth trays



Code: **Sun-MPT40/LD**
Locking door for Vista 40 trolley

product information



Tilt & locking tray runner system



Trays integrate with all other Vista products



100mm castors with front brakes



Optional tamper resistant locking door



Removable top shelf with integrated upstand



Easy clean design and materials



2 year warranty



Supplied fully assembled



British made

accessories * These accessories can only be fitted when ordered with the trolley



Sun-LC3A *
4 Antistatic Castors (front 2 braked)



Sun-BUF2 *
4 Protective Buffers



Sun-MPT/OXY1 */**
Oxygen Bottle Holder (Size D)



Sun-MPT/LAP1 */**
Adjustable Height Laptop Holder



Sun-MPT/CHART1 */**
Adjustable Height Chart Holder



Sun-MPT/DRIP1 */**
Adjustable Height Drip Pole



Sun-MPTS1 *
Folding Shelf



Sun-SBOX1 */**
Sharps Box Bracket (Daniels - All Sizes)



Sun-SBOX2 */**
Sharps Box Bracket (Frontier) Upto 1 Litre



Sun-SBOX3 */**
Sharps Box Bracket (Frontier) 2 & 3 Litre



Sun-MPT/AGH1 */**
Alcohol Gel Holder



Sun-GBH1 */**
Glove Box Holder



Sun-MPT/Rail1 *
Medi-rail



Sun-MPT/BIN *
2.5 Litre Bin (Requires Medi-rail)



Sun-MPT1/LD */**
Tamper Resistant Locking Door



Sun-DIV1
Single Depth Tray Dividers (per set)



Sun-DIV4
Double Depth Tray Dividers (per set)



Sun-MPT1/CV
Trolley Cover

** Includes Medi-rail *** Locking door option includes a secure top and back panel

dimensions



101cm(H) 56cm(W) 52cm(D)

Single Depth Trays
8cm(H) 47cm(W) 42cm(D)
Double Depth Trays
16cm(H) 47cm(W) 42cm(D)

tray colours

please specify when ordering
(Tray labels and covers included)



clear blue green

Vista® 50 Trolleys

Suitable for everyday general healthcare use, Vista 50 provides standard working height with generous working space, optimum tray storage, and is functional, durable and easy to clean.



Code: Sun-MPT1
Vista 50 Trolley - 6 single depth trays



Code: Sun-MPT1/LD
Locking door for Vista 50 trolley



Code: Sun-MPT2
Vista 50 Trolley - 4 single/1 double depth trays



Code: Sun-MPT1/LD
Locking door for Vista 50 trolley



Code: Sun-MPT13
Vista 50 Trolley - 3 double depth trays



Code: Sun-MPT1/LD
Locking door for Vista 50 trolley

product information



Tilt & locking tray runner system



Trays integrate with all other Vista products



100mm castors with front brakes



Optional tamper resistant locking door



Removable top shelf with integrated upstand



Easy clean design and materials



2 year warranty



Supplied fully assembled



British made

accessories * These accessories can only be fitted when ordered with the trolley



Sun-LC3A *
4 Antistatic Castors (front 2 braked)



Sun-BUF2 *
4 Protective Buffers



Sun-MPT/OXY1 */**
Oxygen Bottle Holder (Size D)



Sun-MPT/LAP1 */**
Adjustable Height Laptop Holder



Sun-MPT/CHART1 */**
Adjustable Height Chart Holder



Sun-MPT/DRIP1 */**
Adjustable Height Drip Pole



Sun-MPTS1 *
Folding Shelf



Sun-SBOX1 */**
Sharps Box Bracket (Daniels - All Sizes)



Sun-SBOX2 */**
Sharps Box Bracket (Frontier) Upto 1 Litre



Sun-SBOX3 */**
Sharps Box Bracket (Frontier) 2 & 3 Litre



Sun-MPT/AGH1 */**
Alcohol Gel Holder



Sun-GBH1 */**
Glove Box Holder



Sun-MPT/Rail1 *
Medi-rail



Sun-MPT/BIN *
2.5 Litre Bin (Requires Medi-rail)



Sun-MPT55/LD */***
Tamper Resistant Locking Door



Sun-DIV1
Single Depth Tray Dividers (per set)




Sun-DIV4
Double Depth Tray Dividers (per set)



Sun-MPT55/CV
Trolley Cover

** Includes Medi-rail *** Locking door option includes a secure top and back panel


dimensions



110cm(H) 56cm(W) 52cm(D)

tray colours

please specify when ordering
(Tray labels and covers included)



clear blue green

Single Depth Trays
8cm(H) 47cm(W) 42cm(D)
Double Depth Trays
16cm(H) 47cm(W) 42cm(D)

Vista® 55 Trolleys

Suitable for everyday general healthcare use, Vista 55 provides standard working height with generous working space, smart design optimises storage with upto seven trays, and is functional, durable and easy to clean.



Code: **Sun-MPT55**
Vista 55 Trolley - 7 single depth trays



Code: **Sun-MPT55/LD**
Locking door for Vista 55 trolley



Code: **Sun-MPT56**
Vista 55 Trolley - 3 single/2 double depth trays



Code: **Sun-MPT55/LD**
Locking door for Vista 55 trolley



Code: **Sun-MPT57**
Vista 55 Trolley - 1 single/3 double depth trays



Code: **Sun-MPT55/LD**
Locking door for Vista 55 trolley

product information



Tilt & locking tray runner system



Trays integrate with all other Vista products



100mm castors with front brakes



Optional tamper resistant locking doors



Removable top shelf with integrated upstand



Easy clean design and materials



2 year warranty



Supplied fully assembled



British made

accessories * These accessories can only be fitted when ordered with the trolley



Sun-LC3A *
4 Antistatic Castors (front 2 braked)



Sun-BUF2 *
4 Protective Buffers



Sun-MPT/OXY1 **
Oxygen Bottle Holder (Size D)



Sun-MPT/LAP1 **
Adjustable Height Laptop Holder



Sun-MPT/CHART1 **
Adjustable Height Chart Holder



Sun-MPT/DRIP1 **
Adjustable Height Drip Pole



Sun-MPTS1 *
Folding Shelf



Sun-SBOX1 **
Sharps Box Bracket (Daniels - All Sizes)



Sun-SBOX2 **
Sharps Box Bracket (Frontier) Upto 1 Litre



Sun-SBOX3 **
Sharps Box Bracket (Frontier) 2 & 3 Litre



Sun-MPT/AGH1 **
Alcohol Gel Holder



Sun-GBH1 **
Glove Box Holder



Sun-MPT/Rail1 *
Medi-rail



Sun-MPT/BIN *
2.5 Litre Bin (Requires Medi-rail)



Sun-MPT60/LD ***
2 Tamper Resistant Locking Doors



Sun-DIV2
Double Depth Tray Dividers (per set)




Sun-DIV3
Single Depth Tray Dividers (per set)



Sun-MPT60/CV
Trolley Cover

** Includes Medi-rail *** Locking door option includes a secure top and back panel


dimensions



Single Depth Trays
8cm(H) 31cm(W) 42cm(D)
Double Depth Trays
16cm(H) 31cm(W) 42cm(D)

tray colours

please specify when ordering
(Tray labels and covers included)



clear blue green

Vista® 60 Trolleys

Designed to optimise storage where space is limited, Vista 60 provides additional working space and compact storage with upto 12 trays, and is functional, durable and easy to clean.



Code: **Sun-MPT60**
Vista 60 Trolley - 12 single depth trays

Code: **Sun-MPT60/LD**
2 Locking doors for Vista 60 trolley



Code: **Sun-MPT61**
Vista 60 Trolley - 6 single/3 double depth trays

Code: **Sun-MPT60/LD**
2 Locking doors for Vista 60 trolley



Code: **Sun-MPT62**
Vista 60 Trolley - 6 double depth trays

Code: **Sun-MPT60/LD**
2 Locking doors for Vista 60 trolley


product information



Tilt & locking tray runner system




Trays integrate with all other Vista products



100mm castors with front brakes




Optional tamper resistant locking doors



Removable top shelf with integrated upstand



Easy clean design and materials



2 year warranty



Supplied fully assembled



British made

accessories * These accessories can only be fitted when ordered with the trolley



Sun-LC3A *
4 Antistatic Castors (front 2 braked)



Sun-BUF2 *
4 Protective Buffers



Sun-MPT/OXY1 */**
Oxygen Bottle Holder (Size D)



Sun-MPT/LAP1 */**
Adjustable Height Laptop Holder



Sun-MPT/CHART1 */**
Adjustable Height Chart Holder



Sun-MPT/DRIP1 */**
Adjustable Height Drip Pole



Sun-MPTS1 *
Folding Shelf



Sun-SBOX1 */**
Sharps Box Bracket (Daniels - All Sizes)



Sun-SBOX2 */**
Sharps Box Bracket (Frontier) Upto 1 Litre



Sun-SBOX3 */**
Sharps Box Bracket (Frontier) 2 & 3 Litre



Sun-MPT/AGH1 */**
Alcohol Gel Holder



Sun-GBH1 */**
Glove Box Holder



Sun-MPT/Rail1 *
Medi-rail



Sun-MPT/BIN *
2.5 Litre Bin (Requires Medi-rail)



Sun-MPT4/LD */**
2 Tamper Resistant Locking Doors



Sun-DIV1
Single Depth Tray Dividers (per set)



Sun-DIV4
Double Depth Tray Dividers (per set)



Sun-MPT4/CV
Trolley Cover

** Includes Medi-rail *** Locking door option includes a secure top and back panel


dimensions



Single Depth Trays
8cm(H) 47cm(W) 42cm(D)
Double Depth Trays
16cm(H) 47cm(W) 42cm(D)

tray colours

please specify when ordering
(Tray labels and covers included)



Vista® 100 Trolleys

Designed for the most demanding healthcare environments, Vista 100 combines maximum working space with optimum tray storage, and is functional, durable and easy to clean.



Code: **Sun-MPT16**
Vista 100 Trolley - 12 single depth trays

Code: **Sun-MPT4/LD**
2 Locking doors for Vista 100 trolley



Code: **Sun-MPT4**
Vista 100 Trolley - 8 single/2 double depth trays

Code: **Sun-MPT4/LD**
2 Locking doors for Vista 100 trolley



Code: **Sun-MPT17**
Vista 100 Trolley - 6 single/3 double depth trays

Code: **Sun-MPT4/LD**
2 Locking doors for Vista 100 trolley

product information



Tilt & locking tray runner system



Trays integrate with all other Vista products



100mm castors with front brakes



Removable top shelf



Easy clean design and materials



2 year warranty



Supplied fully assembled



British made

accessories * These accessories can only be fitted when ordered with the trolley



Sun-LC3A *
4 Antistatic Castors
(front 2 braked)



Sun-BUF3 *
4 Protective Buffers



Sun-MPT/OXY1 */**
Oxygen Bottle Holder
(Size D)



Sun-MPTS1 *
Folding Shelf



Sun-SBOX1 */**
Sharps Box Bracket
(Daniels - All Sizes)



Sun-SBOX2 */**
Sharps Box Bracket
(Frontier) Upto 1 Litre



Sun-SBOX3 */**
Sharps Box Bracket
(Frontier) 2 & 3 Litre



Sun-MPT/AGH1 */**
Alcohol Gel Holder



Sun-GBH1 */**
Glove Box Holder



Sun-MPT/Rail1 *
Medi-rail



Sun-MPT/BIN *
2.5 Litre Bin
(Requires Medi-rail)



Sun-DIV1
Single Depth
Tray Dividers (per set)



Sun-DIV4
Double Depth
Tray Dividers (per set)

** Includes Medi-rail

dimensions

Vista 120



146cm(H)
54cm(W)
50cm(D)

Vista 240




146cm(H)
102cm(W)
50cm(D)

Single Depth Trays
8cm(H) 47cm(W) 42cm(D)
Double Depth Trays
16cm(H) 47cm(W) 42cm(D)

tray colours

please specify when ordering
(Tray labels and covers included)



clear

blue

green

Vista® 120 & 240 Trolleys

Suitable for healthcare environments with demanding mobile storage requirements and limited space, Vista 120 & 240 combine high level storage, with multiple tray configurations and are functional, durable and easy to clean.



Code: **Sun-MPT120**
Vista 120 Trolley - 12 single depth trays



Code: **Sun-MPT240**
Vista 240 Trolley - 24 single depth trays



Code: **Sun-MPT121**
Vista 120 Trolley - 6 single/3 double depth trays



Code: **Sun-MPT241**
Vista 240 Trolley - 12 single/6 double depth trays



Code: **Sun-MPT122**
Vista 120 Trolley - 6 double depth trays



Code: **Sun-MPT242**
Vista 240 Trolley - 12 double depth trays

product information



Tilt & locking tray runner system



Trays integrate with all other Vista products



Removable mirror polished recessed stainless steel tray (304 grade)



Narrow width (39cm) compact design



100mm castors with front brakes



4 Protective buffers as standard



Optional tamper resistant locking door



Easy clean design and materials



2 year warranty



Supplied fully assembled



British made

accessories * These accessories can only be fitted when ordered with the trolley



Sun-LC3A *
4 Antistatic Castors (front 2 braked)



Sun-MPTS1 *
Folding Shelf



Sun-SBOX1 */**
Sharps Box Bracket (Daniels - All Sizes)



Sun-SBOX2 */**
Sharps Box Bracket (Frontier) Upto 1 Litre



Sun-SBOX3 */**
Sharps Box Bracket (Frontier) 2 & 3 Litre



Sun-MPT/AGH1 */**
Alcohol Gel Holder



Sun-GBH1 */**
Glove Box Holder



Sun-MPT/Rail1 *
Medi-rail



Sun-MPT/BIN *
2.5 Litre Bin (Requires Medi-rail)



Sun-CVT30/LD */***
Tamper Resistant Locking Door



Sun-DIV2
Double Depth Tray Dividers (per set)




Sun-DIV3
Single Depth Tray Dividers (per set)



Sun-MPT5/CV
Trolley Cover

** Includes Medi-rail *** Locking door option includes a secure top and back panel

dimensions



101cm(H) 39cm(W) 52cm(D)
47cm(W) 60cm(D) with buffers

Single Depth Trays
8cm(H) 31cm(W) 42cm(D)
Double Depth Trays
16cm(H) 31cm(W) 42cm(D)

tray colours

please specify when ordering
(Tray labels and covers included)



clear blue green

Vista® 30 Clinical Trolleys

Designed for environments with limited space, Vista 30 Clinical provides standard working height with compact design, optimum storage, a removable recessed stainless steel tray and is functional, durable and easy to clean.



Code: **Sun-CVT30**
Vista 30 Clinical trolley - 6 single depth trays



Code: **Sun-CVT30/LD**
Locking door for Vista 30 clinical trolley



Code: **Sun-CVT31**
Vista 30 Clinical trolley -
2 single/2 double depth trays



Code: **Sun-CVT30/LD**
Locking door for Vista 30 clinical trolley



Code: **Sun-CVT32**
Vista 30 Clinical trolley - 3 double depth trays



Code: **Sun-CVT30/LD**
Locking door for Vista 30 clinical trolley

product information



Tilt & locking tray runner system



Trays integrate with all other Vista products



Removable mirror polished recessed stainless steel tray (304 grade)




Low level (80cm), ideal for under worktops



100mm castors with front brakes



Optional tamper resistant locking door



Stainless steel handle



4 Protective buffers as standard



Easy clean design and materials



2 year warranty



Supplied fully assembled



British made

accessories * These accessories can only be fitted when ordered with the trolley



Sun-LC3A *
4 Antistatic Castors
(front 2 braked)



Sun-MPT1 *
Folding Shelf



Sun-SBOX1 */**
Sharps Box Bracket
(Daniels - All Sizes)



Sun-SBOX2 */**
Sharps Box Bracket
(Frontier) Upto 1 Litre



Sun-SBOX3 */**
Sharps Box Bracket
(Frontier) 2 & 3 Litre



Sun-MPT/AGH1 */**
Alcohol Gel Holder



Sun-GBH1 */**
Glove Box Holder



Sun-MPT/Rail1 *
Medi-rail



Sun-MPT/BIN *
2.5 Litre Bin
(Requires Medi-rail)



Sun-CVT40/LD */***
Tamper Resistant
Locking Door



Sun-DIV1
Single Depth
Tray Dividers (per set)



Sun-DIV4
Double Depth
Tray Dividers (per set)



Sun-MPT40/CV
Trolley Cover

** Includes Medi-rail *** Locking door option includes a secure top and back panel

dimensions



80cm(H) 56cm(W) 52cm(D)
64cm(W) 60cm(D) with buffers

Single Depth Trays
8cm(H) 47cm(W) 42cm(D)
Double Depth Trays
16cm(H) 47cm(W) 42cm(D)

tray colours

please specify when ordering
(Tray labels and covers included)



clear blue green

Vista® 40 Clinical Trolleys

Designed to fit under a worktop, Vista 40 Clinical can be stored neatly away, while providing optimum number of trays for storage, removable recessed stainless steel tray and is functional, durable and easy to clean.



Code: **Sun-CVT40**
Vista 40 Clinical trolley
5 single depth trays



Code: **Sun-CVT40/LD**
Locking door for Vista 40 clinical trolley



Code: **Sun-CVT41**
Vista 40 Clinical trolley
3 single/1 double depth tray



Code: **Sun-CVT40/LD**
Locking door for Vista 40 clinical trolley



Code: **Sun-CVT42**
Vista 40 Clinical trolley
1 single/2 double depth trays



Code: **Sun-CVT40/LD**
Locking door for Vista 40 clinical trolley

product information



Tilt & locking tray runner system



Trays integrate with all other Vista products



Removable mirror polished recessed stainless steel tray (304 grade)



100mm castors with front brakes



Optional tamper resistant locking door



4 Protective buffers as standard



Easy clean design and materials



2 year warranty



Supplied fully assembled



British made

accessories * These accessories can only be fitted when ordered with the trolley



Sun-LC3A *
4 Antistatic Castors (front 2 braked)



Sun-MPT/OXY1 **
Oxygen Bottle Holder (Size D)



Sun-MPT/LAP1 **
Adjustable Height Laptop Holder



Sun-MPT/CHART1 **
Adjustable Height Chart Holder



Sun-MPT/DRIP1 **
Adjustable Height Drip Pole



Sun-MPTS1 *
Folding Shelf



Sun-SBOX1 **
Sharps Box Bracket (Daniels - All Sizes)



Sun-SBOX2 **
Sharps Box Bracket (Frontier) Upto 1 Litre



Sun-SBOX3 **
Sharps Box Bracket (Frontier) 2 & 3 Litre



Sun-MPT/AGH1 **
Alcohol Gel Holder



Sun-GBH1 **
Glove Box Holder



Sun-MPT/Rail1 *
Medi-rail



Sun-MPT/BIN *
2.5 Litre Bin (Requires Medi-rail)



Sun-CVT50/LD ***
Tamper Resistant Locking Door



Sun-DIV1
Single Depth Tray Dividers (per set)




Sun-DIV4
Double Depth Tray Dividers (per set)



Sun-MPT1/CV
Trolley Cover

** Includes Medi-rail *** Locking door option includes a secure top and back panel


dimensions



101cm(H) 56cm(W) 52cm(D)
64cm(W) 60cm(D) with buffers

tray colours

please specify when ordering
(Tray labels and covers included)



clear blue green

Single Depth Trays
8cm(H) 47cm(W) 42cm(D)
Double Depth Trays
16cm(H) 47cm(W) 42cm(D)

Vista® 50 Clinical Trolleys

Suitable for everyday general healthcare use, Vista 50 Clinical provides standard working height with generous working space, optimum tray storage, removable recessed stainless steel tray and is functional, durable and easy to clean.



Code: **Sun-CVT50**
Vista 50 Clinical trolley - 6 single depth trays



Code: **Sun-CVT50/LD**
Locking door for Vista 50 clinical trolley



Code: **Sun-CVT51**
Vista 50 Clinical trolley
4 single/1 double depth tray



Code: **Sun-CVT50/LD**
Locking door for Vista 50 clinical trolley



Code: **Sun-CVT52**
Vista 50 Clinical Trolley - 3 double depth trays



Code: **Sun-CVT50/LD**
Locking door for Vista 50 clinical trolley

product information



Tilt & locking tray runner system



Trays integrate with all other Vista products



Removable mirror polished recessed stainless steel tray (304 grade)



100mm castors with front brakes



Optional tamper resistant locking doors



4 Protective buffers as standard



Easy clean design and materials



2 year warranty



Supplied fully assembled



British made

accessories * These accessories can only be fitted when ordered with the trolley



Sun-LC3A *
4 Antistatic Castors (front 2 braked)



Sun-MPT/OXY1 **
Oxygen Bottle Holder (Size D)



Sun-MPT/LAP1 **
Adjustable Height Laptop Holder



Sun-MPT/CHART1 **
Adjustable Height Chart Holder



Sun-MPT/DRIP1 **
Adjustable Height Drip Pole



Sun-MPTS1 *
Folding Shelf



Sun-SBOX1 **
Sharps Box Bracket (Daniels - All Sizes)



Sun-SBOX2 **
Sharps Box Bracket (Frontier) Upto 1 Litre



Sun-SBOX3 **
Sharps Box Bracket (Frontier) 2 & 3 Litre



Sun-MPT/AGH1 **
Alcohol Gel Holder



Sun-GBH1 **
Glove Box Holder



Sun-MPT/Rail1 *
Medi-rail



Sun-MPT/BIN *
2.5 Litre Bin (Requires Medi-rail)



Sun-CVT60/LD ***
2 Tamper Resistant Locking Doors



Sun-DIV2
Double Depth Tray Dividers (per set)




Sun-DIV3
Single Depth Tray Dividers (per set)



Sun-MPT60/CV
Trolley Cover

** Includes Medi-rail *** Locking door option includes a secure top and back panel

dimensions




101cm(H) 70cm(W) 52cm(D)
78cm(W) 60cm(D) with buffers

Single Depth Trays
8cm(H) 31cm(W) 42cm(D)
Double Depth Trays
16cm(H) 31cm(W) 42cm(D)

tray colours

please specify when ordering
(Tray labels and covers included)



clear blue green

Vista® 60 Clinical Trolleys

Designed to optimise storage where space is limited, Vista 60 Clinical provides additional working space with compact storage using upto 12 trays, removable recessed stainless steel tray and is functional, durable and easy to clean.



Code: **Sun-CVT60**
Vista 60 Clinical trolley - 12 single depth trays



Code: **Sun-CVT60/LD**
2 Locking doors for Vista 60 clinical trolley



Code: **Sun-CVT61**
Vista 60 Clinical trolley
6 single/3 double depth trays



Code: **Sun-CVT60/LD**
2 Locking doors for Vista 60 clinical trolley




Code: **Sun-CVT62**
Vista 60 Clinical trolley - 6 double depth trays



Code: **Sun-CVT60/LD**
2 Locking doors for Vista 60 clinical trolley

product information



Choice of 5 coloured tray and frame colours



Tilt & locking tray runner system



Trays integrate with all other Vista products



Tamper resistant locking door



Top, bottom and back panels provide added security



Translucent tray lids, tray labels and covers included



Removable top shelf with integrated upstand



Clean & Dirty labels included



100mm castors with front brakes



4 Protective buffers as standard



Front mounted stainless steel handle (Vista 40 only)



Easy clean design and materials



2 year warranty



Supplied fully assembled



British made

accessories * These accessories can only be fitted when ordered with the trolley



Sun-LC3A *
4 Antistatic Castors (front 2 braked)



Sun-MPT/OXY1 **
Oxygen Bottle Holder Size D (Vista 50 Only)



Sun-MPT/LAP1 **
Adjustable Height Laptop Holder (Vista 50 Only)



Sun-MPT/CHART1 **
Adjustable Height Chart Holder (Vista 50 Only)



Sun-MPT/DRIP1 **
Adjustable Height Drip Pole (Vista 50 Only)



Sun-MPTS1 *
Folding Shelf



Sun-SBOX1 **
Sharps Box Holder (Daniels - All Sizes)



Sun-SBOX2 **
Sharps Box Holder (Frontier) Upto 1 Litre



Sun-SBOX3 **
Sharps Box Holder (Frontier) 2 & 3 Litre



Sun-MPT/AGH1 **
Alcohol Gel Holder



Sun-GBH1 **
Glove Box Holder



Sun-MPT/Rail1 *
Medi-rail



Sun-MPT/BIN *
2.5 Litre Bin (Requires Medi-rail)



Sun-DIV1
Single Depth Tray Dividers (per set)



Trolley Covers
Sun-MPT40/CV (Vista 40 Clinical)
Sun-MPT1/CV (Vista 50 Clinical)

** Includes Medi-rail

dimensions



80cm(H) 56cm(W) 52cm(D)
64cm(W) 60cm(D) with buffers


Vista 50 Clinical Trolley



101cm(H) 56cm(W) 52cm(D)
64cm(W) 60cm(D) with buffers

frame & tray colours

please specify when ordering



blue green yellow white/clear red

Designed for clinical & disinfecting environments, Vista Clinical are colour designated for clean and dirty areas. Vista 50 provides standard working height with maximum tray storage, while Vista 40 fits neatly under a worktop, both trolleys are fully secure, functional, durable and easy to clean.



Code: **Sun-CMPT40**
Vista 40 Clinical Trolley - 5 single depth trays with optional sharps box holder and folding shelf



Code: **Sun-CMPT1**
Vista 50 Clinical Trolley - 6 single depth trays with optional sharps box holder, glove box holder and folding shelf

product information

J.I.T. design

Designed for "just in time" working

Laptop or vital pack area

Medical notes storage in trays

Tilt & locking tray runner system

Top locking drawer for secure storage

Combination lock (no keys required)

4 Protective buffers as standard

100mm castors with front brakes

Removable top shelf with integrated upstand

Easy clean design and materials

2 year warranty

Supplied fully assembled

British made

accessories * These accessories can only be fitted when ordered with the trolley

Sun-LC3A *
4 Antistatic Castors (front 2 braked)

Sun-MPT/OXY1 */**
Oxygen Bottle Holder (Size D)

Sun-MPTS1 *
Folding Shelf

Sun-SBOX1 */**
Sharps Box Bracket (Daniels - All Sizes)

Sun-SBOX2 */**
Sharps Box Bracket (Frontier) Upto 1 Litre

Sun-SBOX3 */**
Sharps Box Bracket (Frontier) 2 & 3 Litre

Sun-MPT/AGH1 */**
Alcohol Gel Holder

Sun-GBH1 */**
Glove Box Holder

Sun-MPT/Rail1 *
Medi-rail

Sun-MPT/BIN *
2.5 Litre Bin (Requires Medi-rail)

Sun-MPT/LABEL1
Wipe Clean Labels 2 Per Tray (5 Sets Per Pack)


Sun-DIV1
Single Depth Tray Dividers (per set)

Sun-DIV4
Double Depth Tray Dividers (per set)

Sun-MPT90/LD */***
Tamper Resistant Locking Door

** Includes Medi-rail *** Locking door option includes a secure top and back panel


dimensions



Single Depth Trays
8cm(H) 47cm(W) 42cm(D)
Double Depth Trays
16cm(H) 47cm(W) 42cm(D)

tray colours

please specify when ordering



Designed for "just in time" working during ward rounds, Vista 90 Action Station provides working space for a laptop or "Vital Pack" with a lockable top tray, additional tray storage for patient notes and equipment and is functional, durable and easy to clean.



Single column trolley specifically designed in conjunction with leading Consultant Dr Gordon Caldwell FRCP (London), Consultant Physician, of Worthing Hospital, Diabetes Centre for "just in time" working during ward rounds.

Dr Caldwell comments, "I have found that this trolley significantly improves ward round efficiency, it holds the team of doctors around the trolley and gives them a large surface to write on, the horizontal filing system is also far easier to use than traditional vertical files, we more clearly define team members' roles, therefore more productive doctors".



Code: **Sun-MPT90**
Vista 90 Action Station Trolley
9 single depth trays



Code: **Sun-MPT91**
Vista 90 Action Station Trolley
5 single/ 2 double depth trays

product information

J.I.T. design

Designed for "just in time" working

Laptop or vital pack area

Medical notes storage in trays

Tilt & locking tray runner system

Top locking drawer for secure storage

Combination lock (no keys required)

4 Protective buffers as standard

100mm castors with front brakes

Removable top shelf with integrated upstand

Easy clean design and materials

2 year warranty

Supplied fully assembled

British made

accessories * These accessories can only be fitted when ordered with the trolley

Sun-LC3A *
4 Antistatic Castors (front 2 braked)

Sun-MPT/OXY1 */**
Oxygen Bottle Holder (Size D)

Sun-MPTS1 *
Folding Shelf

Sun-SBOX1 */**
Sharps Box Bracket (Daniels - All Sizes)

Sun-SBOX2 */**
Sharps Box Bracket (Frontier) Upto 1 Litre

Sun-SBOX3 */**
Sharps Box Bracket (Frontier) 2 & 3 Litre

Sun-MPT/AGH1 */**
Alcohol Gel Holder

Sun-GBH1 */**
Glove Box Holder

Sun-MPT/Rail1 *
Medi-rail

Sun-MPT/BIN *
2.5 Litre Bin (Requires Medi-rail)

Sun-MPT/LABEL1
Wipe Clean Labels 2 Per Tray (5 Sets Per Pack)

Sun-DIV1
Single depth Tray Dividers (per set)

Sun-DIV4
Double Depth Tray Dividers (per set)

Sun-MPT180/LD */***
2 Tamper Resistant Locking Doors

** Includes Medi-rail *** Locking door option includes a secure top and back panel

dimensions

122cm(H) 102cm(W) 50cm(D)
110cm(W) 58cm(D) with buffers

Single Depth Trays
8cm(H) 47cm(W) 42cm(D)
Double Depth Trays
16cm(H) 47cm(W) 42cm(D)

tray colours

please specify when ordering

clear blue green

Vista® 180 Action Station Trolley

Designed for "just in time" working during ward rounds, Vista 180 Action Station provides working space for a laptop or "Vital Pack" with a lockable top tray, additional tray storage for patient notes and equipment and is functional, durable and easy to clean.



Double column trolley specifically designed in conjunction with leading Consultant Dr Gordon Caldwell FRCP (London), Consultant Physician, of Worthing Hospital, Diabetes Centre for "just in time" working during ward rounds.

Dr Caldwell comments, "I have found that this trolley significantly improves ward round efficiency, it holds the team of doctors around the trolley and gives them a large surface to write on, the horizontal filing system is also far easier to use than traditional vertical files, we more clearly define team members' roles, therefore more productive doctors".



Code: **Sun-MPT180**
Vista 180 Action Station Trolley
18 single depth trays



Code: **Sun-MPT181**
Vista 180 Action Station Trolley
10 single/4 double depth trays

product information



Tilt & locking tray runner system



Upto eight trays



Trays integrate with all other Vista products



Single or double depth tray combinations



Adjustable feet for levelling



Optional tamper resistant locking door



Removable top shelf with integrated upstand



Modules can be attached together to form larger systems



Optional castor upgrades available



Durable epoxy coated white frame



Easy clean design and materials



2 year warranty



Supplied fully assembled



British made

accessories * These accessories can only be fitted when ordered with the storage module



Sun-SMLC2 *
4 Castors
(front 2 braked)



Sun-SMLC3A *
4 Antistatic Castors
(front 2 braked)



Sun-BUF3 *
4 Protective Buffers



Sun-SMS1/LD */***
Tamper Resistant
Locking Door




Sun-DIV1
Single Depth
Tray Dividers (per set)



Sun-DIV4
Double Depth
Tray Dividers (per set)


*** Locking door option includes a secure top and back panel

dimensions



Single Depth Trays
8cm(H) 47cm(W) 42cm(D)
Double Depth Trays
16cm(H) 47cm(W) 42cm(D)

tray colours
please specify when ordering
(Tray labels and covers included)



Vista® Low Level Storage Modules

Modular, flexible and colourful, Vista Low Level Storage Modules provide a general storage solution for healthcare environments, using premium materials which are functional, durable and easy to clean.



Code: **Sun-SM1/8**
Vista low level module
8 single depth trays

Code: **Sun-SM1/4/2D**
Vista low level module
4 single/2 double depth trays



Code: **Sun-SM1/4D**
Vista low level module
4 double depth trays

Code: **Sun-SM1/LD**
Locking door for Vista
low level module



Modules can be attached together to form larger systems

product information



Tilt & locking tray runner system



Upto sixteen trays



Trays integrate with all other Vista products



Single or double depth tray combinations



Adjustable feet for levelling



Optional tamper resistant locking doors



Removable top shelf with integrated upstand



Modules can be attached together to form larger systems



Optional castor upgrades available



Durable epoxy coated white frame



Easy clean design and materials



2 year warranty



Supplied fully assembled



British made

accessories * These accessories can only be fitted when ordered with the storage module



Sun-SMLC2 *
4 Castors
(front 2 braked)



Sun-SMLC3A *
4 Antistatic Castors
(front 2 braked)



Sun-BUF3 *
4 Protective Buffers



Sun-SMS2/LD * /***
2 Tamper Resistant Locking Doors




Sun-DIV1
Single Depth Tray Dividers (per set)



Sun-DIV4
Double Depth Tray Dividers (per set)


*** Locking door option includes a secure top and back panel

dimensions



Single Depth Trays
8cm(H) 47cm(W) 42cm(D)
Double Depth Trays
16cm(H) 47cm(W) 42cm(D)

tray colours
please specify when ordering
(Tray labels and covers included)



Vista® High Level Storage Modules

Modular, flexible and colourful, Vista High Level Storage Modules provide optimum storage for general healthcare environments, using premium materials which are functional, durable and easy to clean.



Code: **Sun-SM2/16**
Vista high level module
16 single depth trays

Code: **Sun-SM2/8/4D**
Vista high level module
8 single/4 double depth trays


Code: **Sun-SM2/8D**
Vista high level module
8 double depth trays




Code: **Sun-SM2/LD**
2 locking doors for Vista
high level modules

Modules can be attached
together to form larger systems


product information




Tilt & locking tray runner system




Clear, durable, easy to clean ABS plastic trays




Contemporary door and profile with 'high gloss' finish




Trays integrate with all other Vista products




Strong durable 18mm white carcass construction




Locking doors (2 keys supplied)




Durable 38mm White postformed worktop (VBU4W & VBU1W)




One piece removable plinth




Tray labels and covers included




Easy clean design and materials



2 year warranty



Supplied fully assembled



British made

accessories



Sun-VBU/KEY
Extra Keys (each)




Sun-DIV1
Single Depth
Tray Dividers (per set)



Sun-DIV4
Double Depth
Tray Dividers (per set)


dimensions

Sun-VBU1W




91cm(H)
51.5cm(W)
60cm(D)

Sun-VBU4W



91cm(H)
103cm(W)
60cm(D)

Sun-VLUT1W



161cm(H)
51.5cm(W)
60cm(D)

Single Depth Trays
8cm(H) 47cm(W) 42cm(D)

Double Depth Trays
16cm(H) 47cm(W) 42cm(D)

HEAVY ITEM

Designed to be aesthetically pleasing, accessible and secure, Vista base & larder cabinets provide maximum tray storage for equipment and disposable, using premium materials which are functional, durable and easy to clean.



Code: **Sun-VBU4W**
103cm wide Vista base unit with 38mm white postformed worktop
10 single/2 double depth trays



Code: **Sun-VBU1W**
51.5cm wide Vista base unit with 38mm white postformed worktop
5 single/1 double depth trays







Code: **Sun-VLUT1W**
51.5cm wide Vista larder unit with 18mm white MFC top
7 single/3 double depth trays



hygienic
safe
procedures

Dressing & Surgical Trolleys, Theatre Equipment

Designed for clinical procedures in healthcare environments, our range of trolleys, support surfaces and ancillary equipment can be sterilised, whilst being functional, durable and easy to clean.

	PAGE
 Dressing Trolleys	93-96
 Surgical Trolleys	97-100
 Surgical Trolleys (Mirror Polished)	101-104
 Theatre Equipment	105-106

product information



Unique glass effect safety trays



2 or 3 trays



Lift-off top tray with integral handles



Hygienically designed raised upstands



Trays 'locate' around frame to prevent movement



Tray protection buffer ring



Durable non-marking 75mm castors



Durable epoxy coated white frame



Optional drawers available in a choice of 4 finishes



4 Drawer finishes



Maximum individual tray load



Easy clean design and materials



2 year warranty



Supplied fully assembled



British made

accessories * These accessories can only be fitted when ordered with the trolley



Sun-DT4D*
Single Drawer Unit
(small trolleys)



Sun-DTDT1
Drawer insert for small
trolley drawer (each)



Sun-DT6D*
Single Drawer Unit
(medium trolleys)



Sun-DTDT2
Drawer insert for medium
trolley drawer (each)



Sun-DT8D*
Double Drawer Unit
(large trolleys)



Sun-DTDT3
Drawer insert for large
trolley drawer (each)



Sun-BUF1*
4 Protective buffers




Sun-LC1*
2 front
braked castors



Sun-LC1A*
4 Antistatic castors
(front 2 braked)


dimensions

small




86cm(H) 46cm(W) 52cm(D)

medium



86cm(H) 66cm(W) 52cm(D)

large



86cm(H) 86cm(W) 52cm(D)

drawer finishes - please specify when ordering



White (high gloss)



Titanium (high gloss)



Beech effect



Maple effect

Dressing Trolleys - Glass Effect Safety Trays

Innovative in design, and in use of materials, our range of dressing trolleys provide hygienic mobile support and storage, whilst being functional, durable and easy to clean.

small - 46cm wide dressing trolleys



Code: **Sun-DT4G2**
2 glass effect safety trays



Code: **Sun-DT4G3**
3 glass effect safety trays



Sun-DT4G3 with optional
Sun-DT4D single drawer unit

medium - 66cm wide dressing trolleys



Code: **Sun-DT6G2**
2 glass effect safety trays



Code: **Sun-DT6G3**
3 glass effect safety trays



Sun-DT6G3 with optional
Sun-DT6D single drawer unit

large - 86cm wide dressing trolleys



Code: **Sun-DT8G2**
2 glass effect safety trays



Code: **Sun-DT8G3**
3 glass effect safety trays



Sun-DT8G3 with optional
Sun-DT8D double drawer unit



product information



Stainless steel trays
(304 grade/DP2 finish)



2 or 3
trays



Lift-off top tray with
integral handles



Hygienically designed
raised upstands



Trays 'locate' around
frame to prevent
movement



Tray protection
buffer ring



Durable non-marking
75mm castors



Durable epoxy coated white frame



Optional drawers
available in a choice
of 4 finishes



4 Drawer finishes



10Kg



Easy clean design
and materials



2 year warranty




Supplied fully
assembled



British made

accessories * These accessories can only be fitted when ordered with the trolley



Sun-DT4D*
Single Drawer Unit
(small trolleys)



Sun-DTD1
Drawer insert for small
trolley drawer (each)




Sun-DT6D*
Single Drawer Unit
(medium trolleys)



Sun-DTD2
Drawer insert for medium
trolley drawer (each)




Sun-DT8D*
Double Drawer Unit
(large trolleys)



Sun-DTD3
Drawer insert for large
trolley drawer (each)



Sun-BUF1*
4 Protective buffers




Sun-LC1*
2 front
braked castors



Sun-LC1A*
4 Antistatic castors
(front 2 braked)


dimensions

small




86cm(H) 46cm(W) 52cm(D)

medium



86cm(H) 66cm(W) 52cm(D)

large



86cm(H) 86cm(W) 52cm(D)

drawer finishes - please specify when ordering



White (high gloss)



Titanium (high gloss)



Beech effect



Maple effect

Dressing Trolleys - Stainless Steel Trays

Innovative in design, and in use of materials, our range of dressing trolleys provide hygienic mobile support and storage, whilst being functional, durable and easy to clean.

small - 46cm wide dressing trolleys



Code: Sun-DT4S2
2 stainless steel trays



Code: Sun-DT4S3
3 stainless steel trays



Sun-DT4S3 with optional
Sun-DT4D single drawer unit

medium - 66cm wide dressing trolleys



Code: Sun-DT6S2
2 stainless steel trays



Code: Sun-DT6S3
3 stainless steel trays



Sun-DT6S3 with optional
Sun-DT6D single drawer unit

large - 86cm wide dressing trolleys



Code: Sun-DT8S2
2 stainless steel trays



Code: Sun-DT8S3
3 stainless steel trays



Sun-DT8S3 with optional
Sun-DT8D double drawer unit



product information



Stainless steel trays
(304 grade/DP2 finish)



2 or 3 stainless steel
removable trays



Lift-off top tray with
integral handles



Hygienically designed
raised upstands



Trays 'locate' around
frame to prevent
movement



Tray protection
buffer ring



Durable non-marking
75mm castors



Polished stainless
steel frame
(304 grade)



Maximum individual
tray load



Easy clean design
and materials



2 year warranty



Supplied fully
assembled



British made

accessories

* These accessories can only be fitted when ordered with the trolley



Sun-BUF1*
4 Protective buffers




Sun-LC1*
2 front
braked castors




Sun-LC1A*
4 Antistatic castors
(front 2 braked)


dimensions



86cm(H) 46cm(W) 52cm(D)



86cm(H) 66cm(W) 52cm(D)



86cm(H) 86cm(W) 52cm(D)

Surgical Trolleys - Stainless Steel Trays

Innovative in design, and in use of materials, our range of surgical trolleys provide hygienic mobile support for clinical environments, whilst being functional, durable and easy to clean.

small - 46cm wide surgical trolleys



Code: **Sun-ST4S2**
2 stainless steel trays



Code: **Sun-ST4S3**
3 stainless steel trays

medium - 66cm wide surgical trolleys



Code: **Sun-ST6S2**
2 stainless steel trays



Code: **Sun-ST6S3**
3 stainless steel trays

large - 86cm wide surgical trolleys



Code: **Sun-ST8S2**
2 stainless steel trays



Code: **Sun-ST8S3**
3 stainless steel trays

product information



Removable stainless steel shelf (304 grade/DP2 finish))



Lift-off stainless steel trays (304 grade/DP2 finish)



2 or 3 stainless steel shelf/tray combinations



Shelves & trays 'locate' around frame to prevent movement



Shelf/tray protection buffer ring



Durable non-marking 75mm castors



Polished stainless steel frame (304 grade)



Maximum individual tray load



Easy clean design and materials



Supplied fully assembled



Supplied fully assembled



British made

accessories

* These accessories can only be fitted when ordered with the trolley



Sun-BUF1*
4 Protective buffers




Sun-LC1*
2 front braked castors




Sun-LC1A*
4 Antistatic castors (front 2 braked)


dimensions



86cm(H) 46cm(W) 52cm(D)



86cm(H) 66cm(W) 52cm(D)



86cm(H) 86cm(W) 52cm(D)

Surgical Trolleys - Stainless Steel Flat Top Shelf & Trays

Innovative in design, and in use of materials, our range of surgical trolleys provide hygienic mobile support for clinical environments, whilst being functional, durable and easy to clean.

small - 46cm wide surgical trolleys



Code: **Sun-ST4S2SF**
1 stainless steel shelf and 1 tray



Code: **Sun-ST4S3SF**
1 stainless steel shelf and 2 trays

medium - 66cm wide surgical trolleys



Code: **Sun-ST6S2SF**
1 stainless steel shelf and 1 tray



Code: **Sun-ST6S3SF**
1 stainless steel shelf and 2 trays

large - 86cm wide surgical trolleys



Code: **Sun-ST8S2SF**
1 stainless steel shelf and 1 tray



Code: **Sun-ST8S3SF**
1 stainless steel shelf and 2 trays

product information

**premium
quality**

Stainless steel shelves, trays and frame (304 grade)

Mirror polished stainless steel shelves, trays and frame

Fully welded seamless frame

BS 4068 : 1977 Compliant *

**3 shelf
options**

Flat, fixed or removable/reversible shelves

Mirror polished lift off flat shelves (Sun-STFW/RFS2 only)

Removable, reversible folded shelves (Sun-STFW/RRFS2 only)

Fully welded fixed shelves (Sun-STFW/FFS2 only)

Anti-static shelf/tray protection buffers

75mm antistatic castors

4 Protective buffers

10Kg

Maximum individual tray load

eze

Easy clean design and materials

2 Year Warranty

Supplied fully assembled

* Designed to comply with the requirements of the hospital market (BS4068: 1977)

accessories * These accessories can only be fitted when ordered with the trolley

Sun-LC1A*
4 Antistatic castors (front 2 braked)

dimensions

small

86cm(H) 45cm(W) 45cm(D)
53cm(W) 53cm(D) with buffers

medium

86cm(H) 61cm(W) 45cm(D)
69cm(W) 53cm(D) with buffers

Premium use of materials, and finished to the highest standard, our range of 'mirror polished' surgical trolleys provide hygienic mobile support for clinical environments and theatres, whilst being functional, durable and easy to clean.

small - 45cm wide
surgical trolleys (mirror polished)



Code: **Sun-STFW4/RFS2**
2 removable flat shelves



Code: **Sun-STFW4/RRFS2**
2 removable reversible folded shelves (flange up or down)



Code: **Sun-STFW4/FFS2**
2 fully welded fixed shelves (flange down)

medium - 61cm wide
surgical trolleys (mirror polished)



Code: **Sun-STFW6/RFS2**
2 removable flat shelves



Code: **Sun-STFW6/RRFS2**
2 removable reversible folded shelves (flange up or down)



Code: **Sun-STFW6/FFS2**
2 fully welded fixed shelves (flange down)

product information

**premium
quality**

Stainless steel shelves, trays and frame (304 grade)

Mirror polished stainless steel shelves, trays and frame

Fully welded seamless frame

BS 4068 : 1977 Compliant *

**3 shelf
options**

Flat, fixed or removable/reversible shelves

Mirror polished lift off flat shelves (Sun-STFW/RFS2 only)

Removable, reversible folded shelves (Sun-STFW/RRFS2 only)

Fully welded fixed shelves (Sun-STFW/FFS2 only)

Anti-static shelf/tray protection buffers

75mm antistatic castors

4 Protective buffers

10Kg

Maximum individual tray load

eze

Easy clean design and materials

2 Year Warranty

Supplied fully assembled

* Designed to comply with the requirements of the hospital market (BS4068: 1977)

accessories * These accessories can only be fitted when ordered with the trolley

Sun-LC1A*
4 Antistatic castors (front 2 braked)

dimension

large

86cm(H) 75cm(W) 45cm(D)
83cm(W) 53cm(D) with buffers

extra large

86cm(H) 91cm(W) 45cm(D)
99cm(W) 53cm(D) with buffers

Premium use of materials, and finished to the highest standard, our range of 'mirror polished' surgical trolleys provide hygienic mobile support for clinical environments and theatres, whilst being functional, durable and easy to clean.

large - 75cm wide
surgical trolleys (mirror polished)



Code: **Sun-STFW7/RFS2**
2 removable flat shelves



Code: **Sun-STFW7/RRFS2**
2 removable reversible folded shelves (flange up or down)



Code: **Sun-STFW7/FFS2**
2 fully welded fixed shelves (flange down)

extra large - 91cm wide
surgical trolleys (mirror polished)



Code: **Sun-STFW9/RFS2**
2 removable flat shelves



Code: **Sun-STFW9/RRFS2**
2 removable reversible folded shelves (flange up or down)



Code: **Sun-STFW9/FFS2**
2 fully welded fixed shelves (flange down)

product information



Premium quality
(304 grade)
stainless steel



Fully welded
seamless frame



Polished (304 grade)
stainless steel bowls,
trays and buckets



HEIGHT RANGE
72-108cm
Adjustable height
bowl stands
(Sun-TBSV01/02 only)



HEIGHT RANGE
74-112cm
Adjustable height
mayo tables
(Sun-TMAY02 only)



Five castor base
for stability
and manoeuvrability



Removable bowls



Removable spiral
bowl and bucket
buffers



Removable
bucket with handle
(18 litre capacity)



Removable
stainless steel
instrument tray



Height adjustment
locking handle



50mm antistatic
castors



eze
Easy clean design
and materials




2 year warranty
2 Year Warranty




Supplied Assembled
Supplied fully
assembled


dimensions




77cm - 108cm(H)
58cm(W) 58cm(D)
Bowl 31cm Ø - 9cm(H)




77cm - 108cm(H)
74cm(W) 58cm(D)
Bowls 31cm Ø - 9cm(H)




85cm(H) 45cm(W) 45cm(D)
Bowls 31cm Ø - 9cm(H)




85cm(H) 75cm(W) 39cm(D)
Bowls 31cm Ø - 9cm(H)



74cm - 112cm(H)
67cm(W) 47cm(D)
Tray 2cm(H) 47cm(W)
36cm(D)



52cm(H) 38cm(W) 38cm(D)
Bucket 29cm Ø - 28cm(H)
18 Litre capacity



19cm(H) 31cm(W) 31cm(D)
Bowl 31cm Ø - 9cm(H)

Premium use of materials, and finished to the highest standard, our range of polished stainless steel theatre equipment provide hygienic mobile support and solutions for clinical environments and theatres, whilst being functional, durable and easy to clean.



Code: **Sun-TBSV01**
Single variable height
bowl stand with bowl



Code: **Sun-TBSV02**
Double variable height
bowl stand with 2 bowls



Code: **Sun-TBSF01**
Two tier double bowl stand
with 2 bowls



Code: **Sun-TBSF02**
Side by side double bowl stand
with 2 bowls



Code: **Sun-TMAY02**
Mayo table, adjustable
height with 4 castors



Code: **Sun-TBKS01**
Bucket stand with bucket
18 litre capacity







Code: **Sun-TKAB01**
Kickabout stand with bowl

strong
secure
protection

Medicine & Controlled Drug Cabinets

Professionally designed and independently tested for your peace of mind, we provide secure storage solutions for most types of medicines and pharmaceuticals including the latest innovative BIDOSE system.

	PAGE
	Medicine Cabinets 108-111
	Controlled Drug Cabinets 112-115
	Medicine Cabinets with Controlled Drug Inner Cabinets 116-121
	MDS Medicine Cabinets 122-123
	Self Administration Cabinets 122-123
	BIDOSE Medicine Cabinets 124-125

product information



Independently tested



BS 2881:1989
Security Level 1
compliant



High security 8 lever
lock designed to
EN1300 and Vds 2396



Anti-pick levers &
deadlock mechanism



3 keys supplied
(per lock)
Extra keys available



Protective key
escutcheon



Adjustable shelves and
door trays




Internal door
retaining flange



Durable and hygienic
white epoxy coating



Single door cabinets
are right hand hinged
(left hand optional)



Optional warning light
and buzzer alarm



Easy clean design
and materials



2 year warranty



Supplied fully
assembled



British made



Delivery normally within
3-5 working days *
(Sun-MC1, MC2
& MC3 only)

* Please contact us regarding availability for quantities of 10 or more.

accessories * Factory fit - only available at time of ordering cabinet



Sun-MC/LIGHT1 *
(For Sun-MC1/MC2/MC3 Cabinets)
Sun-MC/LIGHT2 *
(For Sun-MC6/MC7 Cabinets)



Sun-MC/BUZZ *
Buzzer alarm (has to
be combined with
Sun-MC/Light)



Sun-MC/LHH *
Left hand hinge option



Sun-MC/KEY *
Additional key
(when ordered
with cabinet)

5V DC external warning light shows green when closed, red when open and flashes red when left open for 2.5 minutes. Nurse call output, (volt free contact) for connection to nurse call system. Internal illumination when door is opened. Must be hard wired to a 220/240V AC mains supply. Buzzer alarm activates when cabinet is left open for 5 minutes (only available when combined with optional warning light)

dimensions - NB: Optional warning light transformer box adds 55mm to overall height



MC1	MC2	MC3	MC6	MC7
60cm(H)	60cm(H)	60cm(H)	85cm(H)	85cm(H)
30cm(W)	50cm(W)	60cm(W)	50cm(W)	60cm(W)
30cm(D)	30cm(D)	30cm(D)	30cm(D)	30cm(D)

Independently tested to conform to or exceed all the required standards and legislation, all our medicine cabinets provide secure storage, which can be adapted internally to your exact requirements.



Code: **Sun-MC1**
Medicine cabinet with 3 shelves, 2 door trays,
shown with optional warning light



Code: **Sun-MC2**
Medicine cabinet with 3 shelves, 2 door trays,
shown with optional warning light



Code: **Sun-MC3**
Medicine cabinet with 3 shelves, 2 door trays,
shown with optional warning light



Code: **Sun-MC6**
Medicine cabinet with 4 shelves, 4 door trays,
shown with optional warning light



Code: **Sun-MC7**
Medicine cabinet with 4 shelves, 4 door trays,
shown with optional warning light



product information



Independently tested



BS 2881:1989
Security Level 1
compliant



High security 8 lever
lock designed to
EN1300 and Vds 2396



Anti-pick levers &
deadlock mechanism



3 keys supplied
(per lock)
Extra keys available



Protective key
escutcheon



Adjustable shelves and
door trays




Internal door
retaining flange



Durable and hygienic
white epoxy coating



Single door cabinets
are right hand hinged
(left hand optional)



Optional warning light
and buzzer alarm



Easy clean design
and materials



2 year warranty



Supplied fully
assembled



British made

accessories * Factory fit - only available at time of ordering cabinet



Sun-MC/LIGHT3 *
(For Sun-MC4/MC5/MC8/MC9 Cabinets)



Sun-MC/BUZZ *
Buzzer alarm (has to
be combined with
Sun-MC/Light)



Sun-MC/KEY *
Additional key
(when ordered
with cabinet)

5V DC external warning light shows green when closed, red when open and flashes red when left open for 2.5 minutes. Nurse call output, (volt free contact) for connection to nurse call system. Internal illumination when door is opened. Must be hard wired to a 220/240V AC mains supply. Buzzer alarm activates when cabinet is left open for 5 minutes (only available when combined with optional warning light)


dimensions - NB: Optional warning light transformer box adds 55mm to overall height



60cm(H)
80cm(W)
30cm(D)



60cm(H)
100cm(W)
30cm(D)



85cm(H)
80cm(W)
30cm(D)



85cm(H)
100cm(W)
30cm(D)

Independently tested to conform to or exceed all the required standards and legislation, all our medicine cabinets provide secure storage, which can be adapted internally to your exact requirements.



Code: **Sun-MC4**
Medicine cabinet with 6 shelves, 4 door trays,
shown with optional warning light



Code: **Sun-MC5**
Medicine cabinet with 6 shelves, 4 door trays,
shown with optional warning light



Code: **Sun-MC8**
Medicine cabinet with 8 shelves, 8 door trays,
shown with optional warning light




Code: **Sun-MC9**
Medicine cabinet with 8 shelves, 8 door trays,
shown with optional warning light



product information



Independently tested



Meets the Misuse of Drugs (Safe Custody Regulation 1973)



BS 2881:1989 Security Level 1 compliant



High security 8 lever lock designed to EN1300 and Vds 2396



Anti-pick levers & deadlock mechanism



3 keys supplied (per lock)
Extra keys available



Protective key escutcheon



Adjustable shelves and door trays



Internal door retaining flange



Durable and hygienic white epoxy coating



Single door cabinets are right hand hinged (left hand optional)



Optional warning light and buzzer alarm



Easy clean design and materials



2 year warranty



Supplied fully assembled




British made



Delivery normally within 3-5 working days *
(Sun-CDC21, 22 & 23 only)

* Please contact us regarding availability for quantities of 10 or more.

accessories * Factory fit - only available at time of ordering cabinet



Sun-CDC/LIGHT1 *
(For Sun-CDC21/CDC22/CDC23 Cabinets)
Sun-CDC/LIGHT2 *
(For Sun-CDC24/CDC25 Cabinets)



Sun-CDC/BUZZ *
Buzzer alarm (has to be combined with Sun-CDC/Light)



Sun-CDC/LHH *
Left hand hinge option



Sun-CDC/KEY *
Additional key (when ordered with cabinet)

5V DC external warning light shows green when closed, red when open and flashes red when left open for 2.5 minutes. Nurse call output, (volt free contact) for connection to nurse call system. Internal illumination when door is opened. Must be hard wired to a 220/240V AC mains supply. Buzzer alarm activates when cabinet is left open for 5 minutes (only available when combined with optional warning light)

dimensions - NB: Optional warning light transformer box adds 55mm to overall height

CDC21	CDC22	CDC23	CDC24	CDC25
30cm(H) 21cm(W) 27cm(D)	30cm(H) 34cm(W) 27cm(D)	55cm(H) 34cm(W) 27cm(D)	55cm(H) 50cm(W) 30cm(D)	60cm(H) 60cm(W) 30cm(D)

Controlled Drug Cabinets

Independently tested to conform to or exceed all the required standards and legislation, all our controlled drug cabinets provide secure storage, which can be adapted internally to your exact requirements.



Code: **Sun-CDC21**
Controlled drug cabinet with 1 shelf



Code: **Sun-CDC22**
Controlled drug cabinet with 1 shelf and 1 door tray, shown with optional warning light



Code: **Sun-CDC23**
Controlled drug cabinet with 2 shelves and 2 door trays, shown with optional warning light



Code: **Sun-CDC24**
Controlled drug cabinet with 2 shelves and 2 door trays, shown with optional warning light



Code: **Sun-CDC25**
Controlled drug cabinet with 2 shelves, 2 door trays and 2 locks, shown with optional warning light

Code: **Sun-CDC23**

product information



Independently tested



Meets the Misuse of Drugs (Safe Custody Regulation 1973)



BS 2881:1989 Security Level 1 compliant



High security 8 lever lock designed to EN1300 and Vds 2396



Anti-pick levers & deadlock mechanism



3 keys supplied (per lock)
Extra keys available



Protective key escutcheon



Adjustable shelves and door trays



Internal door retaining flange



Durable and hygienic white epoxy coating



Single door cabinets are right hand hinged (left hand optional)



Optional warning light and buzzer alarm



Easy clean design and materials



2 year warranty




Supplied fully assembled



British made

accessories * Factory fit - only available at time of ordering cabinet



Sun-CDC/LIGHT2 *
(For Sun-CDC26/CDC28 Cabinets)
Sun-CDC/LIGHT3 *
(For Sun-CDC27 Cabinet)



Sun-CDC/BUZZ *
Buzzer alarm (has to be combined with Sun-CDC/Light)



Sun-CDC/LHH *
Left hand hinge option



Sun-CDC/KEY *
Additional key (when ordered with cabinet)

5V DC external warning light shows green when closed, red when open and flashes red when left open for 2.5 minutes. Nurse call output, (volt free contact) for connection to nurse call system. Internal illumination when door is opened. Must be hard wired to a 220/240V AC mains supply. Buzzer alarm activates when cabinet is left open for 5 minutes (only available when combined with optional warning light)


dimensions - NB: Optional warning light transformer box adds 55mm to overall height

CDC26



85cm(H)
50cm(W)
30cm(D)

CDC27



85cm(H)
100cm(W)
30cm(D)

CDC28



85cm(H)
50cm(W)
45cm(D)

Controlled Drug Cabinets

Independently tested to conform to or exceed all the required standards and legislation, all our controlled drug cabinets provide secure storage, which can be adapted internally to your exact requirements.



Code: **Sun-CDC26**
Controlled drug cabinet with 4 shelves and 4 door trays, shown with optional warning light



Code: **Sun-CDC28**
Controlled drug cabinet with 4 shelves and 4 door trays, shown with optional warning light



Code: **Sun-CDC27**
Controlled drug cabinet with 8 shelves and 8 door trays, shown with optional warning light



product information



Independently tested



Meets the Misuse of Drugs (Safe Custody Regulation 1973)



BS 2881:1989 Security Level 1 compliant



High security 8 lever lock designed to EN1300 and Vds 2396



Anti-pick levers & deadlock mechanism



3 keys supplied (per lock)
Extra keys available



Protective key escutcheon on outer cabinet



5V External warning light is green when closed, red when open - fitted as standard on outer medicine cabinet



Optional buzzer alarm



Adjustable shelves and door trays



Identification coloured key fobs



Identification key escutcheon on inner cabinet



Internal door retaining flange



Durable and hygienic white epoxy coating



Single door cabinets are right hand hinged (left hand optional)



Easy clean design and materials



2 year warranty




Supplied fully assembled




British made


accessories * Factory fit - only available at time of ordering cabinet




Sun-MC/BUZZ *
Buzzer alarm



Sun-MC/LHH *
Left hand hinge option










Sun-MC/KEY *
Additional key for outer cabinet (when ordered with cabinet)



Sun-CDC/KEY *
Additional key for controlled drug inner cabinet (when ordered with cabinet)

Warning Light (fitted as standard)
5V DC external warning light shows green when closed, red when open and flashes red when left open for 2.5 minutes. Nurse call output, (volt free contact) for connection to nurse call system. Internal illumination when door is opened. Must be hard wired to a 220/240V AC mains supply. Buzzer alarm activates when cabinet is left open for 5 minutes (only available when combined with optional warning light)

dimensions - NB: Height measurement for outer cabinets includes transformer box

Controlled Drug Inner Cabinets			Medicine Outer Cabinets			
CDC21	CDC22	CDC23	MC2	MC3	MC6	MC7
						
30cm(H) 21cm(W) 27cm(D)	30cm(H) 34cm(W) 27cm(D)	55cm(H) 34cm(W) 27cm(D)	66cm(H) 50cm(W) 30cm(D)	66cm(H) 60cm(W) 30cm(D)	91cm(H) 50cm(W) 30cm(D)	91cm(H) 60cm(W) 30cm(D)

Independently tested to conform to or exceed all the required standards and legislation. Our range of combined medicine and drug cabinets provide secure storage, which can be adapted internally to your exact requirements. A warning light and buzzer provide additional security.



CODE: **Sun-MCDC221**
MC2 Medicine outer cabinet combined with CDC21 controlled drugs inner cabinet



CODE: **Sun-MCDC322**
MC3 Medicine outer cabinet combined with CDC22 controlled drugs inner cabinet



CODE: **Sun-MCDC321**
MC3 Medicine outer cabinet combined with CDC21 controlled drugs inner cabinet



CODE: **Sun-MCDC621**
MC6 Medicine outer cabinet combined with CDC21 controlled drugs inner cabinet



CODE: **Sun-MCDC723**
MC7 Medicine outer cabinet combined with CDC23 controlled drugs inner cabinet



CODE: **Sun-MCDC722**
MC7 Medicine outer cabinet combined with CDC22 controlled drugs inner cabinet

product information



Independently tested



Meets the Misuse of Drugs (Safe Custody Regulation 1973)



BS 2881:1989 Security Level 1 compliant



High security 8 lever lock designed to EN1300 and Vds 2396



Anti-pick levers & deadlock mechanism



3 keys supplied (per lock)
Extra keys available



Protective key escutcheon on outer cabinet



5V External warning light is green when closed, red when open - fitted as standard on outer medicine cabinet



Optional buzzer alarm



Adjustable shelves and door trays



Identification coloured key fobs



Identification key escutcheon on inner cabinet



Internal door retaining flange



Durable and hygienic white epoxy coating



Easy clean design and materials



2 year warranty




Supplied fully assembled




British made


accessories * Factory fit - only available at time of ordering cabinet



Sun-MC/BUZZ *
Buzzer alarm









Sun-MC/KEY *
Additional key for outer cabinet (when ordered with cabinet)



Sun-CDC/KEY *
Additional key for controlled drug inner cabinet (when ordered with cabinet)

Warning Light (fitted as standard)
5V DC external warning light shows green when closed, red when open and flashes red when left open for 2.5 minutes. Nurse call output, (volt free contact) for connection to nurse call system. Internal illumination when door is opened. Must be hard wired to a 220/240V AC mains supply. Buzzer alarm activates when cabinet is left open for 5 minutes (only available when combined with optional warning light)

dimensions - NB: Height measurement for outer cabinets includes transformer box

Controlled Drug Inner Cabinets			Medicine Outer Cabinets		
CDC21	CDC22	CDC23	MC4	MC5	MC8
					
30cm(H) 21cm(W) 27cm(D)	30cm(H) 34cm(W) 27cm(D)	55cm(H) 34cm(W) 27cm(D)	66cm(H) 80cm(W) 30cm(D)	66cm(H) 100cm(W) 30cm(D)	91cm(H) 80cm(W) 30cm(D)

Independently tested to conform to or exceed all the required standards and legislation. Our range of combined medicine and drug cabinets provide secure storage, which can be adapted internally to your exact requirements. A warning light and buzzer provide additional security.



CODE: **Sun-MCDC421**
MC4 Medicine outer cabinet combined with CDC21 controlled drugs inner cabinet



CODE: **Sun-MCDC521**
MC5 Medicine outer cabinet combined with CDC21 controlled drugs inner cabinet



CODE: **Sun-MCDC522**
MC5 Medicine outer cabinet combined with CDC22 controlled drugs inner cabinet



CODE: **Sun-MCDC821**
MC8 Medicine outer cabinet combined with CDC21 controlled drugs inner cabinet

product information



Independently tested



Meets the Misuse of Drugs (Safe Custody Regulation 1973)



BS 2881:1989 Security Level 1 compliant



High security 8 lever lock designed to EN1300 and Vds 2396



Anti-pick levers & deadlock mechanism



3 keys supplied (per lock)
Extra keys available



Protective key escutcheon on outer cabinet



5V External warning light is green when closed, red when open - fitted as standard on outer medicine cabinet



Optional buzzer alarm



Adjustable shelves and door trays



Identification coloured key fobs



Identification key escutcheon on inner cabinet



Internal door retaining flange



Durable and hygienic white epoxy coating



Easy clean design and materials



2 year warranty




Supplied fully assembled




British made


accessories * Factory fit - only available at time of ordering cabinet



Sun-MC/BUZZ *
Buzzer alarm



Sun-MC/KEY *
Additional key for outer cabinet (when ordered with cabinet)



Sun-CDC/KEY *
Additional key for controlled drug inner cabinet (when ordered with cabinet)

Warning Light (fitted as standard)
5V DC external warning light shows green when closed, red when open and flashes red when left open for 2.5 minutes. Nurse call output, (volt free contact) for connection to nurse call system. Internal illumination when door is opened. Must be hard wired to a 220/240V AC mains supply. Buzzer alarm activates when cabinet is left open for 5 minutes (only available when combined with optional warning light)

dimensions - NB: Height measurement for outer cabinets includes transformer box

Controlled Drug Inner Cabinets

CDC21



30cm(H)
21cm(W)
27cm(D)

CDC22



30cm(H)
34cm(W)
27cm(D)

CDC23



55cm(H)
34cm(W)
27cm(D)

MC9



91cm(H)
100cm(W)
30cm(D)

Independently tested to conform to or exceed all the required standards and legislation. Our range of combined medicine and drug cabinets provide secure storage, which can be adapted internally to your exact requirements. A warning light and buzzer provide additional security.



CODE: **Sun-MCDC921**
MC9 Medicine outer cabinet combined with CDC21 controlled drugs inner cabinet



CODE: **Sun-MCDC922**
MC9 Medicine outer cabinet combined with CDC22 controlled drugs inner cabinet



CODE: **Sun-MCDC923**
MC9 Medicine outer cabinet combined with CDC23 controlled drugs inner cabinet




product information




Independently tested



BS 2881:1989
Security Level 1
compliant




High security 8 lever
lock designed to
EN1300 and Vds 2396




Anti-pick levers &
deadlock mechanism




3 keys supplied
(per lock)
Extra keys available




Protective key
escutcheon



Integral brackets for
hanging MDS racks



Internal door
retaining flange




Durable and hygienic
white epoxy coating



Single door cabinets
are right hand hinged
(left hand optional)



Optional warning light
and buzzer alarm




Easy clean design
and materials



2 Year Warranty



Supplied fully
assembled



British made

accessories * Factory fit - only available at time of ordering cabinet



Sun-MC/LIGHT1 *
(For Sun-MC1/MC2/MC3 Cabinets)
Sun-MC/LIGHT2 *
(For Sun-MC6/MC7 Cabinets)
Sun-MC/LIGHT3 *
(For Sun-MC4/MC5/MC8/MC9 Cabinets)



Sun-MC/BUZZ *
Buzzer alarm (has to
be combined with
Sun-MC/Light)




Sun-MC/LHH *
Left hand hinge option



Sun-MC/KEY *
Additional key
(when ordered
with cabinet)


5V DC external warning light shows green when closed, red when open and flashes red when left open for 2.5 minutes. Nurse call output, (volt free contact) for connection to nurse call system. Internal illumination when door is opened. Must be hard wired to a 220/240V AC mains supply. Buzzer alarm activates when cabinet is left open for 5 minutes (only available when combined with optional warning light)

dimensions - NB: Optional warning light transformer box adds 55mm to overall height




MC1

60cm(H)
30cm(W)
30cm(D)




MC2

60cm(H)
50cm(W)
30cm(D)




MC3

60cm(H)
60cm(W)
30cm(D)



MC4


60cm(H)
80cm(W)
30cm(D)



MC5


60cm(H)
100cm(W)
30cm(D)

Self Administration Cabinets




MC6

85cm(H)
50cm(W)
30cm(D)




MC7

85cm(H)
60cm(W)
30cm(D)




MC8

85cm(H)
80cm(W)
30cm(D)




MC9

85cm(H)
100cm(W)
30cm(D)



MDSCAB1

42cm(H)
30cm(W)
24cm(D)



MDSCAB2

72cm(H)
40cm(W)
24cm(D)

Independently tested to conform to or exceed all the required standards and legislation, our range of Monitored Dosage System (MDS) cabinets provide secure storage, and designed for nursing and care homes using the Manrex (Boots) and Venalink systems.




Code: **Sun-MC1/MDS1** Medicine cabinet with 1 MDS rack
Code: **Sun-MC2/MDS2** Medicine cabinet with 2 MDS racks
Code: **Sun-MC3/MDS2** Medicine cabinet with 2 MDS racks
Code: **Sun-MC6/MDS4** Medicine cabinet with 4 MDS racks
Code: **Sun-MC7/MDS4** Medicine cabinet with 4 MDS racks
(shown with optional warning light & left hand hinge option)



Code: **Sun-MC4/MDS2** Medicine cabinet with 2 MDS racks
Code: **Sun-MC5/MDS4** Medicine cabinet with 4 MDS racks
Code: **Sun-MC8/MDS4** Medicine cabinet with 4 MDS racks
Code: **Sun-MC9/MDS8** Medicine cabinet with 8 MDS racks
(shown with optional warning light)

Self-Administration Cabinets


Designed for storage of self administration medicines, our cabinets provide low level security for the Monitored Dosage System (MDS) and are suitable for storing the Manrex (Boots) and Venalink systems. Optional keyed to differ available.



Code: **Sun-MDSCAB1**
Code: **Sun-MDSCAB1/KD (keyed to differ)**
Self administration MDS beech effect locking wall cabinet with 1 shelf & 1 rack bracket (supplied with 2 keys)



Code: **Sun-MDSCAB2**
Code: **Sun-MDSCAB2/KD (keyed to differ)**
Self administration MDS beech effect locking wall cabinet with 1 shelf & 2 rack brackets (supplied with 2 keys)



* Please contact us regarding availability for quantities of 10 or more.



FSC
www.fsc.org
FSC® C123891
The mark of responsible forestry



Delivery normally within
3-5 working days *

product information - BODOSE medicine cabinets



Independently tested

15 or 30 trays

Cabinets store 15 or 30 BODOSE trays



BS 2881:1989 Security Level 1 compliant



Individual tray slots with easy access



Protective key escutcheon



High security 8 lever lock designed to EN1300 and Vds 2396



Anti-pick levers & deadlock mechanism



3 keys supplied (per lock) Extra keys available



Internal door retaining flange




Mild steel with hygienic white epoxy coating


product information - BODOSE self administration medicine cabinets

6


Holds 6 BODOSE trays




Individual tray slots with easy access



Lock provided with 2 keys




Mild steel with hygienic white epoxy coating (Sun-BIO/SA6 only)




Oak, beech or maple effect finishes


additional product information




Easy clean design and materials



2 year warranty




Supplied fully assembled




British made

accessories (only available for Sun-BIO/MC15 & Sun-BIO/MC30 cabinets)


* Factory fit - only available at time of ordering cabinet




Sun-BIO/MC/LIGHT *
External Warning Light (see below)



Sun-BIO/MC/BUZZ *
Buzzer Alarm (has to be combined with warning light)



Sun-BIO/MC/KEY *
Additional key (when ordered with cabinet)




BIODOSE Drug Trolley Range
(see pages 221-222)

5V DC external warning light shows green when closed, red when open and flashes red when left open for 2.5 minutes. Nurse call output, (volt free contact) for connection to nurse call system. Internal illumination when door is opened. Must be hard wired to a 220/240V AC mains supply. Buzzer alarm activates when cabinet is left open for 5 minutes (only available when combined with optional warning light)


dimensions - NB: Optional warning light transformer box adds 55mm to overall height

Sun-BIO/MC15




73cm(H)
32cm(W)
39cm(D)

Sun-BIO/MC30




73cm(H)
55cm(W)
39cm(D)

Sun-BIO/SA6



32cm(H)
43cm(W)
25cm(D)

**Sun-BIO/WSA6/OAK
Sun-BIO/WSA6/MAPLE
Sun-BIO/WSA6/BEECH**



40cm(H)
39cm(W)
26cm(D)

Designed and independently tested to conform to or exceed all the required standards and legislation. These cabinets provide secure storage with easy access to BODOSE medicine trays.



Code: **Sun-BIO/MC15**
BIODOSE medicine cabinet, 15 tray capacity



Code: **Sun-BIO/MC30**
BIODOSE medicine cabinet, 30 tray capacity

BIODOSE® Self Administration Medicine Cabinets

Designed specifically for the self-administration and storage of BODOSE medicine trays. These cabinets provide low level security, whilst being aesthetically pleasing and easy to clean.



Code: **Sun-BIO/SA6**
BIODOSE self administration mild steel medicine cabinet, 6 tray capacity



Code: **Sun-BIO/WSA6/MAPLE**
BIODOSE self administration medicine cabinet, 6 tray capacity - maple effect finish



Code: **Sun-BIO/WSA6/OAK**
BIODOSE self administration medicine cabinet, 6 tray capacity - oak effect finish



Code: **Sun-BIO/WSA6/BEECH**
BIODOSE self administration medicine cabinet, 6 tray capacity - beech finish



simple
hygienic
disposal

Bins & Sack Holders

Ergonomic and hygienic design, combined with strength and functionality is what makes our bins unique. Compliant to the latest Safe Management of Waste guidelines with fire retardancy, hands free operation and silent closing all as standard.



Bins

PAGE

127-128



Sack Holders

129-130

product information

3

sizes

20, 50 & 70 litre capacity

hands free

Hygienic hands free operation, foot pedal must be used to open

HTM07-01 compliant

Safe management of healthcare waste compliant

HTM05-03 tested

Designed & tested to HTM05-03 Part A

Silent, dampened closing lid as standard

Corrosion resistant epoxy coated white galvanised bodies

Silicone bag retaining cord with nylon finger loop

Stainless steel foot pedal with pedal label (where applicable)

Sealed base prevents leakages

wall mounting bracket

Integrated wall mounting bracket

Off the wall design to help prevent damage

Flame retardant **

Easy clean design and materials

2 year warranty

Delivery normally within 3-5 working days *

** Designed and tested to HTM05-03 Part A * Please contact us regarding availability for quantities of 10 or more.

designated usage - please specify lid colour when ordering

<div>WHITE</div> <div></div> <div>General Use</div>	<div>YELLOW</div> <div></div> <div>Waste for Incineration</div>	<div>BLACK</div> <div></div> <div>Domestic (Municipal) Waste</div>	<div>YELLOW/BLK</div> <div></div> <div>Offensive/ Hygiene Waste</div>	<div>ORANGE</div> <div></div> <div>Waste Which May Be Treated</div>
<div>PURPLE</div> <div></div> <div>Cytotoxic & Cytostatic Waste</div>	<div>AMALGAM</div> <div></div> <div>Amalgam Waste</div>	<div>BLUE</div> <div></div> <div>Medical Waste for Incineration</div>	<div>RED</div> <div></div> <div>Anatomical Waste for Incineration</div>	

dimensions

20 Litre

51cm(H)
32cm(W)
33cm(D)

50 Litre

59cm(H)
42cm(W)
43cm(D)

70 Litre

78cm(H)
42cm(W)
43cm(D)

We are "The Strong Silent Type" and we don't like hands either! Designed for the multi function disposal of waste (HTM07-01), with usage designated by an appropriate coloured lid and pedal label.



Code: Sun-BIN20 - 20 Litre bin



Code: Sun-BIN50 - 50 Litre bin



Code: Sun-BIN70 - 70 Litre bin

product information - sack holders

3
sizes

20, 70 & 90
litre capacity

Hygienic operation,
foot pedal should be
used to open

Elasticated bag
retaining cord with
nylon finger loop

Corrosion resistant
materials

HTM07-01
compliant

Safe management
of healthcare
waste compliant

Easy clean design
and materials

1 year warranty

Supplied fully
assembled

British made

Delivery normally within
3-5 working days *

* Please contact us regarding availability for quantities of 10 or more.

designated usage - please specify lid colour when ordering

General Use	Waste for Incineration	Domestic (Municipal) Waste	Offensive/ Hygiene Waste	Waste Which May Be Treated
Cytotoxic & Cytostatic Waste	Amalgam Waste	Medical Waste for Incineration	Anatomical Waste for Incineration	

dimensions

20 Litre

58cm(H)
22cm(W)
35cm(D)

70 Litre

78cm(H)
43cm(W)
44cm(D)

90 Litre

86cm(H)
43cm(W)
44cm(D)

Hands-free operation and sturdy design (HTM07-01) is what makes our sack holders hygienic and simple to use. They are ideal for areas where the risk of fire is minimal. Usage is designated by an appropriate coloured lid.



Code: **Sun-BINSH20** - 20 Litre sack holder



Code: **Sun-BINSH70** - 70 Litre sack holder




Code: **Sun-BINSH90** - 90 Litre sack holder

hygienic
mobile
privacy

Mobile Screens

Smart design combined with ease of use is what makes our hygienic range of mobile screens unique. Choose from either our solid panel system or disposable one-piece curtain range, for ease of cleaning or to assist in reducing cross infection.



new

Mobile Screens - Solid Panels

PAGE

132-133



Mobile Screens - Disposable Curtains

134-135

product information

3/4/5
panel

Available in 3, 4 and 5 panel sizes

8
designs

Unique range of 8 patterns and colours

Hygienic PVC panel system

Panels can be removed or replaced

Flame retardant panels - BS476 part 7

Multiple leg supports for stability

Swivel castors for ease of movement

Durable epoxy coated white frame

Folds flat for storage

Easy clean design and materials

Minimal assembly required

2 year warranty

British made

dimensions

3 Panel
172cm(H) 188cm(W) 3cm(D)
38cm(D) Across Legs

4 Panel
172cm(H) 250cm(W) 3cm(D)
38cm(D) Across Legs

5 Panel
172cm(H) 312cm(W) 3cm(D)
38cm(D) Across Legs

screen panel patterns and colours - please specify when ordering

white

beige

silver

sky blue

blue patchwork

bubble design

wave design

jungle design

Designed to provide a temporary or mobile privacy solution, our range of folding screens incorporate our unique hygienic panel system in 8 patterns and colours and are functional, durable and easy to clean.




Code: **Sun-MFS3** - Three panel screen - Wave Design



Code: **Sun-MFS4** - Four panel screen - Jungle Design



Code: **Sun-MFS5** - Five panel screen - Bubble Design



product information

3/4/5
sections

Available in 3, 4 and 5 sections

4
colours

Choice of 4 curtain colours

££££
SAVING

Cost effective disposable curtain system

unique

Unique design

No laundry costs

Curtain labelling system for monitoring change dates

One piece pleated curtain - easy to fit or replace

Flame retardant curtains BS5867 part 2 type C

Infection control friendly

Recyclable polypropylene material

Multiple leg supports for stability

Durable epoxy coated white frame

Swivel castors for ease of movement

Folds flat for storage

Easy clean design and materials

minimal

Minimal assembly required

2 year warranty on frame

British made

Environmentally Friendly

Designed in partnership with Marlux Medical

accessories - replacement curtains

pastel blue

Sun-CUR3/DCPB - 3 Panel
Sun-CUR4/DCPB - 4 Panel
Sun-CUR5/DCPB - 5 Panel

summer blue

Sun-CUR3/DCSB - 3 Panel
Sun-CUR4/DCSB - 4 Panel
Sun-CUR5/DCSB - 5 Panel

forest green

Sun-CUR3/DCFG - 3 Panel
Sun-CUR4/DCFG - 4 Panel
Sun-CUR5/DCFG - 5 Panel

beige

Sun-CUR3/DCB - 3 Panel
Sun-CUR4/DCB - 4 Panel
Sun-CUR5/DCB - 5 Panel

dimensions

3 Panel
172cm(H) 188cm(W) 3cm(D)
38cm(D) Across Legs

4 Panel
172cm(H) 250cm(W) 3cm(D)
38cm(D) Across Legs

5 Panel
172cm(H) 312cm(W) 3cm(D)
38cm(D) Across Legs

curtain colours - please specify when ordering

pastel blue

summer blue

forest green

beige

Designed to provide a temporary or mobile privacy solution, and assist in reducing cross infection, our range of folding screens incorporate our unique disposable curtain system in 4 colours and are functional, durable and easy to clean.



Code: **Sun-MFS3/DC**
Three section screen with disposable curtain



Code: **Sun-MFS4/DC**
Four section screen with disposable curtain



Code: **Sun-MFS5/DC**
Five section screen with disposable curtain



individual
practical
solutions

Mirrors, Height Measure & Mobile Drip Stands

Practical solutions for general medical environments which are designed to assist both the patient and healthcare professional. Durable, functional and easy to clean and intended to provide a safer, more efficient and hygienic environment.



Medical Safety Mirrors

PAGE

137-138



Mobile Drip Stands

139-140


product information



Durable shatter proof design



5mm aircraft grade "plaskolite" acrylic



Small (A5) to extra large



Wall mounting tape supplied



Smooth rounded corners



Easy clean design and materials



2 year warranty



British made



Delivery normally within 3-5 working days *

* Please contact us regarding availability for quantities of 10 or more.

dimensions



Mirror 1
210mm(H)
148mm(W)
5mm(D)
(A5 Size)




Mirror 2
297mm(H)
210mm(W)
5mm(D)
(A4 Size)



Mirror 3
420mm(H)
297mm(W)
5mm(D)
(A3 Size)



Mirror 4
600mm(H)
350mm(W)
5mm(D)




Mirror 5
1200mm(H)
350mm(W)
5mm(D)

Medical Safety Mirrors


Smart design and utilisation of modern materials is what makes our medical safety mirrors unique. Using shatter proof material as a modern alternative to glass makes them an ideal solution for medical rooms, public areas or high risk environments.




Code: **Sun-MIRROR1** (A5 Size)




Code: **Sun-MIRROR2** (A4 Size)



Code: **Sun-MIRROR3** (A3 Size)



Code: **Sun-MIRROR4**



Code: **Sun-MIRROR5**

product information

<p>HEIGHT RANGE 160-245cm</p> <p>Sun-IV01, IV02</p>	<p>HEIGHT RANGE 156-267cm</p> <p>Sun-IV08 & IV09</p>	<p>Simple locking height adjustment</p>	<p>Complies with BS3619:1976</p>	<p>Chrome or stainless steel poles</p>
<p>2kg hook load or 4kg total</p>	<p>Lightweight design approx 3kg (Sun-IV01, IV02)</p>	<p>Stable 5 star Plastic castor base (Sun-IV01 & IV02)</p>	<p>5 Star Weighted base (Sun-IV08 & IV09)</p>	<p>Antistatic castors (Sun-IV08, IV09)</p>
<p>Simple self assembly required</p>	<p>Easy clean design and materials</p>	<p>2 year warranty</p>	<p>British made</p>	<p>Delivery normally within 3-5 working days *</p>

* Please contact us regarding availability for quantities of 10 or more.

dimensions

<p>Sun-IV01 160cm - 245cm(H) 51cm(W)</p>	<p>Sun-IV02 160cm - 245cm(H) 51cm(W)</p>	<p>Sun-IV08 156cm - 267cm(H) 58cm(W)</p>	<p>Sun-IV09 156cm - 267cm(H) 58cm(W)</p>
---	---	---	---


Simple, practical, easy to adjust and manoeuvre, our drip stands are suitable for most healthcare environments. We can provide a lightweight design through to ones with heavy duty weighted bases, available in a choice of materials to suit your budget and hygiene requirements.




Code: **Sun-IV01**
Chrome steel drip stand, plastic base, 2 plastic hooks



Code: **Sun-IV02**
Chrome steel drip stand, plastic base, 4 chrome hooks



Code: **Sun-IV08**
Brushed stainless steel drip stand and weighted base, 2 chrome hooks, antistatic castors



Code: **Sun-IV09**
Brushed stainless steel drip stand and weighted base, 4 chrome hooks, antistatic castors

individual
flexible
solutions

Furniture Packages

A choice of furniture packages to suit your budget, individual requirements and environment.



Furniture Packages

PAGE

142-143

furniture package one

Colenso examination couch (navy upholstery*) with paper roll holder
Single tier step (rectangular)
3 panel mobile folding screen - white panels
46cm wide dressing trolley with two glass effect safety trays
Medium first aid cabinet with 2 shelves
20 Litre bin - white



All additional room contents are for display purposes only

product information - furniture package two

Practitioner deluxe examination couch (navy upholstery*) with paper roll holder
Single couch drawer unit in white
Single tier step (rectangular)
4 Panel mobile folding screen - silver panels
46cm wide dressing trolley with three glass effect safety trays and drawer unit in white
Medium first aid cabinet with 2 shelves
20 Litre bin - white



All additional room contents are for display purposes only

product information - furniture package three

Practitioner deluxe examination couch (navy upholstery*) with paper roll holder
Three couch drawer units in white
Single tier step (oval)
4 Panel mobile folding screen - silver panels
Vista 30 trolley with 2 single and 2 double clear trays
50cm Wide locking base cabinet - white
50cm Wide locking wall cabinet - white
20 Litre bin - white



All additional room contents are for display purposes only

product information - furniture package four

Specialist treatment couch (grey upholstery*) with paper roll holder and one drawerline unit & one drawer pack in white
Double tier step (oval)
4 Panel mobile folding screen - silver panels
Vista 50 trolley with 4 single and 1 double clear tray
100cm Wide locking base cabinet - white
100cm Wide locking wall cabinet - white
50 Litre bin - white



All additional room contents are for display purposes only

product information - furniture package five
(school/restroom)

Alberti rest couch (sky blue upholstery*) with paper roll holder
4 Panel mobile folding screen - bubble design panels
VISTA 50 trolley with 6 single blue trays
50cm Wide locking base cabinet & 50cm Wide locking wall cabinet - white
20 Litre bin - white



All additional room contents are for display purposes only
* Other upholstery colours available, please see inside back cover
** Please specify at time of ordering



FIRA Standard 6250 2005
Specification:
Adhesion Performance
General/Severe



Anti-bacterial Vinyl Strength Test
BS EN ISO 13934-1: 2013
BS EN ISO 13973-3: 2000
BS EN ISO 12947-2: 1998



Fire (Flammability) Testing
Flammability BS 7176 Medium Hazard
BS 5852 : 2006 Clause 11
(composites)



Sun-INST2**
Assembly Service &
Packaging Recycling

Code: Sun-FP1
Furniture package one



Code: Sun-FP2
Furniture package two



Code: Sun-FP3
Furniture package three



Code: Sun-FP4
Furniture package four



Code: Sun-FP5
Furniture package five



stylish
modular
storage

Cabinet Ranges

Aesthetically designed, contemporary range of cabinetry incorporating wall, base, corner and sink units all available in five finishes. Modular sizes, adjustable shelving and locks as standard.



Beech Cabinet Range

PAGE

145-146



White (High Gloss) Cabinet Range

147-148

product information



Quality, 18mm MFC carcass with 2mm edging



Locking doors supplied with 2 keys



Contemporary door and drawer profile



Stylish flush fitting handles



38mm deep dark grey laminate work top



Base cabinets have one adjustable shelf, wall cabinets have two



Easy clean design and materials



FIRA Standard 6250 2005 Specification: Adhesion Performance General/Severe

- Base cabinets can be combined and supplied with a one piece work top in lengths of up to 4.1m
- Cut to length worktop can be supplied upon request

accessories



Sun-DT1 Drawer Insert for 50cm Wide Drawers



Integrated Plumbing System (IPS) (see pages 150-151)



Sinks (HTM64) (see pages 152-155)



Taps (HTM64) (see pages 156-159)




Sinks & Taps (General) (see pages 160-161)


dimensions

All single door base & wall cabinets are supplied right handed (unless otherwise specified)


Sun-BU8B 30cm Base Cabinet




Sun-BU1B 50cm Base Cabinet




Sun-BU2B 50cm Drawerline Cabinet




Sun-BU3B 50cm 5 Drawer Pack




Sun-BU4B 100cm Base Cabinet




Sun-BU5B 100cm Drawerline Cabinet



Sun-BU6B 100cm Sink Cabinet (excludes sink)




Sun-BU7B L-Shaped Corner Base Cabinet




Also available in 40cm, 60cm and 80cm widths


Sun-WU1B 30cm Wall Cabinet




Sun-WU2B 50cm Wall Cabinet



Sun-WU3B 100cm Wall Cabinet



Sun-WU4B Corner Wall Cabinet



Also available in 40cm, 60cm and 80cm widths

Aesthetically designed, contemporary range of modular cabinetry incorporating wall, base, corner and sink units finished in beech effect. Easy to clean, adjustable shelves and locks as standard.



White (High Gloss) Cabinet Range


product information



Quality, 18mm MFC carcass with 2mm edging




Locking doors supplied with 2 keys



Contemporary door and drawer profile



Stylish flush fitting handles



38mm deep crystal white laminate work top



Base cabinets have one adjustable shelf, wall cabinets have two



Easy clean design and materials



FIRA Standard 6250 2005 Specification: Adhesion Performance General/Severe

- Base cabinets can be combined and supplied with a one piece work top in lengths of up to 4.1m
- Cut to length worktop can be supplied upon request

accessories



Sun-DT1
Drawer Insert for
50cm Wide Drawers



Integrated Plumbing
System (IPS)
(see pages 150-151)



Sinks (HTM64)
(see pages 152-155)



Taps (HTM64)
(see pages 156-159)




Sinks & Taps (General)
(see pages 160-161)


dimensions

All single door base & wall cabinets are supplied right handed (unless otherwise specified)


Sun-BU8W 30cm
Base Cabinet




Sun-BU1W 50cm
Base Cabinet




Sun-BU2W 50cm
Drawerline Cabinet




Sun-BU3W 50cm
5 Drawer Pack




Sun-BU4W 100cm
Base Cabinet




Sun-BU5W 100cm
Drawerline Cabinet



Sun-BU6W 100cm
Sink Cabinet (excludes sink)




Sun-BU7W L-Shaped
Corner Base Cabinet




Also available in 40cm, 60cm and 80cm widths


Sun-WU1W 30cm
Wall Cabinet




Sun-WU2W 50cm
Wall Cabinet



Sun-WU3W 100cm
Wall Cabinet



Sun-WU4W
Corner Wall Cabinet



Also available in 40cm, 60cm and 80cm widths



White (High Gloss) Cabinet Range

Aesthetically designed, contemporary range of modular cabinetry incorporating wall, base, corner and sink units finished in high gloss white. Easy to clean, adjustable shelves and locks as standard.



Integrated Plumbing Systems, Sinks & Taps

Designed to assist in helping reduce cross infection and improve general hygiene, we can provide a range of HTM64 compliant clinical solutions from Integrated Plumbing Systems (IPS), sinks and taps for clean and dirty areas, as well as general hand washing applications.

	PAGE
 Integrated Plumbing Systems	150-151
 Sinks (HTM64 Compliant)	152-155
 Taps (HTM64 Compliant)	156-159
 Sinks & Taps (General Healthcare)	160-161

product information

HTM64
compliant

For use in clinical areas

Hygienic and durable design

Removable postformed panels

18mm
LFMDf

Laminate faced MDF panels

Optional sloping top for extra high ceilings

4 Widths
3 Heights

4 standard widths in 3 different heights

eze

Easy clean design and materials

FIRA
Independently tested by **FIRA**

FIRA Standard 6250 2005 Specification: Adhesion Performance General/Severe

installation

Professional installation available

accessories * Required for sink installation and for HTM64 compliance

Sun-SNK14
Washbasin - medium (see pages 154-155)

FIXING
kit

SINK
adapter

Sun-DISPENSER9
GOJO touch free dispenser for foaming soap

Sun-DISPENSER11
Purell touch free dispenser for hand rub gel

Sun-SNK15
Washbasin - large (see pages 154-155)

FIXING
kit

SINK
adapter

Sun-DISPENSER4
Paper towel dispenser

Medical Safety Mirrors
(see pages 137-138)

Sun-TAP15
Electronic thermostatic mixer tap with time flow sensor

Sun-TAP16
Electronic thermostatic mixer tap with proximity sensor

Sun-TAP17
Thermostatic sequential mixer tap

dimensions & codes

240cm(H)
IPS Panel System

240cm(H) 70cm(W) 32cm(D)

Sun-IPS1/LFMDf/240/70/WHITE

240cm(H) 80cm(W) 32cm(D)

Sun-IPS2/LFMDf/240/80/WHITE

240cm(H) 90cm(W) 32cm(D)

Sun-IPS3/LFMDf/240/90/WHITE

240cm(H) 100cm(W) 32cm(D)

Sun-IPS4/LFMDf/240/100/WHITE

250cm(H) 70cm(W) 32cm(D)

Sun-IPS1/LFMDf/250/70/WHITE

250cm(H) 80cm(W) 32cm(D)

Sun-IPS2/LFMDf/250/80/WHITE

250cm(H) 90cm(W) 32cm(D)

Sun-IPS3/LFMDf/250/90/WHITE

250cm(H) 100cm(W) 32cm(D)

Sun-IPS4/LFMDf/250/100/WHITE

260cm(H) 70cm(W) 32cm(D)

Sun-IPS1/LFMDf/260/70/WHITE

260cm(H) 80cm(W) 32cm(D)

Sun-IPS2/LFMDf/260/80/WHITE

260cm(H) 90cm(W) 32cm(D)

Sun-IPS3/LFMDf/260/90/WHITE

260cm(H) 100cm(W) 32cm(D)

Sun-IPS4/LFMDf/260/100/WHITE

240cm(H) 70cm(W) 32cm(D)

Sun-IPS8/LFMDf/240/70/WHITE

240cm(H) 80cm(W) 32cm(D)

Sun-IPS2/LFMDf/240/80/WHITE

240cm(H) 90cm(W) 32cm(D)

Sun-IPS3/LFMDf/240/90/WHITE

240cm(H) 100cm(W) 32cm(D)

Sun-IPS4/LFMDf/240/100/WHITE

Our range of HTM64 compliant integrated plumbing systems (IPS), combine a range of healthcare sanitary ware fitted into a specially designed pre-assembled panel arrangement. They are suitable for use in both clinical and personal washing areas and provide a cost-effective, functional and easy to maintain solution.




2 year warranty on panels




Sinks (HTM64 Compliant)


product information - Sun-SNK24 & 25




For use in clinical areas




No plug, chain or overflow




Available in left or right hand drainer




Premium quality stainless steel




200mm tap centres




Fixing kit supplied




Bottle trap and waste kit needed (see accessories)*



Easy clean design and materials



1 year warranty



British made

* Required for sink installation and for HTM64 compliance

product information - Sun-SNK28



For use in clinical areas



Premium quality stainless steel



No plug, chain or overflow



Fixing kit supplied



Bottle trap & waste kit needed (see accessories)*



Easy clean design and materials



1 year warranty



British made

* Required for sink installation and for HTM64 compliance

accessories - (please see pages 167-172 for our full range & specification of taps)



Sun-SNK24 & 25/WK*
Waste kit



Sun-SNK24 & 25/BT*
Bottle trap



Sun-SNK28/WK*
Waste kit



Sun-SNK28/BT*
Bottle trap




Sun-TAP13
Deck mounted sequential thermostatic mixer tap


* Required for sink installation and for HTM64 compliance

dimensions

Sun-SNK24 & Sun-SNK25



Sun-SNK28



Sinks (HTM64 Compliant)

Our range of HTM64 compliant inset stainless steel sinks and bowls are suitable for clean and dirty areas as well as clinical healthcare applications.

Inset Sinks

HTM64 compliant inset stainless steel sink, available with left or right hand drainer.



Code: Sun-SNK24 - Inset stainless steel sink with right hand drainer



Code: Sun-SNK25 - Inset stainless steel sink with left hand drainer



Code: Sun-SNK28 - Inset stainless steel hemispherical bowl

Sinks (HTM64 Compliant)

product information - Sun-SNK14 & 15

HTM64
compliant

For use in
clinical areas

No plug, chain or
overflow

Suitable for
IPS systems

Back
Outlet
Adapter

Back outlet adapter
needed (see accessories)*

1 Year Warranty

1 year warranty

product information - Sun-SNK16, 17 & 18

HTM64
compliant

Suitable for personal
washing only

No plug, chain or
overflow

Single tap hole
(right hand side)

Bottle Trap
&
Waste Kit

Bottle trap & waste kit
needed (see accessories)*

1 Year Warranty

1 year warranty

accessories

Panel
Fixing Kit

Sun-SNK14/KIT*
(for Sun-SNK14)
Sun-SNK15/KIT*
(for Sun-SNK15)

Back
Outlet
Adapter

Sun-SNK14/KIT1*
(for Sun-SNK14)
Sun-SNK15/KIT1*
(for Sun-SNK15)

Wall
Fixing Kit

Sun-SNK16/KIT1*
(for Sun-SNK16)
Sun-SNK17/KIT1*
(for Sun-SNK17)
Sun-SNK18/KIT1*
(for Sun-SNK18)

Panel
Fixing Kit

Sun-SNK16/KIT2*
(for Sun-SNK16)
Sun-SNK17/KIT2*
(for Sun-SNK17)
Sun-SNK18/KIT2*
(for Sun-SNK18)

Waste Kit

Sun-SNK16/WK*
(for Sun-SNK16)
Sun-SNK17/WK*
(for Sun-SNK17)
Sun-SNK18/WK*
(for Sun-SNK18)

Bottle Trap

Sun-SNK16/BT*
(for Sun-SNK16)
Sun-SNK17/BT*
(for Sun-SNK17)
Sun-SNK18/BT*
(for Sun-SNK18)

Sun-TAP15
Electronic thermostatic
mixer tap with time
flow sensor
(for Sun-SNK14 & 15)

Sun-TAP16
Electronic
thermostatic mixer tap
with proximity sensor
(for Sun-SNK14 & 15)


Sun-TAP17
Thermostatic
sequential mixer tap
(for Sun-SNK14 & 15)

Sun-TAP18
Washbasin mixer tap
(for Sun-SNK16,
17 & 18)

* Required for sink installation and for HTM64 compliance


dimensions

Sun-SNK14 & Sun-SNK15



	A	B	C	D	E	G	H	J
Sun-SNK14 (M)	500	400	270	425	360	187	53	7
Sun-SNK15 (L)	600	480	342	500	450	196	55	9

Sun-SNK16, Sun-SNK17 & Sun-SNK18



	A	B	C	D	E	G	H	M
Sun-SNK16 (S)	400	365	333	216	190-210	62	190	150
Sun-SNK17 (M)	500	420	420	265	215-265	62	195	150
Sun-SNK18 (L)	600	479	503	315	270-290	63	195	200

Sinks (HTM64 Compliant)

Washbasins

Our range of HTM64 compliant washbasins are suitable for use in clinical areas and for duct panels on integrated plumbing systems (IPS). Made from vitreous china to BS 3402.



Code: **Sun-SNK14**
Washbasin (rear draining) - medium



Code: **Sun-SNK15**
Washbasin (rear draining) - large



Our range of HTM64 compliant washbasins are suitable for personal washing (face, forearms and hands). Made from vitreous china to BS 3402.



Code: **Sun-SNK16**
Washbasin (bottom draining) - small



Code: **Sun-SNK17**
Washbasin (bottom draining) - medium



Code: **Sun-SNK18**
Washbasin (bottom draining) - large

product information - Sun-TAP13



Suitable for clinical use



200mm tap centres



Integrated thermostat regulates temperature



180° Span single lever operation



Built in thermostat cleansing feature (requires bridging hose)



Durable attractive chrome finish



Easy clean design and materials



1 year warranty



British made

product information - Sun-TAP15



Suitable for clinical use



Pre-determined water running length time (set up on installation)



Suitable for scrub up applications



200mm tap centres



Integral thermostat



Built in thermostat cleansing feature (requires bridging hose)



Durable attractive chrome finish



Easy clean design and materials



1 year warranty



British made

product information - Sun-TAP16



Suitable for clinical use



For use in general areas where water saving benefits are required



Activates and stops by movement and proximity



200mm tap centres



Integral thermostat



Built in thermostat cleansing feature (requires bridging hose)



Durable attractive chrome finish



Easy clean design and materials




1 year warranty




British made


Our range of HTM64 compliant taps are suitable for use in clinical areas and in duct panels on integrated plumbing systems (IPS).



Code: **Sun-TAP13**
Deck mounted sequential thermostatic mixer tap



Code: **Sun-TAP15**
Electronic thermostatic mixer tap with time flow sensor



Code: **Sun-TAP16**
Electronic thermostatic mixer tap with proximity sensor

product information - Sun-TAP17



Suitable for clinical use



For use in general areas



200mm tap centres



Integral thermostat



180° Span single lever operation



Built in thermostat cleansing feature (requires bridging hose)



Durable attractive chrome finish



Easy clean design and materials



1 year warranty



British made

product information - Sun-TAP18



Suitable for personal washing only



Sequential operation means water starts at cold, removing risk of scalding



Integral thermostat



Single lever operation providing ease of use for the less able



Suitable for use with Sun-SNK17 & 18 washbasins



Durable attractive chrome finish



Easy clean design and materials






1 year warranty





British made

Our range of HTM64 compliant taps are suitable for use in clinical areas and personal cleansing.



Code: **Sun-TAP17**
Thermostatic sequential mixer tap



Code: **Sun-TAP18**
Washbasin mixer tap

product information - Sun-SNK30



Suitable for personal washing only



200mm tap centres



Premium quality vitreous china



Counter/worktop inset design



No plug, chain or overflow



Easy clean design and materials



1 year warranty



British made

product information - Sun-TAP5



Suitable for personal washing only



200mm projection swivel spout



200mm tap centres



Durable attractive chrome finish



Traditional twin lever "elbow operation" design



Easy clean design and materials




1 year warranty




British made

Suitable for general healthcare and personal washing, our range of sinks and taps are designed to improve hygiene in healthcare environments.



Code: **Sun-SNK30**
Vitreous china inset washbasin




Code: **Sun-TAP5**
Traditional twin lever mixer tap

relaxing
contemporary
support

Visitor, Patient & Waiting Room Seating

Contemporary, durable and hygienic range of visitor and patient seating designed for healthcare environments. Available in a wide range of colours, patterns and fabrics, including anti-bacterial.

		PAGE
	 Neptune Visitor 3 Seat Modular Seating	163-164
	 Neptune Visitor 4 Seat Modular Seating	165-166
	 Neptune Visitor 5 Seat Modular Seating	167-168
	 Neptune Visitor Chairs	169-170
	Venus Visitor 3 Seat Modular Seating	171-172
	Venus Visitor 4 Seat Modular Seating	173-174
	Venus Visitor 5 Seat Modular Seating	175-176
	Visitor & Patient Chairs	177-182
	Atlas Seating	183-186
	Bariatric Chairs	187-190

product information

127kg

Maximum user weight (20st.) per seat

One piece durable moulded polypropylene seats

Comfortable foam upholstered seat pads

Strong ergonomic chrome arms

Stable ergonomic chrome legs

Adjustable feet for uneven floors

7

Seat Colours

Available in 7 moulded seat colours

4

Seat Pads

Choice of 4 upholstered seat pads

Vinyl or inter/vene Anti-bacterial upholstery

easy assembly

Minimal assembly required (tools supplied)

installation

Optional installation service available

5

5 year upholstery warranty

2

2 year warranty on vinyl, frame and moulded seats

Easy clean design and materials

FIRA

Independently tested by FIRA


Anti-bacterial Vinyl Strength Test
BS EN ISO 13934-1: 2013
BS EN ISO 13973-3: 2000
BS EN ISO 12947-2: 1998

Anti-bacterial Vinyl

Fire (Flammability) Testing
Flammability BS 7176 Medium Hazard
BS 5852 : 2006 Clause 11 (composites)


British made

dimensions




moulded seat and upholstery colours - please specify when ordering
please see inside back cover for full specification of our upholstery range

moulded seat colours



anti-bacterial vinyl seat pad upholstery colours



anti-bacterial inter/vene™ seat pad upholstery colours



Neptune offers a contemporary range of seating modules and chairs, with great design and stability whilst being durable, hygienic and easy to clean. Available in 7 modern and vibrant colours with the option of upholstered seat pads in 2 anti-bacterial materials, vinyl and inter/vene™.



Code: **Sun-SEAT73/3/BLACK** - Neptune visitor 3 seat module in black - 2 chrome arms



Code: **Sun-SEAT73/3/GREEN/VYL/BLK**
Neptune visitor 3 seat module in green with 3 black seat pads in anti-bacterial vinyl - 2 chrome arms



Code: **Sun-SEAT73/3/YELLOW/IV/GRY**
Neptune visitor 3 seat module in yellow with 3 grey seat pads in anti-bacterial inter/vene - 2 chrome arms

product information

Maximum user weight
(20st.) per seat

One piece
durable moulded
polypropylene seats

Comfortable foam
upholstered
seat pads

Strong ergonomic
chrome arms

Stable ergonomic
chrome legs

Adjustable feet for
uneven floors

7
Seat Colours

4
Seat Pads

Vinyl or inter/vene
Anti-bacterial upholstery

easy assembly
Minimal assembly
required
(tools supplied)

installation
Optional installation
service available

5
5 year upholstery
wearability warranty
(inter/vene only)

2
2 year warranty
on vinyl, frame and
moulded seats

eze
Easy clean design
and materials

FIRA
Independently
tested by FIRA
Anti-bacterial Vinyl Strength Test
BS EN ISO 13934-1: 2013
BS EN ISO 13973-3: 2000
BS EN ISO 12947-2: 1998

West Yorkshire
Materials
Testing Service
Anti-bacterial Vinyl

British made


Fire (Flammability) Testing
Flammability BS 7176 Medium Hazard
BS 5852 : 2006 Clause 11 (composites)

dimensions




moulded seat and upholstery colours - please specify when ordering
please see inside back cover for full specification of our upholstery range

moulded seat colours



anti-bacterial vinyl
seat pad upholstery colours



anti-bacterial inter/vene™
seat pad upholstery colours



Neptune offers a contemporary range of seating modules and chairs, with great design and stability whilst being durable, hygienic and easy to clean. Available in 7 modern and vibrant colours with the option of upholstered seat pads in 2 anti-bacterial materials, vinyl and inter/vene™.



Code: **Sun-SEAT74/4/ORANGE** - Neptune visitor 4 seat module in orange - 2 chrome arms



Code: **Sun-SEAT74/4/IVORY/VYL/GRY**
Neptune visitor 4 seat module in ivory with 4 grey seat pads in anti-bacterial vinyl - 2 chrome arms



Code: **Sun-SEAT74/4/BBLUE/IV/BLK**
Neptune visitor 4 seat module in blue with 4 black seat pads in anti-bacterial inter/vene - 2 chrome arms

product information

127kg

Maximum user weight (20st.) per seat

One piece durable moulded polypropylene seats

Comfortable foam upholstered seat pads

Strong ergonomic chrome arms

Stable ergonomic chrome legs

Adjustable feet for uneven floors

7 Seat Colours

Available in 7 moulded seat colours

4 Seat Pads

Choice of 4 upholstered seat pads

Vinyl or inter/vene Anti-bacterial upholstery

easy assembly

Minimal assembly required (tools supplied)

installation

Optional installation service available

5 Year Upholstery Warranty

5 year upholstery wearability warranty (inter/vene only)

2 Year Warranty

2 year warranty on vinyl, frame and moulded seats

eze

Easy clean design and materials

FIRA

Independently tested by FIRA

Anti-bacterial Vinyl Strength Test
BS EN ISO 13934-1: 2013
BS EN ISO 13973-3: 2000
BS EN ISO 12947-2: 1998

West Yorkshire Materials Testing Service

Anti-bacterial Vinyl

Fire (Flammability) Testing
Flammability BS 7176 Medium Hazard
BS 5852 : 2006 Clause 11 (composites)

British Made


Quality Assured

dimensions




moulded seat and upholstery colours - please specify when ordering
please see inside back cover for full specification of our upholstery range

moulded seat colours



anti-bacterial vinyl seat pad upholstery colours



anti-bacterial inter/vene™ seat pad upholstery colours



Neptune offers a contemporary range of seating modules and chairs, with great design and stability whilst being durable, hygienic and easy to clean. Available in 7 modern and vibrant colours with the option of upholstered seat pads in 2 anti-bacterial materials, vinyl and inter/vene™.



Code: Sun-SEAT75/5/YELLOW - Neptune visitor 5 seat module in yellow - 2 chrome arms



Code: Sun-SEAT75/5/GREY/VYL/BLK
Neptune visitor 5 seat module in grey with 5 black seat pads in anti-bacterial vinyl - 2 chrome arms



Code: Sun-SEAT75/5/ORANGE/IV/GRY
Neptune visitor 5 seat module in orange with 5 grey seat pads in anti-bacterial inter/vene - 2 chrome arms

product information

Maximum user weight (20st.) per seat

One piece durable moulded polypropylene seat

Comfortable foam upholstered seat pad

Comfortable and durable moulded arm rests (Sun-SEAT71 only)

Optional chair linking kit

7 Seat Colours

Available in 7 moulded seat colours

4 Seat Pads

Choice of 4 upholstered seat pads

Anti-bacterial vinyl and inter/vene upholstery

Supplied fully assembled

5 year upholstery wearability warranty (inter/vene only)

2 year warranty on vinyl, frame and moulded seats

Easy clean design and materials

Anti-bacterial Vinyl Strength Test BS EN ISO 13934-1: 2013 BS EN ISO 13973-3: 2000 BS EN ISO 12947-2: 1998

Fire (Flammability) Testing Flammability BS 7176 Medium Hazard BS 5852 : 2006 Clause 11 (composites)

British made

dimensions

Sun-SEAT70

78cm(H)
47cm(W)
47cm(D)
Seat Height 47cm

Sun-SEAT71

78cm(H)
61cm(W)
47cm(D)
Seat Height 47cm

accessories

Sun-SEAT7071/LINK
Neptune Visitor Chair Linking Kit
(2 sets per chair for front & back legs)

moulded seat and upholstery colours - please specify when ordering
please see inside back cover for full specification of our upholstery range

moulded seat colours

black

grey

blue

green

yellow

orange

ivory

anti-bacterial vinyl seat pad upholstery colours

black

grey

anti-bacterial inter/vene™ seat pad upholstery colours

black

grey

Neptune offers a contemporary range of chairs and seating modules, with great design and stability whilst being durable, hygienic and easy to clean. Available in 7 modern and vibrant colours with the option of upholstered seat pads in 2 anti-bacterial materials, vinyl and inter/vene™.



Code: **Sun-SEAT70/ORANGE/VYL/BLK**
Neptune visitor chair in orange with black vinyl upholstered seat pad



Code: **Sun-SEAT71/GREEN/VYL/GRY**
Neptune visitor chair with arms in green and grey vinyl upholstered seat pad



Sun-SEAT7071/LINK
Neptune Visitor Chair Linking Kit
(2 sets per chair for front & back legs)



Code: **Sun-SEAT71/BLUE/IV/GRY**
Neptune visitor chair with arms in blue and grey inter/vene upholstered seat pad



Code: **Sun-SEAT70/YELLOW/IV/BLK**
Neptune visitor chair in yellow with black inter/vene upholstered seat pad



product information

Maximum user weight
(17st.) per seat

Durable moulded
polypropylene seats

Comfortable deep
foam upholstered seats

Rounded/upholstered
seat corners

Vinyl or inter/vene
Anti-bacterial upholstery

Adjustable feet for
uneven floors

Durable silver epoxy
coated frame

Magazine table
option

easy assembly
Minimal assembly
required (tools supplied)

installation
Optional installation
service available

5 year upholstery
wearability warranty
(inter/vene only)

2 year warranty
on frame and
moulded seats

Easy clean design
and materials

FIRA
Independently
tested by FIRA

West Yorkshire
Materials
Testing Service
Anti-bacterial Vinyl

British made

dimensions

Sun-SEAT3
78cm(H) 153cm(W) 66cm(D)

Sun-SEAT3T
78cm(H) 153cm(W) 66cm(D)

Magazine Table
45cm(W) 45cm(D)

upholstery colours - please specify when ordering
please see inside back cover for full specification of our upholstery range

polypropylene seats

blackbluegreybeige

anti-bacterial vinyl upholstery colours

blackgreybeigeprimrosewalnutred winesalmonlavendernavymid bluesky bluecool bluegreenmintwhite

anti-bacterial inter/vene™ upholstery colours

blackgreybeigeyellowplumrednavysky bluegreenpastel green

Contemporary range of 3 seat seating modules, which are durable, hygienic and easy to clean. Available in a wide range of colours and anti-bacterial fabrics.



Code: **Sun-SEAT3** - Venus visitor 3 seat module with moulded polypropylene seats
Code: **Sun-SEAT3T** - Venus visitor 2 seat & 1 table module with moulded polypropylene seats



Code: **Sun-SEAT3VYL** - Venus visitor 3 seat module with anti-bacterial vinyl upholstery
Code: **Sun-SEAT3TVYL** - Venus visitor 2 seat & 1 table module with anti-bacterial vinyl upholstery



Code: **Sun-SEAT3IV** - Venus visitor 3 seat module with anti-bacterial inter/vene upholstery
Code: **Sun-SEAT3TIV** - Venus visitor 2 seat & 1 table module with anti-bacterial inter/vene upholstery

product information

Maximum user weight (17st.) per seat

Durable moulded polypropylene seats

Comfortable deep foam upholstered seats

Rounded/upholstered seat corners

Vinyl or inter/vene
Anti-bacterial upholstery

Adjustable feet for uneven floors

Durable silver epoxy coated frame

Magazine table option

easy assembly
Minimal assembly required (tools supplied)

installation
Optional installation service available

5 year upholstery wearability warranty (inter/vene only)

2 year warranty on frame and moulded seats

Easy clean design and materials

FIRA
Independently tested by FIRA

West Yorkshire Materials Testing Service
Anti-bacterial Vinyl

British made

Anti-bacterial Vinyl Strength Test
BS EN ISO 13934-1: 2013
BS EN ISO 13973-3: 2000
BS EN ISO 12947-2: 1998

Fire (Flammability) Testing
Flammability BS 7176 Medium Hazard
BS 5852 : 2006 Clause 11 (composites)

dimensions

Sun-SEAT4 78cm(H) 203cm(W) 66cm(D)

Sun-SEAT4T 78cm(H) 203cm(W) 66cm(D)

Magazine Table
45cm(W) 45cm(D)

upholstery colours - please specify when ordering
please see inside back cover for full specification of our upholstery range

polypropylene seats

black

blue

grey

beige

anti-bacterial vinyl upholstery colours

black

grey

beige

primrose

walnut

red wine

salmon

lavender

navy

mid blue

sky blue

cool blue

green

mint

white

anti-bacterial inter/vene™ upholstery colours

black

grey

beige

yellow

plum

red

navy

sky blue

green

pastel green

Contemporary range of 4 seat seating modules, which are durable, hygienic and easy to clean. Available in a wide range of colours and anti-bacterial fabrics.

Code: **Sun-SEAT4** - Venus visitor 4 seat module with moulded polypropylene seats
Code: **Sun-SEAT4T** - Venus visitor 3 seat & 1 table module with moulded polypropylene seats

FSC
www.fsc.org
FSC® C122891
The mark of responsible forestry

Code: **Sun-SEAT4VYL** - Venus visitor 4 seat module with anti-bacterial vinyl upholstery
Code: **Sun-SEAT4TVYL** - Venus visitor 3 seat & 1 table module with anti-bacterial vinyl upholstery

Code: **Sun-SEAT4IV** - Venus visitor 4 seat module with anti-bacterial inter/vene upholstery
Code: **Sun-SEAT4TIV** - Venus visitor 3 seat & 1 table module with anti-bacterial inter/vene upholstery

product information

Maximum user weight
(17st.) per seat

Durable moulded
polypropylene seats

Comfortable deep
foam upholstered seats

Rounded/upholstered
seat corners

Vinyl or inter/vene
Anti-bacterial upholstery

Adjustable feet for
uneven floors

Durable silver epoxy
coated frame

Magazine table
option

easy assembly

installation

5 year upholstery
wearability warranty
(inter/vene only)

2 year warranty
on frame and
moulded seats

Easy clean design
and materials

FIRA
Independently
tested by FIRA

West Yorkshire
Materials
Testing Service


Anti-bacterial Vinyl

British made

Anti-bacterial Vinyl Strength Test
BS EN ISO 13934-1: 2013
BS EN ISO 13973-3: 2000
BS EN ISO 12947-2: 1998

Fire (Flammability) Testing
Flammability BS 7176 Medium Hazard
BS 5852 : 2006 Clause 11
(composites)

dimensions




upholstery colours - please specify when ordering
please see inside back cover for full specification of our upholstery range

polypropylene seats



anti-bacterial vinyl upholstery colours



anti-bacterial inter/vene™ upholstery colours



Contemporary range of 5 seat seating modules, which are durable, hygienic and easy to clean. Available in a wide range of colours and anti-bacterial fabrics.



Code: **Sun-SEAT5** - Venus visitor 5 seat module with moulded polypropylene seats
Code: **Sun-SEAT5T** - Venus visitor 4 seat & 1 table module with moulded polypropylene seats



Code: **Sun-SEAT5VYL** - Venus visitor 5 seat module with anti-bacterial vinyl upholstery
Code: **Sun-SEAT5TVYL** - Venus visitor 4 seat & 1 table module with anti-bacterial vinyl upholstery



Code: **Sun-SEAT5IV** - Venus visitor 5 seat module with anti-bacterial inter/vene upholstery
Code: **Sun-SEAT5TIV** - Venus visitor 4 seat & 1 table module with anti-bacterial inter/vene upholstery

product information - venus visitor chairs

Maximum user weight (17st.)

Durable moulded polypropylene seats

Plastic feet prevent movement/protect floors

Safe and easy to stack (max. 10 high)

Supplied fully assembled

1 year warranty on frame & moulded seats

Easy clean design and materials

British made

Delivery normally within 3-5 working days *

Chair linking kits available (please see below)

* Please contact us regarding availability for quantities of 10 or more.

product information - galaxy visitor chairs

Maximum user weight (17st.)

Comfortable deep foam upholstery

Rounded/upholstered seat corners

Safe and easy to stack (max. 5 high)

Plastic feet prevent movement/protect floors

Vinyl or inter/vene Anti-bacterial upholstery

5 year upholstery wearability warranty (inter/vene only)

1 year warranty on frame

Supplied fully assembled

Supplied fully assembled

BS EN 16139: 2016 Furniture Strength, Durability and Safety Requirements

Anti-bacterial Vinyl Strength Test BS EN ISO 13934-1: 2013 BS EN ISO 13973-3: 2000 BS EN ISO 12947-2: 1998

Fire (Flammability) Testing BS 5852 : 2006 Clause 11 (composites)

British made

Delivery normally within 3-5 working days *

* Please contact us regarding availability for quantities of 10 or more.

upholstery colours - please specify when ordering

please see inside back cover for full specification of our upholstery range

anti-bacterial vinyl upholstery colours

anti-bacterial inter/vene™ upholstery colours

Contemporary range of chairs, which are durable, hygienic and easy to clean, designed for waiting, reception and public areas. Available in a wide range of colours and anti-bacterial fabrics.

Venus Visitor Chairs

Multi-purpose economy chairs, ergonomically designed moulded seats and backs utilise modern materials, combining smooth surfaces, rounded corners with comfort, durability and ease of cleaning.

Code: **Sun-SEAT1** - Venus visitor chair with moulded polypropylene seat

Code: **Sun-SEAT2** - Venus visitor chair and arms with moulded polypropylene seat

Sun-SEAT/LINK1 Chair Linking Kit for Chairs with No Arms

Sun-SEAT/LINK2 Chair Linking Kit for Chairs with Arms

save when buying quantities of 4+ see price list for details

Sun-SEAT1

79cm(H) 54cm(W) 56cm(D)
Seat Height 44cm(H)

moulded polypropylene

black blue grey beige

Sun-SEAT2

79cm(H) 58cm(W) 56cm(D)
Seat Height 44cm(H)

moulded polypropylene

black blue grey beige

Galaxy Visitor Chairs

Multi-purpose upholstered economy chairs with comfortable deep foam ergonomically designed seat & back.

Code: **Sun-SEAT32VYL** - Galaxy visitor chair in anti-bacterial vinyl upholstery

Code: **Sun-SEAT33VYL** - Galaxy visitor chair with arms in anti-bacterial vinyl upholstery

Code: **Sun-SEAT32IV** - Galaxy visitor chair in anti-bacterial inter/vene upholstery

Code: **Sun-SEAT33IV** - Galaxy visitor chair with arms in anti-bacterial inter/vene upholstery

Sun-SEAT/LINK1 Chair Linking Kit for Chairs with No Arms

Sun-SEAT/LINK2 Chair Linking Kit for Chairs with Arms

save when buying quantities of 3+ see price list for details

product information - aurora visitor chairs



110kg


Maximum user weight (17st.)



Comfortable deep foam upholstery



Vinyl or inter/vene Anti-bacterial upholstery



Solid lacquered beech wood frame

3

stackable

Safe and easy to stack (max. 3 high)

product information - aurora easy access chairs



110kg

Maximum user weight (17st.)



Comfortable deep foam upholstery



Vinyl or inter/vene Anti-bacterial upholstery



Solid lacquered beech frame

3

stackable

Safe and easy to stack (max. 3 high)

product information - aurora bariatric chair



222kg

Maximum user weight (34st.)



Comfortable deep foam upholstery



Solid beech arms



Solid beech arms



Ergonomic design

additional product information



5 year upholstery wearability warranty (inter/vene only)



2 year warranty on frame



Supplied fully assembled



Supplied fully assembled



BS EN 16139: 2016 Furniture Strength, Durability and Safety Requirements (Sun-SEAT35 & Sun-SEAT36 only)



Anti-bacterial Vinyl Strength Test
BS EN ISO 13934-1: 2013
BS EN ISO 13973-3: 2000
BS EN ISO 12947-2: 1998



Fire (Flammability) Testing
Flammability BS 7176 Medium Hazard
BS 5852 : 2006 Clause 11 (composites)



British made



Delivery normally within 3-5 working days */**


* Please contact us regarding availability for quantities of 10 or more.

** Sun-SEAT35, Sun-SEAT36 & Sun-SEAT37 only


upholstery colours - please specify when ordering

please see inside back cover for full specification of our upholstery range


anti-bacterial vinyl upholstery colours




black




grey




beige




primrose




walnut




red wine




salmon




lilac




navy




mid blue




sky blue




cool blue



green




mint




white


anti-bacterial inter/vene™ upholstery colours




black




grey




beige




yellow




plum




red




navy



sky blue



green



pastel green


Aurora Visitor Chairs

Deluxe beech frame upholstered chairs with comfortable deep foam padded seat & back, with or without arms.




- Code: **Sun-SEAT35VYL** - Aurora visitor chair with anti-bacterial vinyl upholstery
- Code: **Sun-SEAT36VYL** - Aurora visitor chair and arms with anti-bacterial vinyl upholstery
- Code: **Sun-SEAT35IV** - Aurora visitor chair with anti-bacterial inter/vene upholstery
- Code: **Sun-SEAT36IV** - Aurora visitor chair and arms with anti-bacterial inter/vene upholstery

Sun-SEAT35



84cm(H)
51cm(W)
57cm(D)
Seat Height 45cm

Sun-SEAT36



84cm(H)
53cm(W)
57cm(D)
Seat Height 45cm

save when buying quantities of 3+

see price list for details

Aurora Easy Access Chairs

The Aurora Easy Access Chair is a range extension to the existing and successful Aurora seating. With an extended squared arm rest to aid in the entry and exit from the chair seat. This is ideal for all users but particularly of use to older or less able visitors.



Sun-SEAT37



84cm(H)
56cm(W)
57cm(D)
Seat Height 48cm

save when buying quantities of 3+

see price list for details

- Code: **Sun-SEAT37VYL** - Aurora easy access chair with arms in anti-bacterial vinyl upholstery
- Code: **Sun-SEAT37IV** - Aurora easy access chair with arms in anti-bacterial inter/vene upholstery

Aurora Bariatric Chair

Our Aurora bariatric chair combines strength, stability and comfort with ease of access and is easy to clean.



Sun-CHA75



88cm(H)
84cm(W)
58cm(D)
Seat Width 72cm
Seat Height 45cm

- Code: **Sun-CHA75VYL** - Aurora Bariatric 222kg (34st.) arm chair with anti-bacterial vinyl upholstery
- Code: **Sun-CHA75IV** - Aurora Bariatric 222kg (34st.) arm chair with anti-bacterial inter/vene upholstery

product information

Maximum user weight (17st.)

Comfortable deep foam upholstery

Vinyl or inter/vene Anti-bacterial upholstery

Safe and easy to stack (max. 3 high)

5 year upholstery wearability warranty (inter/vene only)

2 year warranty on frame

Supplied fully assembled

Easy clean design and materials

BS EN 16139: 2016 Furniture Strength, Durability and Safety Requirements

Anti-bacterial Vinyl Strength Test
BS EN ISO 13934-1: 2013
BS EN ISO 13973-3: 2000
BS EN ISO 12947-2: 1998

Fire (Flammability) Testing
Flammability BS 7176 Medium Hazard
BS 5852: 2006 Clause 11 (composites)

British made

Delivery normally within 3-5 working days *

* Please contact us regarding availability for quantities of 10 or more.

upholstery colours - please specify when ordering
please see inside back cover for full specification of our upholstery range

anti-bacterial vinyl upholstery colours

black	grey	beige	primrose	walnut	red wine	salmon	lilac	navy	mid blue	sky blue	cool blue	green	mint	white
-------	------	-------	----------	--------	----------	--------	-------	------	----------	----------	-----------	-------	------	-------

anti-bacterial inter/vene™ upholstery colours

black	grey	beige	yellow	plum	red	navy	sky blue	green	pastel green
-------	------	-------	--------	------	-----	------	----------	-------	--------------

Genesis Chairs

The Genesis Easy Access Chair provides an extended squared arm rest to aid in the entry and exit from the chair seat. This is ideal for all users but particularly of use to older or less able visitors. This easy access feature is combined with an attractive, modern, slightly deeper, square section seat back. The chair is available as a side chair to compliment the arm chair.



Code: **Sun-SEAT41VYL** - Genesis easy access chair in anti-bacterial vinyl upholstery

Code: **Sun-SEAT41IV** - Genesis easy access chair in anti-bacterial inter/vene upholstery

Sun-SEAT41

86cm(H)
56cm(W)
57cm(D)
Seat Height 48cm

save when buying quantities of 3+ see price list for details



Code: **Sun-SEAT42VYL** - Genesis easy access chair with arms in anti-bacterial vinyl upholstery

Code: **Sun-SEAT42IV** - Genesis easy access chair with arms in anti-bacterial inter/vene upholstery

Sun-SEAT42

86cm(H)
56cm(W)
57cm(D)
Seat Height 48cm

save when buying quantities of 3+ see price list for details

product information

Maximum user weight (17st.)

Comfortable deep foam upholstery

Solid beech frame

Solid beech arms

Removable seat pad for cleaning

Optional pressure relieving foam upgrade

Optional drop-arms to aid lateral transfer (mid & high back chairs)

Wings for head support (Sun-CHA52 only)

Optional adjustable height legs (10cm extension)

Optional house keeping wheels (fitted to rear legs only)

Vinyl or inter/vene Anti-bacterial upholstery

5 year upholstery wearability warranty (inter/vene only)

2 year warranty

Supplied fully assembled

Easy clean design and materials

BS EN 16139: 2016 Furniture Strength, Durability and Safety Requirements

Anti-bacterial Vinyl Strength Test BS EN ISO 13934-1: 2013 BS EN ISO 13973-3: 2000 BS EN ISO 12947-2: 1998

Fire (Flammability) Testing Flammability BS 7176 Medium Hazard BS 5852 : 2006 Clause 11 (composites)

British made

accessories * these accessories can only be fitted when ordered with the chair

Sun-CHA48/PRF Pressure Relieving Foam Upgrade for Seat & Back

Sun-CHA49/PRF Pressure Relieving Foam Upgrade for Seat & Back

Sun-CHA50/PRF Pressure Relieving Foam Upgrade for Seat & Back

Sun-CHA50/AHL* Adjustable Height Legs (10cm extension)

Sun-CHA50/HKW* House Keeping Wheels (fitted to rear legs only)

upholstery colours - please specify when ordering
please see inside back cover for full specification of our upholstery range

anti-bacterial vinyl upholstery colours

black	grey	beige	primrose	walnut	red wine	salmon	lilac	navy	mid blue	sky blue	cool blue	green	mint	white
-------	------	-------	----------	--------	----------	--------	-------	------	----------	----------	-----------	-------	------	-------

anti-bacterial inter/vene™ upholstery colours

black	grey	beige	yellow	plum	red	navy	sky blue	green	pastel green
-------	------	-------	--------	------	-----	------	----------	-------	--------------

Atlas Patient Arm Chair Range
Deluxe solid beech frame upholstered seating range with deep foam padded seat and back for optimum comfort.



Code: **Sun-CHA48VYL** - Low back side chair with anti-bacterial vinyl upholstery
Code: **Sun-CHA48IV** - Low back side chair with anti-bacterial inter/vene upholstery



Code: **Sun-CHA49VYL** - Low back arm chair with anti-bacterial vinyl upholstery
Code: **Sun-CHA49IV** - Low back arm chair with anti-bacterial inter/vene upholstery



Code: **Sun-CHA50VYL** - Mid back arm chair with anti-bacterial vinyl upholstery
Code: **Sun-CHA50IV** - Mid back arm chair with anti-bacterial inter/vene upholstery



Code: **Sun-CHA50DAVYL** - Mid back arm chair with drop arms in anti-bacterial vinyl upholstery
Code: **Sun-CHA50DAIV** - Mid back arm chair with drop arms in anti-bacterial inter/vene upholstery

save when buying quantities of 2+ see price list for details

dimensions

Sun-CHA48	Sun-CHA49	Sun-CHA50	Sun-CHA50DA
<p>94cm(H) 54cm(W) 63cm(D)</p>	<p>94cm(H) 60cm(W) 60cm(D)</p>	<p>104cm(H) 60cm(W) 60cm(D)</p>	<p>104cm(H) 60cm(W) 60cm(D)</p>
Seat Heights 48cm - Internal Seat Widths 49cm			

product information

Maximum user weight (17st.)

Comfortable deep foam upholstery

Solid beech frame

Solid beech arms

Removable seat pad for cleaning

Optional pressure relieving foam upgrade

Optional drop-arms to aid lateral transfer (mid & high back chairs)

Wings for head support (Sun-CHA52 only)

Optional adjustable height legs (10cm extension)

Optional house keeping wheels (fitted to rear legs only)

Vinyl or inter/vene Anti-bacterial upholstery

5 year upholstery wearability warranty (inter/vene only)

2 year warranty

Supplied fully assembled

Easy clean design and materials

BS EN 16139: 2016 Furniture Strength, Durability and Safety Requirements

Anti-bacterial Vinyl Strength Test BS EN ISO 13934-1: 2013 BS EN ISO 13973-3: 2000 BS EN ISO 12947-2: 1998

Fire (Flammability) Testing BS 5852 : 2006 Clause 11 (composites)

British made

accessories * these accessories can only be fitted when ordered with the chair

Sun-CHA51/PRF Pressure Relieving Foam Upgrade for Seat & Back

Sun-CHA52/PRF Pressure Relieving Foam Upgrade for Seat & Back

Sun-CHA50/AHL* Adjustable Height Legs (10cm extension)

Sun-CHA50/HKW* House Keeping Wheels (fitted to rear legs only)

upholstery colours - please specify when ordering
please see inside back cover for full specification of our upholstery range

anti-bacterial vinyl upholstery colours

black	grey	beige	primrose	walnut	red wine	salmon	lilac	navy	mid blue	sky blue	cool blue	green	mint	white
-------	------	-------	----------	--------	----------	--------	-------	------	----------	----------	-----------	-------	------	-------

anti-bacterial inter/vene™ upholstery colours

black	grey	beige	yellow	plum	red	navy	sky blue	green	pastel green
-------	------	-------	--------	------	-----	------	----------	-------	--------------

Atlas Patient Arm Chair Range
Deluxe solid beech frame upholstered seating range with deep foam padded seat and back for optimum comfort.



Code: **Sun-CHA51VYL** - High back arm chair with anti-bacterial vinyl upholstery

Code: **Sun-CHA51IV** - High back arm chair with anti-bacterial inter/vene upholstery

Code: **Sun-CHA51DAVYL** - High back arm chair with drop arms in anti-bacterial vinyl upholstery

Code: **Sun-CHA51DAIV** - High back arm chair with drop arms in anti-bacterial inter/vene upholstery



Code: **Sun-CHA52VYL** - High back arm chair with wings in anti-bacterial vinyl upholstery

Code: **Sun-CHA52IV** - High back arm chair with wings in anti-bacterial inter/vene upholstery

Code: **Sun-CHA52DAVYL** - High back arm chair with wings and drop arms in anti-bacterial vinyl upholstery

Code: **Sun-CHA52DAIV** - High back arm chair with wings and drop arms in anti-bacterial inter/vene upholstery

dimensions

Sun-CHA51	Sun-CHA51DA	Sun-CHA52	Sun-CHA52DA
Seat Heights 48cm - Internal Seat Widths 49cm			

product information

254kg

Maximum user weight (40st.)

Comfortable deep foam upholstery

Solid beech frame

Solid beech arms

Removable seat pad for cleaning

Optional pressure relieving foam upgrade

Optional drop-arms to aid lateral transfer (mid & high back chairs)

Wings for head support (Sun-CHA56 only)

Optional adjustable height legs (10cm extension)

Optional house keeping wheels (fitted to rear legs only)

Vinyl or inter/vene Anti-bacterial upholstery

5 year upholstery wearability warranty (inter/vene only)

2 year warranty

Supplied fully assembled

Easy clean design and materials

BS EN 16139: 2016 Furniture Strength, Durability and Safety Requirements

Anti-bacterial Vinyl Strength Test BS EN ISO 13934-1: 2013 BS EN ISO 13973-3: 2000 BS EN ISO 12947-2: 1998

Fire (Flammability) Testing BS 5852 : 2006 Clause 11 (composites)

British made

accessories * these accessories can only be fitted when ordered with the chair

Sun-CHA55/PRF Pressure Relieving Foam Upgrade for Seat & Back

Sun-CHA56/PRF Pressure Relieving Foam Upgrade for Seat & Back

Sun-CHA50/AHL* Adjustable Height Legs (10cm extension)

Sun-CHA50/HKW* House Keeping Wheels (fitted to rear legs only)

upholstery colours - please specify when ordering
please see inside back cover for full specification of our upholstery range

anti-bacterial vinyl upholstery colours

black	grey	beige	primrose	walnut	red wine	salmon	lilac	navy	mid blue	sky blue	cool blue	green	mint	white
-------	------	-------	----------	--------	----------	--------	-------	------	----------	----------	-----------	-------	------	-------

anti-bacterial inter/vene™ upholstery colours

black	grey	beige	yellow	plum	red	navy	sky blue	green	pastel green
-------	------	-------	--------	------	-----	------	----------	-------	--------------

Bariatric Chairs

Atlas Bariatric High Back Arm Chair Range
Deluxe solid beech frame upholstered seating range with deep foam padded seat and back for optimum comfort.



Code: **Sun-CHA55VYL** - High back bariatric arm chair with anti-bacterial vinyl upholstery
Code: **Sun-CHA55IV** - High back bariatric arm chair with anti-bacterial inter/vene upholstery
Code: **Sun-CHA55DAVYL** - High back bariatric arm chair with drop arms in anti-bacterial vinyl upholstery
Code: **Sun-CHA55DAIV** - High back bariatric arm chair with drop arms in anti-bacterial inter/vene upholstery




Code: **Sun-CHA56VYL** - High back bariatric arm chair with wings in anti-bacterial vinyl upholstery
Code: **Sun-CHA56IV** - High back bariatric arm chair with wings in anti-bacterial inter/vene upholstery
Code: **Sun-CHA56DAVYL** - High back bariatric arm chair with wings and drop arms in anti-bacterial vinyl upholstery
Code: **Sun-CHA56DAIV** - High back bariatric arm chair with wings and drop arms in anti-bacterial inter/vene upholstery

dimensions


Sun-CHA55	Sun-CHA55DA	Sun-CHA56	Sun-CHA56DA
110cm(H) 93cm(W) 73cm(D)	110cm(H) 93cm(W) 73cm(D)	126cm(H) 93cm(W) 73cm(D)	126cm(H) 93cm(W) 73cm(D)

Seat Heights 51cm - Internal Seat Widths 81cm


product information - aurora bariatric chair (Sun-CHA75)




Maximum user weight (34st.)




Comfortable deep foam upholstery




Solid beech arms




Solid beech arms




Ergonomic design




Vinyl or inter/vene Anti-bacterial upholstery




5 year upholstery wearability warranty (inter/vene only)




2 year warranty




Supplied fully assembled




Easy clean design and materials



Anti-bacterial Vinyl Strength Test
BS EN ISO 13934-1: 2013
BS EN ISO 13973-3: 2000
BS EN ISO 12947-2: 1998



Fire (Flammability) Testing
Flammability BS 7176 Medium Hazard
BS 5852 : 2006 Clause 11 (composites)



British made

Aurora Bariatric Chair
Our Aurora bariatric chair combines strength, stability and comfort with ease of access and is easy to clean.



Code: **Sun-CHA75VYL** - Aurora Bariatric 222kg (34st.) arm chair with anti-bacterial vinyl upholstery
Code: **Sun-CHA75IV** - Aurora Bariatric 222kg (34st.) arm chair with anti-bacterial inter/vene upholstery



88cm(H)
84cm(W)
58cm(D)
Seat Width 72cm
Seat Height 45cm

upholstery colours - please specify when ordering
please see inside back cover for full specification of our upholstery range

anti-bacterial vinyl upholstery colours



black grey beige primrose walnut red wine salmon lilac navy mid blue sky blue cool blue green mint white

anti-bacterial inter/vene™ upholstery colours








black grey beige yellow plum red navy sky blue green pastel green

functional
durable
workspace

Desks, Chairs, Filing, Bookcases & Administration Storage

Durable, contemporary and stylish, our premium range of administration furniture, chairs and storage solutions are ideal for today's busy consultation environments. Available in a choice of finishes, all of which are designed to be easy to clean.

	PAGE
 Desks	192-197
 Desk Pedestals & Storage	198-199
 Tambour Units & Storage	200-201
 Filing Cabinets, Bookcases & Storage	202-203
 Consultation Chairs	204-207

product information

6
widths

Available in six widths
(80 - 180cm)

ten
finishes

Choice of 10 attractive
laminate finishes

Strong & durable
18mm thick desk top
(also available in 25mm)

Stylish silver legs &
modesty back panel

Integrated cable
management system

minimal

Minimal assembly
required

eze

Easy clean design
and materials

5
Year Warranty

5 year warranty

British Made
Quality Assured

British made

accessories - rectangular desks

Sun-FPED1
Suspended under
desk 2 drawer unit
(80cm(D) rectangular
desks only)

Sun-FPED2
Suspended under
desk 3 drawer unit
(80cm(D) rectangular
desks only)

Sun-UPED1
2 drawer under
desk pedestal
(80cm(D) rectangular
desks only)

Sun-UPED2
3 drawer under
desk pedestal
(80cm(D) rectangular
desks only)

Sun-UPED3
4 drawer under
desk pedestal
(80cm(D) rectangular
desks only)

Sun-HPED8
80cm(D) 2 drawer
desk height pedestal

Sun-HPED9
80cm(D) 3 drawer
desk height pedestal

Sun-HPED10
80cm(D) 4 drawer
desk height pedestal

Sun-DESKWS
Wave desktop
screens
(please enquire)

Sun-DESKRS
Rectangular
desktop screens
(please enquire)

accessories - rectangular (narrow) desks

Sun-HPED5
60cm(D) 2 drawer
desk height pedestal

Sun-HPED6
60cm(D) 3 drawer
desk height pedestal

Sun-HPED7
60cm(D) 4 drawer
desk height pedestal

Sun-DESKWS
Wave desktop
screens
(please enquire)

Sun-DESKRS
Rectangular
desktop screens
(please enquire)

finishes - please specify when ordering

beech

maple

cherry

japanese
ash

light oak

modern
oak

walnut

white

silver

plain
grey

Premium range of desks which are stylish & contemporary in design, strong and durable in construction, suitable for consultation and administration environments. Available in a range of 10 attractive laminate finishes and guaranteed for 5 years.

Rectangular Desks - 74cm(H) x 80cm(D)

dimensions
| 80cm - 180cm(W) |
| 80cm(D) |

Code: **Sun-DESK1/80**
80cm wide rectangular desk

Code: **Sun-DESK1/100**
100cm wide rectangular desk

Code: **Sun-DESK1/120**
120cm wide rectangular desk

Code: **Sun-DESK1/140**
140cm wide rectangular desk

Code: **Sun-DESK1/160**
160cm wide rectangular desk

Code: **Sun-DESK1/180**
180cm wide rectangular desk

Rectangular (Narrow) Desks - 74cm(H) x 60cm(D)

dimensions
| 80cm - 180cm(W) |
| 60cm(D) |

Code: **Sun-DESK2/80**
80cm wide rectangular
(narrow) desk

Code: **Sun-DESK2/100**
100cm wide rectangular
(narrow) desk

Code: **Sun-DESK2/120**
120cm wide rectangular
(narrow) desk

Code: **Sun-DESK2/140**
140cm wide rectangular
(narrow) desk

Code: **Sun-DESK2/160**
160cm wide rectangular
(narrow) desk

Code: **Sun-DESK2/180**
180cm wide rectangular
(narrow) desk

product information



Available in three widths
(140 - 180cm)



Choice of 10 attractive
laminate finishes




Strong & durable
18mm thick desk top
(also available in 25mm)



Stylish silver legs &
modesty back panel



Integrated cable
management system



Available in left
or right hand



Minimal assembly
required



Easy clean design
and materials



5 year warranty



British made

accessories



Sun-FPED1
Suspended under
desk 2 drawer unit



Sun-FPED2
Suspended under
desk 3 drawer unit



Sun-UPED1
2 drawer under
desk pedestal



Sun-UPED2
3 drawer under
desk pedestal



Sun-UPED3
4 drawer under
desk pedestal



Sun-HPED5
60cm(D) 2 drawer
desk height pedestal



Sun-HPED6
60cm(D) 3 drawer
desk height pedestal



Sun-HPED7
60cm(D) 4 drawer
desk height pedestal




Sun-HPED8
80cm(D) 2 drawer
desk height pedestal




Sun-HPED9
80cm(D) 3 drawer
desk height pedestal



Sun-HPED10
80cm(D) 4 drawer
desk height pedestal



Sun-DESKWS
Wave desktop
screens
(please enquire)



Sun-DESKRS
Rectangular
desktop screens
(please enquire)

finishes - please specify when ordering



Premium range of desks which are stylish & contemporary in design, strong and durable in construction, suitable for consultation and administration environments. Available in a range of 10 attractive laminate finishes and guaranteed for 5 years.

J Shaped Desks (left hand) - 74cm(H) x 120cm(D)



dimensions

| 140cm - 180cm(W) |



Code: **Sun-DESK5/140**
140cm wide J shaped desk
(left hand)



Code: **Sun-DESK5/160**
160cm wide J shaped desk
(left hand)




Code: **Sun-DESK5/180**
180cm wide J shaped desk
(left hand)

J Shaped Desks (right hand) - 74cm(H) x 120cm(D)



dimensions

| 140cm - 180cm(W) |



Code: **Sun-DESK6/140**
140cm wide J shaped desk
(right hand)



Code: **Sun-DESK6/160**
160cm wide J shaped desk
(right hand)



Code: **Sun-DESK6/180**
180cm wide J shaped desk
(right hand)

product information - Wave & Symmetrical Desks

4

widths

Wave Desks
Available in four widths
(120 - 180cm)

2

widths

Symmetrical Desks
Available in two widths
(140 & 160cm)

ten

finishes

Choice of 10 attractive
laminate finishes

Strong & durable
18mm thick desk top
(also available in 25mm)

Stylish silver legs &
modesty back panel

Integrated cable
management system

Available in left
or right hand

minimal

Minimal assembly
required

Easy clean design
and materials

5 year warranty

British made

accessories - Wave & Symmetrical Desks

Sun-FPED1
Suspended under
desk 2 drawer unit
(Wave desks only)

Sun-FPED2
Suspended under
desk 3 drawer unit
(Wave desks only)

Sun-UPED1
2 drawer under
desk pedestal
(Wave desks only)

Sun-UPED2
3 drawer under
desk pedestal
(Wave desks only)

Sun-UPED3
4 drawer under
desk pedestal
(Wave desks only)

Sun-HPED5
60cm(D) 2 drawer
desk height pedestal

Sun-HPED6
60cm(D) 3 drawer
desk height pedestal

Sun-HPED7
60cm(D) 4 drawer
desk height pedestal

Sun-HPED8
80cm(D) 2 drawer
desk height pedestal

Sun-HPED9
80cm(D) 3 drawer
desk height pedestal

Sun-HPED10
80cm(D) 4 drawer
desk height pedestal

Sun-DESKWS
Wave desktop
screens
(please enquire)

Sun-DESKRS
Rectangular
desktop screens
(please enquire)

finishes - please specify when ordering

beech

maple

cherry

japanese
ash

light oak

modern
oak

walnut

white

silver

plain
grey

Premium range of desks which are stylish & contemporary in design, strong and durable in construction, suitable for consultation and administration environments. Available in a range of 10 attractive laminate finishes and guaranteed for 5 years.

Wave Desks (left or right hand) - 74cm(H) x 100cm(D)

dimensions (left Hand)
| 120cm - 180cm(W) |
| 100cm(D) |
| 80cm(D) |

dimensions (right Hand)
| 120cm - 180cm(W) |
| 80cm(D) |
| 100cm(D) |

Code: **Sun-DESK7/120** (left hand)
Code: **Sun-DESK8/120** (right hand)
120cm wide wave desk

Code: **Sun-DESK7/140** (left hand)
Code: **Sun-DESK8/140** (right hand)
140cm wide wave desk

Code: **Sun-DESK7/160** (left hand)
Code: **Sun-DESK8/160** (right hand)
160cm wide wave desk

Code: **Sun-DESK7/180** (left hand)
Code: **Sun-DESK8/180** (right hand)
180cm wide wave desk

Symmetrical Desks - 74cm(H)

dimensions
| 140cm(W) |
| 140cm(W) |
| 80cm(D) |
| 80cm(D) |

Code: **Sun-DESK9/140**
140cm wide symmetrical desk

160cm(W)
| 160cm(W) |
| 80cm(D) |
| 80cm(D) |

Code: **Sun-DESK9/160**
160cm wide symmetrical desk

product information - fixed under desk pedestals



Suspends underneath desk

2 & 3
drawer
options



Single lock, locks all drawers (2 keys supplies)

ten
finishes

Choice of 10 finishes

product information - under desk pedestals



Fits underneath desk

2, 3 & 4
drawer
options



Single lock, locks all drawers (2 keys supplies)

ten
finishes

Choice of 10 finishes

product information - desk height pedestals



Sits along side desk



Two depths 60cm & 80cm

2, 3 & 4
drawer
options



Single lock, locks all drawers (2 keys supplies)

ten
finishes

Choice of 10 finishes

additional product information



5 year warranty



5 year warranty



Supplied fully assembled



British made

dimensions

Suspended under desk drawer units



51cm(H)
43cm(W)
58cm(D)

Under desk pedestals



69cm(H)
44cm(W)
55cm(D)

Desk height pedestals



74cm(H)
44cm(W)
60cm(D)

Desk height pedestals



74cm(H)
44cm(W)
80cm(D)

finishes - please specify when ordering



beech maple cherry japanese ash light oak modern oak walnut white silver plain grey

Premium range of storage and filing solutions which are stylish & contemporary in design, strong and durable in construction. Available in a range of 10 attractive laminate finishes and guaranteed for 5 years.

Under desk suspended pedestals - 80cm(D)



Code: **Sun-FPED1**
2 Drawer fixed under desk pedestal



Code: **Sun-FPED2**
3 Drawer fixed under desk pedestal



Under desk pedestals - 80cm(D)



Code: **Sun-UPED1**
2 Drawer under desk pedestal



Code: **Sun-UPED2**
3 Drawer under desk pedestal



Code: **Sun-UPED3**
4 Drawer under desk pedestal

Desk height pedestals - 60cm(D)



Code: **Sun-HPED5**
2 Drawer desk height pedestal



Code: **Sun-HPED6**
3 Drawer desk height pedestal



Code: **Sun-HPED7**
4 Drawer desk height pedestal

Desk height pedestals - 80cm(D)



Code: **Sun-HPED8**
2 Drawer desk height pedestal



Code: **Sun-HPED9**
3 Drawer desk height pedestal



Code: **Sun-HPED10**
4 Drawer desk height pedestal

product information - tambour units



Adjustable height shelves



Locking roller doors (2 keys supplied)



Choice of 10 finishes

product information - desk height side filing cabinet



Holds files and folders



Two locking drawers (2 keys supplied)



Anti-tilt castor feature



Sits along side desk



Choice of 10 finishes

product information - desk height base unit



One adjustable height shelf



Locking door (2 keys supplied)



Choice of 10 finishes

additional product information



5 year warranty



5 year warranty




Supplied fully assembled



British made


dimensions

Tall Tambour Units




110cm(W) 53cm(D)

Desk height tambour unit




74cm(H) 80cm(W) 60cm(D)

Corner tambour unit




74cm(H) 80cm(W) 60cm(D)

Desk height side filing cabinet



74cm(H) 80cm(W) 60cm(D)

Desk height base unit



74cm(H) 80cm(W) 60cm(D)

finishes - please specify when ordering



beech maple cherry japanese ash light oak modern oak walnut white silver plain grey

Premium range of storage and filing solutions which are stylish & contemporary in design, strong and durable in construction. Available in a range of 10 attractive laminate finishes and guaranteed for 5 years.

Tall tambour units



Code: **Sun-TAMBOUR2**
120cm(H) Tambour unit - 2 shelves
Code: **Sun-TAMBOUR3**
160cm(H) Tambour unit - 3 shelves
Code: **Sun-TAMBOUR4**
180cm(H) Tambour unit - 3 shelves

Desk height tambour unit



Code: **Sun-TAMBOUR1**
Desk height tambour unit - 1 shelf

Corner tambour unit



Code: **Sun-TAMBOUR7**
Corner tambour unit - 1 shelf

Desk height side filing cabinet



Code: **Sun-FCAB1** - Side filing cabinet

Desk height base unit



Code: **Sun-SCAB1** - Desk height base unit

product information - document file organiser/sorter



15 individual compartments
(Document Organiser)



9 or 18 individual compartments
(Document sorter)



Accommodates A4 size paper



FIRA Standard 6250 2005
Specification:
Adhesion Performance
General/Severe



5 year warranty



Supplied fully assembled



British made

product information - filing cabinets



British made



Available in 3 sizes



Single lock, locks all drawers
(2 keys supplies)



5 year warranty



Supplied fully assembled

product information - bookcases



All with adjustable height shelves



FIRA Standard 6250 2005
Specification:
Adhesion Performance
General/Severe



5 year warranty



Supplied fully assembled



British made

dimensions

40cm(w) 50cm(w) 80cm(w)

Document File Organiser
110cm(H) 80cm(W) 33cm(D)
Individual Compartments
19.5cm(H) 24cm(W) 30.5cm(D)

Document Sorters
107cm(H) 30cm(D)
Individual Compartments
8cm(H) 36cm(W) 28cm(D)

Filing Cabinets
48cm(W) 66cm(D)

Bookcases
80cm(W) 41cm(D)

finishes - please specify when ordering



Premium range of storage and filing solutions which are stylish & contemporary in design, strong and durable in construction. Available in a range of 10 attractive laminate finishes and guaranteed for 5 years.

Document file organiser/sorters



Code: **Sun-DOC1**
Document organiser

Code: **Sun-WUSHMDH40-LH**
40cm wide document sorter with sloping top (Left Hand)

Code: **Sun-WUSHMDH50-LH**
50cm wide document sorter with sloping top (Left Hand)

Code: **Sun-WUSHMDH80**
80cm wide document sorter, sloping top with 2 doors

Code: **Sun-WUSHMDH40-RH**
40cm wide Document sorter with sloping top (Right Hand)

Code: **Sun-WUSHMDH50-RH**
50cm wide Document sorter with sloping top (Right Hand)

Filing cabinets



Code: **Sun-FCAB2**
74cm(H) Filing cabinet 2 drawers

Code: **Sun-FCAB3**
105cm(H) Filing cabinet 3 drawers

Code: **Sun-FCAB4**
136cm(H) Filing cabinet 4 drawers

Bookcases



Code: **Sun-BOOK5**
200cm(H) Bookcase 4 shelves

Code: **Sun-BOOK4**
180cm(H) Bookcase 3 shelves

Code: **Sun-BOOK3**
160cm(H) Bookcase 3 shelves

Code: **Sun-BOOK2**
120cm(H) Bookcase 2 shelves

Code: **Sun-BOOK1**
74cm(H) Bookcase 1 shelf

Consultation Chairs

product information

127kg

Maximum user weight (20st.)

Vinyl or inter/vene
Anti-bacterial upholstery

Gas-lift seat height adjustment

HEIGHT RANGE
44-58cm

Adjustable height range

Mid or high back chairs available

Twin levers adjust seat height & back angle

Optional fixed arms

Five star chrome base

Five star black nylon base

accreditation from
FIRA

Tested by FIRA to comply with BS EN 1335 2&3: 2009

FIRA
Independently tested by **FIRA**

Anti-bacterial Vinyl Strength Test
BS EN ISO 13934-1: 2013
BS EN ISO 13973-3: 2000
BS EN ISO 12947-2: 1998

West Yorkshire
Materials Testing Service

Anti-bacterial Vinyl

Fire (Flammability) Testing
Flammability BS 7176 Medium Hazard
BS 5852 : 2006 Clause 11 (composites)

5 Year Warranty
Upholstery Warranty

5 year upholstery wearability warranty (inter/vene only)

eze

Easy clean design and materials

5 Year Warranty

5 year warranty

British Made
Quality Assured

British made

dimensions

Mid Back

44 - 58cm(H)
Seat Width 48cm
Base 60cm Ø

Mid Back (with arms)

44 - 58cm(H)
Seat Width 48cm
Base 60cm Ø

High Back

44 - 58cm(H)
Seat Width 48cm
Base 60cm Ø

High Back (with arms)

44 - 58cm(H)
Seat Width 48cm
Base 60cm Ø

upholstery colours - please specify when ordering

please see inside back cover for full specification of our upholstery range

anti-bacterial vinyl upholstery colours

black grey beige primrose walnut red wine salmon lilac navy mid blue sky blue cool blue green mint white

anti-bacterial inter/vene™ upholstery colours

black grey beige yellow plum red navy sky blue green pastel green

Consultation Chairs

Premium range of fully adjustable consultation chairs, suitable for daily use in healthcare environments. Available in vinyl or inter/vene anti-bacterial upholstery.

Solitaire Chairs (Mid Back)

Mid back, twin lever consultation chairs with five star castor base.



Code: **Sun-CHA20VYL/M**
Consultation chair in anti-bacterial vinyl upholstery

Code: **Sun-CHA20IV/M**
Consultation chair in anti-bacterial inter/vene upholstery

Code: **Sun-CHA20VYL/MS**
Consultation chair with chrome base in anti-bacterial vinyl upholstery

Code: **Sun-CHA20IV/MS**
Consultation chair with chrome base in anti-bacterial inter/vene upholstery

Code: **Sun-CHA20VYL/MA**
Consultation chair with arms in anti-bacterial vinyl upholstery

Code: **Sun-CHA20IV/MA**
Consultation chair with arms in anti-bacterial inter/vene upholstery

Code: **Sun-CHA20VYL/MSA**
Consultation chair with arms and chrome base in anti-bacterial vinyl upholstery

Code: **Sun-CHA20IV/MSA**
Consultation chair with arms and chrome base in anti-bacterial inter/vene upholstery

Solitaire Chairs (High Back)

High back, twin lever consultation chairs with five star castor base.



Code: **Sun-CHA21VYL/H**
Consultation chair in anti-bacterial vinyl upholstery

Code: **Sun-CHA21IV/H**
Consultation chair in anti-bacterial inter/vene upholstery

Code: **Sun-CHA21VYL/HS**
Consultation chair with chrome base in anti-bacterial vinyl upholstery

Code: **Sun-CHA21IV/HS**
Consultation chair with chrome base in anti-bacterial inter/vene upholstery

Code: **Sun-CHA21VYL/HA**
Consultation chair with arms in anti-bacterial vinyl upholstery

Code: **Sun-CHA21IV/HA**
Consultation chair with arms in anti-bacterial inter/vene upholstery

Code: **Sun-CHA21VYL/HSA**
Consultation chair with arms and chrome base in anti-bacterial vinyl upholstery

Code: **Sun-CHA21IV/HSA**
Consultation chair with arms and chrome base in anti-bacterial inter/vene upholstery

Consultation Chairs

product information

127kg

Maximum user weight (20st.)

Vinyl or inter/vene
Anti-bacterial upholstery

Gas-lift seat height adjustment

HEIGHT RANGE
46-58cm

Gas assisted adjustable height

Three levers adjust seat height, seat tilt & back angle

Independent controlled seat slide

Adjustable lumber support (Sun-CHA23/EAALC only)

Fixed arms

Adjustable arms

Five star chrome base

24hr

Designed for 24hr usage

accreditation from
FIRA

Tested by FIRA to comply with BS EN 1335 2&3: 2009

FIRA
Independently tested by **FIRA**

Anti-bacterial Vinyl Strength Test
BS EN ISO 13934-1: 2013
BS EN ISO 13973-3: 2000
BS EN ISO 12947-2: 1998

West Yorkshire Materials Testing Service
Anti-bacterial Vinyl

Fire (Flammability) Testing
Flammability BS 7176 Medium Hazard
BS 5852 : 2006 Clause 11 (composites)

5 Year Upholstery Warranty

5 year upholstery wearability warranty (inter/vene only)

eze

Easy clean design and materials

5 Year Warranty

5 year warranty

British Made Quality Assured

British made

dimensions

No Arms

46 - 58cm(H)
Seat Width 50cm
Base 60cm Ø

Fixed Arms

46 - 58cm(H)
Seat Width 50cm
Base 60cm Ø

Adjustable Arms

46 - 58cm(H)
Seat Width 50cm
Base 60cm Ø

Adjustable Arms & Lumber

46 - 58cm(H)
Seat Width 50cm
Base 60cm Ø

upholstery colours - please specify when ordering
please see inside back cover for full specification of our upholstery range

anti-bacterial vinyl upholstery colours

black	grey	beige	primrose	walnut	red wine	salmon	lavender	navy	mid blue	sky blue	cool blue	green	mint	white

anti-bacterial inter/vene™ upholstery colours

black	grey	beige	yellow	plum	red	navy	sky blue	green	pastel green

Consultation Chairs

Premium range of fully adjustable consultation chairs, suitable for daily use in healthcare environments. Available in vinyl or inter/vene anti-bacterial upholstery.

Quasar Deluxe Chairs

Executive high back, three lever consultation chairs with chrome five star castor base, fixed or adjustable arms with optional lumbar support.



Code: **Sun-CHA23VYL/EC**
Deluxe consultation chair in anti-bacterial vinyl upholstery

Code: **Sun-CHA23IV/EC**
Deluxe consultation chair in anti-bacterial inter/vene upholstery



Code: **Sun-CHA23VYL/EAC**
Deluxe consultation chair with fixed arms in anti-bacterial vinyl upholstery

Code: **Sun-CHA23IV/EAC**
Deluxe consultation chair with fixed arms in anti-bacterial inter/vene upholstery



Code: **Sun-CHA23VYL/EAALC**
Deluxe consultation chair with adjustable arms in anti-bacterial vinyl upholstery

Code: **Sun-CHA23IV/EAALC**
Deluxe consultation chair with adjustable arms in anti-bacterial inter/vene upholstery



Code: **Sun-CHA23VYL/EAALC**
Deluxe consultation chair with adjustable arms and lumbar support in anti-bacterial vinyl upholstery

Code: **Sun-CHA23IV/EAALC**
Deluxe consultation chair with adjustable arms and lumbar support in anti-bacterial inter/vene upholstery

effective.
ergonomic
design

Drug, Medicine, Records Trolleys & Pharmacy Fridges

Designed to provide a complete solution for securely storing, distributing and dispensing medicines and controlled drugs. Our range complies with all the required standards and guidelines and are independently tested for your peace of mind.

		PAGE
	Unit Dosage System (UDS) Trolleys	209-210
	Monitored Dosage System (MDS) Trolleys	211-212
	new EPMA Ward Drug & Medicine Dispensing Trolleys	213-214
	Ward Drug & Medicine Dispensing Trolleys	215-216
	Medical Notes Trolleys	217-218
	Record Transfer Trolleys	219-220
	BIODOSE Medicine Dispensing Trolleys	221-222
	Pharmacy Fridges	223-224

product information - (DT8/UDS Trolley)



9 wide trays and 9 tray dividers



18 narrow trays and 18 tray dividers



Trays provides 54 compartments



6 adjustable height door shelves



Pull-out clear A4 storage tray

product information - (DT7/UDS Trolley)



6 wide trays and 6 tray dividers



12 narrow trays and 12 tray dividers



Trays provides 36 compartments



3 adjustable height door shelves

additional product information



6 pull-out shelves with finger recess



Fold down work surface



Slam lock with 2 keys & lock escutcheon



Integrated cup, pill cup, water jug and pen recess



Integral shelf for patient records file



Push and pull handles, bin & medi-rail



All-round protective rail & front braked castors



Easy clean design and materials



FIRA Standard 6250 2005 Specification: Adhesion Performance General/Severe



FSC The mark of responsible forestry



2 year warranty



British made

accessories * these accessories can only be fitted when ordered with the trolley



Sun-DT/KEY* Additional Door Key



Sun-DT/SL1 Security Leash & Wall Fixing Bracket



Sun-DT/B25* Additional 2.5 Litre Bin and Bracket



Sun-DT/AGH1* Alcohol Gel Holder



Sun-DT/WT Additional Wide Tray



Sun-DT/WTD Additional Wide Tray Divider (each)



Sun-DT/WTL Additional Wide Tray Labels (10 pack)



Sun-DT/NT Additional Narrow Tray



Sun-DT/NTD Additional Narrow Tray Divider (each)



Sun-DT/NTL Additional Narrow Tray Labels (10 pack)

Unit Dosage System (UDS) Trolleys

Stylish & contemporary in design, strong, secure and durable in construction, unit dosage drug trolleys are the ideal solution for storing, distributing and dispensing medicines and patient prescriptions. Hygienically designed and finished in attractive beech effect laminate which is easy to clean.



Code: **Sun-DT8/UDS**
Unit dosage system (UDS) trolley - large



Sun-DT7/UDS
112cm(H)
62cm(W)
60cm(D)




Code: **Sun-DT7/UDS**
Unit dosage system (UDS) trolley - small




Sun-DT8/UDS
112cm(H)
82cm(W)
60cm(D)


product information




Pull-out shelves with finger recess




Location channel for holding racks




Fold down work surface




3 adjustable height door shelves (DT1/MDS4 DT1/MDS6)




6 adjustable height door shelves (DT2/MDS9)




Pull-out clear A4 storage tray (DT2/MDS9)




Slam lock with 2 keys & lock escutcheon



Integrated cup, pill cup & water jug recesses




Integrated pen recess




Push and pull handles, bin & medi-rail




All-round protective rail




Front push and pull handle



100mm Castors with front brakes




Integral shelf for patients record file




Easy clean design and materials




FIRA
Independently tested by FIRA




FSC
The mark of responsible forestry



2 Year Warranty



Supplied fully assembled



British made

FIRA Standard 6250 2005 Specification: Adhesion Performance General/Severe

accessories * these accessories can only be fitted when ordered with the trolley



Sun-DT/KEY*
Additional Door Key



Sun-DT/SL1
Security Leash & Wall Fixing Bracket



Sun-DT/B25*
Additional 2.5 Litre Bin and Bracket



Sun-DT/AGH1*
Alcohol Gel Holder



MDS Self Administration Cabinets
(see pages 122-123)




MDS Medicine Cabinets (see pages 122-123)

- Cards, racks and blister packs are supplied by the pharmacist


dimensions

Sun-DT1/MDS4




112cm(H)
62cm(W)
60cm(D)

Sun-DT1/MDS6



112cm(H)
62cm(W)
60cm(D)

Sun-DT2/MDS9



112cm(H)
82cm(W)
60cm(D)

Stylish & contemporary in design, strong, secure and durable in construction, monitored dosage drug trolleys are the ideal solution for storing, distributing and dispensing medicines and patient prescriptions stored in blister packs, as supplied by Boots & Venalink. Hygienically designed and finished in attractive beech effect laminate which is easy to clean.



Code: Sun-DT2/MDS9
Monitored dosage system (MDS) trolley - large 9 racks



Code: Sun-DT1/MDS6
Monitored dosage system (MDS) trolley - small 6 racks



Code: Sun-DT1/MDS4
Monitored dosage system (MDS) trolley - small 4 racks



product information

Fully adjustable and configurable laptop/tablet mount allowing comfortable working

Laptop or tablet adjustable securing brackets

Solid and stable laptop plate for error free typing

Dual securing/leashing positions

Optional moulded "Leash/Tether Here" wording on dispensing doors (as standard) helping agency staff to secure the trolley correctly

Lock supplied with 2 keys, top and bottom boxes share same keys (keyed to differ)

Suited key solutions including pharmacy "master" keys (keyed to differ)

Fold down work surface with Hi-Vis safety panel, for protecting hands/fingers from being trapped

Hinged lid with slam lock and protective lock escutcheon

Optional integrated lockable lower storage/charging cabinet

Optional integrated recessed storage tray

Pull out storage trays for waste and consumables

Lift up side table

Gas assisted dispensing door (opening and retention)

Adjustable shelf divider system for drug storage

Infection control friendly

Choice of 9 KYDEX® anti-bacterial door colours

KYDEX® antibacterial moulded construction helping infection control

Free running and easy to manoeuvre with multiple push/pull handles

Robust and secure design with a durable and hygienic silver powder coated frame

Three braked castors and one anti-static castor with protective buffers

Integrated mouse cradle (mouse not included)

Supplied fully assembled

FIRA Standard 6250 2005 Specification:
Adhesion Performance General/Severe

FSC The mark of responsible forestry

*** Keyed to differ as specified by NHS England - (keyed to differ means each trolley has its own unique set of keys)

accessories

* these accessories can only be fitted when ordered with the trolley ** includes medi-rail

Sun-EPMA/KD/KEY *
Additional Keys (Suited)*

Sun-EPMA/STL *
Suited Locks

Sun-EPMA/LC3A *
Four 100mm Antistatic Braked Castors

Sun-EPMA/SBOX1 */**
Sharps Box Bracket (Daniels - All Sizes)

Sun-EPMA/SBOX2 */**
Sharps Box Bracket (Frontier) Upto 1 Litre

Sun-EPMA/SBOX3 */**
Sharps Box Bracket (Frontier) 2 & 3 Litre

Sun-EPMA/AGH1 ***
Alcohol Gel Holder

Sun-EPMA/GBH1 */**
Glove Box Holder

Sun-EPMA/DIV1
Additional Divider Kit (5 Shelf Dividers 100 Base Dividers/10 Labels /5 Clips)

Sun-EPMA/BC *
Bar Code Reader Bracket

Sun-EPMA/BW *
Bespoke Wording on Dispensing Cupboard Door

Sun-EPMA/EC *
Electric Cable/Laptop Transformer Cable Management

The NHS is targeted to be paperless by 2020...

...the Sunflower Medical Ward Drug Trolley Range forms part of the EPMA Solution.



The trolleys are designed to provide a complete solution for storing, distributing & dispensing medicine & drugs. Aesthetically pleasing, hygienically designed and combining smart adjustable internal storage with easy access, all trolleys are functional, strong, secure & easy to clean.

Our EPMA Ward Drug Trolleys are available in three sizes, small, medium and large, each size has four options to choose from giving you greater flexibility for your specific needs.

Code: **Sun-EPMA1/LAPTOP/BOX/WHITE** (Shown)
SMALL EPMA Ward Drug Trolley with Laptop Arm and Bracket, Tray and 1 Secure Bottom Storage Cabinet - White Dispensing Cupboard Doors

Code: **Sun-EPMA5/LAPTOP/BOX/WHITE** (Shown)
MEDIUM EPMA Ward Drug Trolley with Laptop Arm and Bracket, Trays and 1 Secure Bottom Storage Cabinet - White Dispensing Cupboard Doors

Code: **Sun-EPMA11/LAPTOP/TRAYS/WHITE** (Shown)
LARGE EPMA Ward Drug Trolley with Laptop Arm and Bracket, Trays and Bottom Shelf - White Dispensing Cupboard Doors



dispensing door colours

(please specify when ordering)
see inside back cover for full specification of our KYDEX® range

Code: **Sun-EPMA11/LAPTOP/TRAYS/WHITE** (Shown)
LARGE EPMA Ward Drug Trolley with Laptop Arm and Bracket, Trays and Bottom Shelf - White Dispensing Cupboard Doors

arctic blue

frost blue

frost

riva

steel green

yellow

sun coral

orange peel

linen

product information



Lightweight easy to manoeuvre design



Adjustable shelf divider system provides upto 5 compartments per shelf



Adjustable cabinet base divider system provides upto 16 compartments



Hinged lid with slam lock & protective lock escutcheon



Lock supplied with 2 keys (Keyed to differ)***



Shelf labels supplied



Gas assisted lid opening & retention



Hygienic stainless steel handle and fittings



Fold down work surface



Lift up side table



4 removable tray positions (2 trays supplied)



Patient notes/file retention strap



Push pull handle (to 3 sides)



Four 100mm braked castors with four protective buffers



Security Leash & Wall Fixing Bracket



Laminate faced MDF with 2mm edging



Easy clean design and materials



FIRA Standard 6250 2005 Specification: Adhesion Performance General/Severe



FSC The mark of responsible forestry



2 year warranty



Supplied fully assembled



British made

*** Keyed to differ as specified by NHS England - (keyed to differ means each trolley has its own unique set of keys)

accessories

* these accessories can only be fitted when ordered with the trolley ** includes medi-rail



Sun-WDT/T2B Additional Double Depth Wide Tray (blue)



Sun-WDT/4Y Additional Double Depth Narrow Tray (yellow)



Sun-WDT/KD/KEY* Additional Keys



Sun-WDT/SBOX1 */** Sharps Box Bracket (Daniels - All Sizes)



Sun-WDT/SBOX2 */** Sharps Box Bracket (Frontier) Upto 1 Litre



Sun-WDT/SBOX3 */** Sharps Box Bracket (Frontier) 2 & 3 Litre



Sun-WDT/DIV1 Additional Divider Kit (5 Shelf Dividers, 10 Base Dividers, 10 Labels and 5 Clips)



Sun-WDT/AGH1 */** Alcohol Gel Holder



Sun-WDT/GBH1 */** Glove Box Holder



Sun-LC3A * Antistatic Braked Castors

Ward Drug & Medicine Dispensing Trolleys

Designed to provide a complete solution for storing, distributing and dispensing medicine and drugs. Ward drug & medicine dispensing trolleys are aesthetically pleasing, hygienically designed and combine smart adjustable internal storage with easy access, while being functional, strong, secure and easy to clean.



Code: **Sun-WDT20/KD**
Ward drug & medicine dispensing trolley (medium capacity) with divider system and 2 storage trays



Sun-WDT20/KD
Wide/Double Depth Trays
16cm(H) 47cm(W) 42cm(D)
Narrow/Double Depth Trays
16cm(H) 31cm(W) 42cm(D)
NB: Maximum tray load 7.5kg
110cm(H) 90cm(W) 58cm(D)
94cm(W) 60cm(D) with buffers



Adjustable storage divider system on shelves and cabinet base



Keyed to Differ as required by NHS England



Code: **Sun-WDT50/KD**
Ward drug & medicine dispensing trolley (large capacity) with divider system and 2 storage trays



Sun-WDT50/KD
Wide/Double Depth Trays
16cm(H) 47cm(W) 42cm(D)
Narrow/Double Depth Trays
16cm(H) 31cm(W) 42cm(D)
NB: Maximum tray load 7.5kg
128cm(H) 90cm(W) 58cm(D)
94cm(W) 60cm(D) with buffers

product information



Lightweight easy to manoeuvre design



Secure locking top 2 keys per lock (keyed to differ)***



Digital locking system



Secure solid sides and base



Accepts foolscap or X-Ray size filing pockets (not included)**



Accepts 470mm & 360mm filing pockets (not included)**



Filing pocket identification label for easy re-ordering



Wipe clean filing pocket tabs (supplied with pockets)



Hinged top(s) open to create additional working space



Hygienic stainless steel handles



2 or 4 removable storage trays



Side storage for patients record file



Push pull handles



Four 100mm braked castors



Four protective buffers



Durable laminate faced MDF with 2mm edging



Easy clean design and materials



FIRA Standard 6250 2005 Specification: Adhesion Performance General/Severe



FSC The mark of responsible forestry



2 year warranty



Supplied fully assembled



British made

*** Keyed to differ as specified by NHS England - (keyed to differ means each trolley has its own unique set of keys)

accessories * these accessories can only be fitted when ordered with the trolley



Sun-LC3A * Antistatic Braked Castors



Sun-MNT/DCL/KEY* Overriding Key for Digital Locking System



Sun-MNT/FPF/PK10 Filing Pocket Foolscap 360mm(W) x 270mm(D) (Pack of 10)



Sun-MNT/WFPF/PK10 Wide Filing Pocket Foolscap 450mm(W) x 270mm(D) (Pack of 10)



Sun-MNT/FPX/PK10 Filing Pocket X-Ray 360mm(W) x 470mm(D) (Pack of 10)



Sun-MNT/WFPX/PK10 Wide Filing Pocket X-Ray 450mm(W) x 470mm(D) (Pack of 10)



Sun-MNT/FPXF/PK10 Twin Filing Pocket X-Ray & Foolscap 360mm(W) x 270/470mm(D) (Pack of 10)



Sun-MNT/WFPXF/PK10 Twin Wide Filing Pocket X-Ray & Foolscap 450mm(W) x 270/470mm(D) (Pack of 10)

** Sun-MNT10 & Sun-MNT10/DCL Trolley can accept all filing pocket sizes but cannot be combined together.

** Sun-MNT20 & Sun-MNT20/DCL Trolley can accept all filing pocket sizes which can be used together but only in opposite sides of the trolley.

Designed to provide a complete solution for storing, distributing and accessing medical notes and X-Rays. Medical notes trolleys are aesthetically pleasing, hygienically designed and combine smart adjustable filing pockets with easy access, while being functional, strong, secure and easy to clean.



Code: **Sun-MNT10**
Medical Notes Trolley, enclosed sides with hinged locking top (Shown)

Code: **Sun-MNT10/DCL**
Medical Notes Trolley, enclosed sides with hinged top and 1 Digital Combination Lock



Optional Digital locking system with optional overriding key (override key(s) must be ordered with trolley)

Sun-MNT10 & 10/DCL



Wide/Single Depth Trays
8cm(H) 47cm(W) 42cm(D)

Wide/Double Depth Trays
16cm(H) 47cm(W) 42cm(D)

NB: Maximum tray load 7.5kg

106cm(H) 64cm(W) 52cm(D)
66cm(W) 57cm(D) with buffers




Keyed to Differ as required by NHS England



Code: **Sun-MNT20**
Medical Notes Trolley, enclosed sides with hinged locking tops (Shown)

Code: **Sun-MNT20/DCL**
Medical Notes Trolley, enclosed sides with hinged top and 2 Digital Combination Locks

Sun-MNT20 & 20/DCL




Wide/Single Depth Trays
8cm(H) 47cm(W) 42cm(D)

Wide/Double Depth Trays
16cm(H) 47cm(W) 42cm(D)


NB: Maximum tray load 7.5kg

106cm(H) 121cm(W) 52cm(D)
121cm(W) 57cm(D) with buffers


product information




Lightweight easy to manoeuvre design




Secure locking top (Sun-MRTTL1 only)




Lock is provided with 2 keys (keyed to differ)***




Locking top folds neatly back




Hygienic stainless steel handle




Maximum load 60kg




Strong deep cabinet for optimum file storage




Easy access to files and documents




2 removable storage trays




Secure solid sides and base




Push pull handles




Durable powder coated blue frame




100mm castors (2 braked)




Four protective buffers




Durable laminate faced MDF with 2mm edging




Easy clean design and materials




FIRA Standard 6250 2005 Specification: Adhesion Performance General/Severe




FSC The mark of responsible forestry



2 year warranty




Supplied fully assembled




British made

*** Keyed to differ as specified by NHS England - (keyed to differ means each trolley has its own unique set of keys)

accessories * these accessories can only be fitted when ordered with the trolley




Sun-LC3A *
Antistatic Braked Castors




Sun-MRTT/KEY*
Additional Key

dimensions



Sun-MRTT1
98cm(H) 81cm(W) 52cm(D)
83cm(W) 57cm(D) with buffers



Sun-MRTTL1
Narrow/Single Depth Trays
8cm(H) 31cm(W) 42cm(D)
NB: Maximum tray load 7.5kg
Storage Box (internal)
25cm(H) 62cm(W) 40cm(D)
98cm(H) 81cm(W) 52cm(D)
83cm(W) 57cm(D) with buffers

Designed to provide a complete solution for storing, distributing and transferring medical notes and files quickly to different locations. Records transfer trolleys are aesthetically pleasing, hygienically designed and easy to access, while being functional, strong, secure and easy to clean.



Code: Sun-MRTT1
Records Transfer Trolley - open



Code: Sun-MRTTL1
Records Transfer Trolley - with folding locking top

product information

72

Holds 72 BODOSE trays

Easy clean design and materials

FIRA Standard 6250 2005 Specification: Adhesion Performance General/Severe

FSC The mark of responsible forestry

3 year warranty

Supplied fully assembled

British made

accessories

dimensions

112cm(H)
92cm(W)
55cm(D)

Stylish & contemporary in design, strong, secure and durable in construction, BODOSE medicine dispensing trolleys are the ideal solution for storing, distributing and dispensing medicines in BODOSE trays. Aesthetically pleasing design which is functional, durable and easy to clean.



Code: **Sun-BIO/MK2/DT69**
BIODOSE medicine dispensing trolley - holds upto 72 trays



product information



5
sizes

5 sizes/capacities available



Wall mounting, worktop mounted, under worktop or freestanding



Lockable doors (2 keys supplied)



Hi-Low temperature & door open alarm

mains failure alarm

Via battery backup



Off cycle automatic defrost

CFC HCFC Ammonia FREE

ECO friendly refrigerant

DH Green Book & RPSGB/GPhC compliant

DH Green Book & RPSGB/GPhC compliant



Fan assisted air cooling for rapid temperature recovery



Internal light fitted as standard on all models

2°C 8°C

Temperature range 2°C to 8°C



Glass doors for visible stock control



Real time temperature display

Min/Max Recording

Min/max temperature recording



Accurate temperature control



Solid doors provide extra security

13 month

Temperature data storage



Adjustable shelves



White painted corrosion proof aluminium interior



Low energy consumption



Easy clean design and materials

5 year warranty

5 year warranty *




Supplied fully assembled



British made

* 5 year warranty upon registration by the end user (UK only) 1 year warranty for export orders

dimensions



Pharmacy Fridges

Designed specifically to meet the demands of health care professionals and for the storage of temperature sensitive vaccines and pharmaceuticals in accordance with manufacturers licensed temperature ranges and guidance. Low energy use, fan assisted cooling, flexible shelving and a 5 year warranty all come as standard.



Code: **Sun-FRIDGE11**
36 Litre wall/worktop mounted fridge - 2 shelves



Code: **Sun-FRIDGE12**
66 Litre wall/worktop mounted fridge - 3 shelves



Code: **Sun-FRIDGE13**
150 Litre under worktop fridge - 3 shelves



Code: **Sun-FRIDGE14**
150 Litre under worktop, glass door fridge - 3 shelves



Code: **Sun-FRIDGE15**
300 Litre free standing fridge - 6 shelves
Code: **Sun-FRIDGE17** (not shown)
430 Litre free standing fridge - 6 shelves











Code: **Sun-FRIDGE16**
300 Litre free standing, glass door fridge - 6 shelves
Code: **Sun-FRIDGE18** (not shown)
430 Litre free standing, glass door fridge - 6 shelves

co-ordinating
personal
storage

Overbed Tables, Bedside Furniture
& Patient Seating

Stylish and contemporary in design, strong and durable in construction and designed to co-ordinate together. Our range of bedside furniture provides the ideal solution for storage and securing personal items. Available in a choice of finishes, all of which are designed to be easy to clean.

	PAGE
 Over Bed Tables	226-231
 new Axis & Axis Plus Bedside Locker Range	232-239
 Bedside Combination Units (Mobile)	240-241
 Bedside Cabinets	242-243
 Drawer Chests	242-243
 Wardrobes	244-245
 Atlas Patient Chairs & Seating	246-249
 Atlas Bariatric Chairs	250-251

NEW

KEYLESS Locker Solutions...



product information

MFC is strong, durable, and resistant to stains and heat

Compact grade laminate is extremely strong & durable, resistant to stains, moisture, heat & scratching

Compact grade laminate has a smooth finish, rounded corners & infection control friendly

C-Shaped Base for overbed or chair use

Parallel leg base for overbed or chair use

1 Piece flat top

Quick and easy 35cm height adjustment (between 74cm-109cm)

Anti-trapping safety feature prevents trapping of patient between the bed & table

Sealed tubular steel base - infection control friendly

4 High quality swivel castors - for ease of movement

Durable and hygienic white epoxy coating base & frame

Safe working load (Spread evenly)

Product Weight

3 finishes
Beech, Blue and White

FIRA
Independently tested by FIRA
FIRA Standard 6250 2005
Specification:
Adhesion Performance
General/Severe

FSC the mark of responsible forestry

Easy clean design and materials

2 year warranty


Supplied fully assembled

British made

dimensions



finishes



COMING SOON

A new range of colours to match our new range of bedside lockers

Astro Overbed/Overchair Tables
(MFC & Compact Grade Laminate Top)



Stylish & contemporary in design, strong and durable in construction, our overbed tables provide the ideal daily solution for multiple uses. Available in 5 specifically designed table tops in a choice of 3 materials, (MFC, Compact Grade Laminate and High Impact PVC/Kydex), 3 finishes and parallel or C-shaped base.

MFC is strong, durable and resistant to stains and heat.



Code: **Sun-OBT1C**
Overbed/Overchair Table with
C-Shaped Base & MFC Flat Top

Code: **Sun-OBT1P**
Overbed/Overchair Table with
Parallel Base & MFC Flat Top

Compact Grade Laminate is a high-density fibreboard with a high strength melamine coating. It is extremely strong and durable, resistant to stains, moisture, heat and scratching. Its smooth finish and rounded edges make it exceptionally easy to clean.




Code: **Sun-OBT2C**
Overbed/Overchair Table with
C-Shaped Base & Flat Top

Code: **Sun-OBT2P**
Overbed/Overchair Table with
Parallel Base & Flat Top

product information



Compact grade laminate is extremely strong & durable, resistant to stains, moisture, heat & scratching



Compact grade laminate has a smooth finish, rounded corners & infection control friendly



C-Shaped Base for overbed or chair use



Parallel leg base for overbed or chair use



1 Piece flat top, with raised lip (Sun-OBT3)



2 Piece tilting top with retaining lip (Sun-OBT4)




Retaining table top lip allows for activities such as reading



Tilting top raises from 0 - 45° (Sun-OBT4)



Quick and easy 35cm height adjustment (between 74cm-109cm)



Anti-trapping safety feature prevents trapping of patient between the bed & table



Sealed tubular steel base - infection control friendly



4 High quality swivel castors - for ease of movement




Durable and hygienic white epoxy coating base & frame



Safe working load



Product Weight



3 finishes

beech, blue and white



FIRA Standard 6250 2005 Specification: Adhesion Performance General/Severe



FSC the mark of responsible forestry



Easy clean design and materials



2 year warranty




Supplied fully assembled




British made

dimensions



109cm
Height Adjustment 35cm
74cm

Overbed/Overchair Table
109cm(H)
95cm(W)
45cm(D)




Tilting Top Raises From 0 - 45°



Table Top
95cm(W) 45cm(D)

finishes



beech blue white



A new range of colours to match our new range of bedside lockers



Astro Overbed/Overchair Tables
(Compact Grade Laminate Top)



Stylish & contemporary in design, strong and durable in construction, our overbed tables provide the ideal daily solution for multiple uses. Available in 5 specifically designed table tops in a choice of 3 materials, (MFC, Compact Grade Laminate and High Impact PVC/Kydex), 3 finishes and parallel or C-shaped base.



Code: **Sun-OBT3C**
Overbed/Overchair Table with C-Shaped Base & Flat Top with 1 Raised Lip



Code: **Sun-OBT3P**
Overbed/Overchair Table with Parallel Base & Flat Top with 1 Raised Lip

Compact Grade Laminate is a high-density fibreboard with a high strength melamine coating. It is extremely strong and durable, resistant to stains, moisture, heat and scratching. Its smooth finish and rounded edges make it exceptionally easy to clean.



Code: **Sun-OBT4C**
Overbed/Overchair Table with C-Shaped Base & 2 Piece Tilting Top with Retaining Lip



Code: **Sun-OBT4P**
Overbed/Overchair Table with Parallel Base & 2 Piece Tilting Top with Retaining Lip

product information



High impact PVC is strong & durable, resistant to stains & heat



High impact PVC has a smooth finish, rounded corners & infection control friendly



C-Shaped Base for overbed or chair use



Parallel leg base for overbed or chair use



1 Piece recessed top (Sun-OBT5)



2 Piece recessed tilting top (Sun-OBT6)




Recessed table top helps to retain liquids and other items



Tilting top raises from 0 - 45°



Quick and easy 35cm height adjustment (between 74cm-109cm)



Anti-trapping safety feature prevents trapping of patient between the bed & table



Sealed tubular steel base - infection control friendly



4 High quality swivel castors - for ease of movement




Durable and hygienic white epoxy coating base & frame



Safe working load



Product Weight



3 finishes
beech, blue and white



FIRA Standard 6250 2005 Specification: Adhesion Performance General/Severe



FSC the mark of responsible forestry



Easy clean design and materials



2 year warranty




Supplied fully assembled




British made

dimensions



finishes



beech blue white



A new range of colours to match our new range of bedside lockers



Astro Overbed/Overchair Tables
(High Impact PVC Top)

Stylish & contemporary in design, strong and durable in construction, our overbed tables provide the ideal daily solution for multiple uses. Available in 5 specifically designed table tops in a choice of 3 materials, (MFC, Compact Grade Laminate and High Impact PVC/Kydex), 3 finishes and parallel or C-shaped base.



Code: **Sun-OBT5C**
Overbed/Overchair Table with C-Shaped Base & Flat Recessed Top



Code: **Sun-OBT5P**
Overbed/Overchair Table with Parallel Base & Flat Recessed Top

High Impact PVC is a heavy duty material which can be formed into recesses, providing a sealed area which can retain liquids and other items. It is strong and durable, resistant to stains and heat. Its smooth finish and rounded edges make it easy to clean and is infection control friendly.



Code: **Sun-OBT6C**
Overbed/Overchair Table with C-Shaped Base & 2 Piece Recessed Tilting Top



Code: **Sun-OBT6P**
Overbed/Overchair Table with Parallel Base & 2 Piece Recessed Tilting Top

product information - (Sun-BSL1)




product information - (Sun-BSL2)




product information - (Sun-BSL3)



product information - All Models



dimensions



Stylish and contemporary in design, strong and durable in construction our Axis Hospital Bedside Locker range provide the ideal solution for storing and securing medicines and personal belongings on the ward. Manufactured from robust MFC material and designed to BS1765-1:1990. Our Axis bedside lockers come with rounded vinyl wrapped matching tops while the Axis Plus range has white tops with upstands on 3 sides in a rounded vinyl wrap finish. They have been designed to be cost effective, strong, easy to clean and offers good infection control performance.

Axis Bedside Lockers



Code: **Sun-BSL1/COLOUR/KA**
Axis Bedside Locker, keyed alike, large locking top drawer, locking small middle drawer, saloon bottom doors

Axis Plus Bedside Lockers



Code: **Sun-BSL1/COLOUR/KA/WT**
Axis Plus Bedside Locker, keyed alike, white upstand top, large locking top drawer, locking small middle drawer, saloon bottom doors

Lock Options



Code: **Sun-BSL2/COLOUR/KA**
Axis Bedside Locker, keyed alike, large locking fold down shelf, open access shelf, saloon bottom doors



Code: **Sun-BSL2/COLOUR/KA/WT**
Axis Plus Bedside Locker, keyed alike, white upstand top, large locking fold down shelf, open access shelf, saloon bottom doors

Finishes



Code: **Sun-BSL3/COLOUR/KA**
Axis Bedside Locker, keyed alike, large locking top drawer, large open access shelf, right or left hand hinged locking bottom cupboard



Code: **Sun-BSL3/COLOUR/KA/WT**
Axis Plus Bedside Locker, keyed alike, white upstand top, large locking top drawer, large open access shelf, right or left hand hinged locking bottom cupboard



product information - (Sun-BSL4)




product information - (Sun-BSL5)




product information - (Sun-BSL6)



product information - All Models



dimensions



Stylish and contemporary in design, strong and durable in construction our Axis Hospital Bedside Locker range provide the ideal solution for storing and securing medicines and personal belongings on the ward. Manufactured from robust MFC material and designed to BS1765-1:1990. Our Axis bedside lockers come with rounded vinyl wrapped matching tops while the Axis Plus range has white tops with upstands on 3 sides in a rounded vinyl wrap finish. They have been designed to be cost effective, strong, easy to clean and offers good infection control performance.

Axis Bedside Lockers



Code: **Sun-BSL4/COLOUR/KA**
Axis Bedside Locker, keyed alike, large locking top drawer, locking small middle drawer, right or left hand hinged locking bottom cupboard

Axis Plus Bedside Lockers




Code: **Sun-BSL4/COLOUR/KA/WT**
Axis Plus Bedside Locker, keyed alike, white upstand top, large locking top drawer, locking small middle drawer, right or left hand hinged locking bottom cupboard

Lock Options



Finishes



Code: **Sun-BSL5/COLOUR/KA**
Axis Bedside Locker, keyed alike, large locking fold down shelf, large open access shelf, right or left hand hinged locking bottom cupboard



Code: **Sun-BSL5/COLOUR/KA/WT**
Axis Plus Bedside Locker, keyed alike, white upstand top, large locking fold down shelf, large open access shelf, right or left hand hinged locking bottom cupboard



Code: **Sun-BSL6/COLOUR/KA**
Axis Bedside Locker, keyed alike, large locking top drawer, open access shelf, saloon bottom doors



Code: **Sun-BSL6/COLOUR/KA/WT**
Axis Plus Bedside Locker, keyed alike, white upstand top, large locking top drawer, open access shelf, saloon bottom doors



Right or left hand hinged cupboard doors (Sun-BSL4 & 5 only)

product information - (Sun-BSL7)




product information - (Sun-BSL8)




product information - (Sun-BSL9/LSA)



product information - All Models



dimensions



Stylish and contemporary in design, strong and durable in construction our Axis Hospital Bedside Locker range provide the ideal solution for storing and securing medicines and personal belongings on the ward. Manufactured from robust MFC material and designed to BS1765-1:1990. Our Axis bedside lockers come with rounded vinyl wrapped matching tops while the Axis Plus range has white tops with upstands on 3 sides in a rounded vinyl wrap finish. They have been designed to be cost effective, strong, easy to clean and offers good infection control performance.

Axis Bedside Lockers




Code: **Sun-BSL7/COLOUR/KA**
Axis Bedside Locker, keyed alike, 2 locking top drawers, open access shelf, saloon bottom doors

Axis Plus Bedside Lockers



Code: **Sun-BSL7/COLOUR/KA/WT**
Axis Plus Bedside Locker, keyed alike, white upstand top, 2 locking top drawers, open access shelf, saloon bottom doors

Lock Options




Code: **Sun-BSL8/COLOUR/KA**
Axis Bedside Locker, keyed alike, 2 locking top drawers, large open access shelf, right or left hand hinged locking bottom cupboard



Code: **Sun-BSL8/COLOUR/KA/WT**
Axis Plus Bedside Locker, keyed alike, white upstand top, 2 locking top drawers, large open access shelf, right or left hand hinged locking bottom cupboard

Finishes



Code: **Sun-BSL9/LSA/COLOUR/KA**
Axis Bedside Locker, keyed alike, locking top drawer, locking fold down shelf, left side access, open access shelf, saloon bottom doors



Code: **Sun-BSL9/LSA/COLOUR/KA/WT**
Axis Plus Bedside Locker, keyed alike, white upstand top, locking top drawer, locking fold down shelf, left side access, open access shelf, saloon bottom doors



product information - (Sun-BSL9/RSA)




product information - (Sun-BSL10/LSA)




product information - (Sun-BSL10/RSA)



product information - All Models



dimensions



Stylish and contemporary in design, strong and durable in construction our Axis Hospital Bedside Locker range provide the ideal solution for storing and securing medicines and personal belongings on the ward. Manufactured from robust MFC material and designed to BS1765-1:1990. Our Axis bedside lockers come with rounded vinyl wrapped matching tops while the Axis Plus range has white tops with upstands on 3 sides in a rounded vinyl wrap finish. They have been designed to be cost effective, strong, easy to clean and offers good infection control performance.

Axis Bedside Lockers



Code: **Sun-BSL9/RSA/COLOUR/KA**
Axis Bedside Locker, keyed alike, locking top drawer, locking fold down shelf, right side access, open access shelf, saloon bottom doors

Axis Plus Bedside Lockers



Code: **Sun-BSL9/RSA/COLOUR/KA/WT**
Axis Plus Bedside Locker, keyed alike, white upstand top, locking top drawer, locking fold down shelf, right side access, open access shelf, saloon bottom doors

Lock Options



Keyed Alike Locks (KA)



Keyed to Differ Locks (KD)



Digital Combination Locks (DCL)



RFID Locks

Finishes



White (WHT) Oak (OAK)



Walnut (WAL) Beech (BCH)



Code: **Sun-BSL10/LSA/COLOUR/KA**
Axis Bedside Locker, keyed alike, locking top drawer, locking fold down shelf, left side access, large open access shelf, right hand hinged locking bottom cupboard



Code: **Sun-BSL10/LSA/COLOUR/KA/WT**
Axis Plus Bedside Locker, keyed alike, white upstand top, locking top drawer, locking fold down shelf, left side access, large open access shelf, right hand hinged locking bottom cupboard



Code: **Sun-BLS10/RSA/COLOUR/KA**
Axis Bedside Locker, keyed alike, locking top drawer, locking fold down shelf, right side access, large open access shelf, left hand hinge locking bottom cupboard



Code: **Sun-BLS10/RSA/COLOUR/KA/WT**
Axis Plus Bedside Locker, keyed alike, white upstand top, locking top drawer, locking fold down shelf, right side access, large open access shelf, left hand hinge locking bottom cupboard



Right or left hand hinged cupboard doors

product information



Available with or without locks (2 keys per lock)

18/25mm MFC

Manufactured chip board construction

18/25mm LFMDf

Manufactured from laminate faced MDF to BS1765-1 : 1990



Conforms to British Standard BS1765-1 : 1990 (LFMDf models only)



Available in left or right hand



Durable 2mm panel edging all round



Chrome hanging rail in wardrobe



60mm castors (front 2 braked)



Stylish flush fitting handles

4 finishes

Maple, beech, oak and white



FIRA Standard 6250 2005 Specification: Adhesion Performance General/Severe



FSC the mark of responsible forestry



Easy clean design and materials



2 year warranty




Supplied fully assembled



British made

dimensions



finishes



Stylish & contemporary in design, strong and durable in construction, our bedside combination units provide the ideal solution for storing and securing personal items. Choose from standard MFC material or LFMDf for compliance with BS1765-1:1990. Available in four finishes, all of which are designed to be easy to clean.



Left or right hand bedside combination unit with single locking wardrobe & hanging rail, locking top drawer, middle open storage recess & large locking cupboard with built in shelf.

Code: **Sun-CBHBC4/MFC/LOCKS**
Left hand bedside combination unit with locks.
Manufactured from MFC material.

Code: **Sun-CBHBC5/MFC/LOCKS**
Right hand bedside combination unit with locks.
Manufactured from MFC material.

Code: **Sun-CBHBC4/LFMDf/LOCKS**
Left hand bedside combination unit with locks.
Manufactured from laminate faced MDF material.

Code: **Sun-CBHBC5/LFMDf/LOCKS**
Right hand bedside combination unit with locks.
Manufactured from laminate faced MDF material.



Left or right hand bedside combination unit with single wardrobe & hanging rail, top drawer, middle open storage recess, & large cupboard with built in shelf.

Code: **Sun-CBHBC4/MFC**
Left hand bedside combination unit.
Manufactured from MFC material.

Code: **Sun-CBHBC5/MFC**
Right hand bedside combination unit.
Manufactured from MFC material.

Code: **Sun-CBHBC4/LFMDf**
Left hand bedside combination unit.
Manufactured from laminate faced MDF material.

Code: **Sun-CBHBC5/LFMDf**
Right hand bedside combination unit.
Manufactured from laminate faced MDF material.

product information



Locking top drawer supplied with 2 keys (Sun-CBBC1, 2 & 3 only)



18mm MFC



One drawer with large open storage recess (Sun-CBBC1 only)



One drawer and one cupboard (Sun-CBBC2 only)



Three drawers (Sun-CBBC3 only)



Five drawers (Sun-CBCD5 only)



Four drawers (Sun-CBCD4 only)



Three drawers (Sun-CBCD3 only)



Durable 2mm panel edging all round



Stylish flush fitting handles



4 finishes



Easy clean design and materials



2 year warranty




Supplied fully assembled




FSC the mark of responsible forestry




dimensions




Bedside cabinet 1 drawer/open
65cm(H)
44cm(W)
56cm(D)




Bedside cabinet 1 drawer/1 door
65cm(H)
44cm(W)
56cm(D)




Bedside cabinet 3 drawers
65cm(H)
44cm(W)
56cm(D)



5 Drawer chest
118cm(H)
80cm(W)
59cm(D)



4 Drawer chest
95cm(H)
80cm(W)
59cm(D)



3 Drawer chest
72cm(H)
80cm(W)
59cm(D)

finishes



Stylish & contemporary in design, strong and durable in construction, our bedside cabinets and drawer chests provide the ideal solution for storing and securing personal items. Manufactured from MFC material. Available in four finishes, all of which are designed to be easy to clean.



Code: **Sun-CBBC1/MFC/LOCK**
Bedside cabinet with locking top drawer & large open storage recess.



Code: **Sun-CBBC2/MFC/LOCK**
Bedside cabinet with locking top drawer & large cupboard



Code: **Sun-CBBC3/MFC/LOCK**
Bedside cabinet with locking top drawer 2 additional drawers



Code: **Sun-CBDC5/MFC** - 5 drawer chest



Code: **Sun-CBDC4/MFC** - 4 drawer chest



Code: **Sun-CBDC3/MFC** - 3 drawer chest

product information

18mm
MFC

Manufactured chip board construction

4
finishes

Maple, beech, oak and white

Easy clean design and materials

2 year warranty

Supplied fully assembled

FSC the mark of responsible forestry

British made

dimensions

HEAVY
ITEM

Single wardrobe
185cm(H)
43cm(W)
60cm(D)

Gents single wardrobe
185cm(H)
43cm(W)
60cm(D)

Double wardrobe
185cm(H)
80cm(W)
60cm(D)

Gents double wardrobe
185cm(H)
80cm(W)
60cm(D)

finishes

maple

beech

oak

white

Stylish & contemporary in design, strong and durable in construction, our range of wardrobes provide the ideal solution for storing personal items. Available in four finishes, all of which are designed to be easy to clean.



Code: **Sun-CBSW1/MFC**
Single wardrobe with hanging rail and 2 shelves



Code: **Sun-CBSW2/MFC**
Gents single wardrobe with hanging rail, 2 shelves and 2 drawers



Code: **Sun-CBDW1/MFC**
Double wardrobe with hanging rail and 2 shelves



Code: **Sun-CBDW2/MFC**
Gents double wardrobe with hanging rail, 2 shelves and 2 large drawers

product information

Maximum user weight (17st.)

Comfortable deep foam upholstery

Solid beech frame

Solid beech arms

Removable seat pad for cleaning

Optional pressure relieving foam upgrade

Optional drop-arms to aid lateral transfer (mid & high back chairs)

Wings for head support (Sun-CHA52 only)

Optional adjustable height legs (10cm extension)

Optional house keeping wheels (fitted to rear legs only)

Vinyl or inter/vene Anti-bacterial upholstery

5 year upholstery wearability warranty (inter/vene only)

2 year warranty

Supplied fully assembled

Easy clean design and materials

BS EN 16139: 2016 Furniture Strength, Durability and Safety Requirements

Anti-bacterial Vinyl Strength Test BS EN ISO 13934-1: 2013 BS EN ISO 13973-3: 2000 BS EN ISO 12947-2: 1998

Fire (Flammability) Testing BS 5852 : 2006 Clause 11 (composites)

British made

accessories * these accessories can only be fitted when ordered with the chair

Sun-CHA48/PRF Pressure Relieving Foam Upgrade for Seat & Back

Sun-CHA49/PRF Pressure Relieving Foam Upgrade for Seat & Back

Sun-CHA50/PRF Pressure Relieving Foam Upgrade for Seat & Back

Sun-CHA50/AHL* Adjustable Height Legs (10cm extension)

Sun-CHA50/HKW* House Keeping Wheels (fitted to rear legs only)

upholstery colours - please specify when ordering
please see inside back cover for full specification of our upholstery range

anti-bacterial vinyl upholstery colours

black	grey	beige	primrose	walnut	red wine	salmon	lilac	navy	mid blue	sky blue	cool blue	green	mint	white
-------	------	-------	----------	--------	----------	--------	-------	------	----------	----------	-----------	-------	------	-------

anti-bacterial inter/vene™ upholstery colours

black	grey	beige	yellow	plum	red	navy	sky blue	green	pastel green
-------	------	-------	--------	------	-----	------	----------	-------	--------------

Atlas Patient Arm Chair Range
Deluxe solid beech frame upholstered seating range with deep foam padded seat and back for optimum comfort.



Code: **Sun-CHA48VYL** - Low back side chair with anti-bacterial vinyl upholstery
Code: **Sun-CHA48IV** - Low back side chair with anti-bacterial inter/vene upholstery



Code: **Sun-CHA49VYL** - Low back arm chair with anti-bacterial vinyl upholstery
Code: **Sun-CHA49IV** - Low back arm chair with anti-bacterial inter/vene upholstery



Code: **Sun-CHA50VYL** - Mid back arm chair with anti-bacterial vinyl upholstery
Code: **Sun-CHA50IV** - Mid back arm chair with anti-bacterial inter/vene upholstery



Code: **Sun-CHA50DAVYL** - Mid back arm chair with drop arms in anti-bacterial vinyl upholstery
Code: **Sun-CHA50DAIV** - Mid back arm chair with drop arms in anti-bacterial inter/vene upholstery

save when buying quantities of 2+ see price list for details

dimensions

Sun-CHA48	Sun-CHA49	Sun-CHA50	Sun-CHA50DA
<p>94cm(H) 54cm(W) 63cm(D)</p>	<p>94cm(H) 60cm(W) 60cm(D)</p>	<p>104cm(H) 60cm(W) 60cm(D)</p>	<p>104cm(H) 60cm(W) 60cm(D)</p>
Seat Heights 48cm - Internal Seat Widths 49cm			

product information

Maximum user weight (17st.)

Comfortable deep foam upholstery

Solid beech frame

Solid beech arms

Removable seat pad for cleaning

Optional pressure relieving foam upgrade

Optional drop-arms to aid lateral transfer (mid & high back chairs)

Wings for head support (Sun-CHA52 only)

Optional adjustable height legs (10cm extension)

Optional house keeping wheels (fitted to rear legs only)

Vinyl or inter/vene Anti-bacterial upholstery

5 year upholstery wearability warranty (inter/vene only)

2 year warranty

Supplied fully assembled

Easy clean design and materials

BS EN 16139: 2016 Furniture Strength, Durability and Safety Requirements

Anti-bacterial Vinyl Strength Test BS EN ISO 13934-1: 2013 BS EN ISO 13973-3: 2000 BS EN ISO 12947-2: 1998

Fire (Flammability) Testing BS 5852 : 2006 Clause 11 (composites)

British made

accessories * these accessories can only be fitted when ordered with the chair

Sun-CHA51/PRF Pressure Relieving Foam Upgrade for Seat & Back

Sun-CHA52/PRF Pressure Relieving Foam Upgrade for Seat & Back

Sun-CHA50/AHL* Adjustable Height Legs (10cm extension)

Sun-CHA50/HKW* House Keeping Wheels (fitted to rear legs only)

upholstery colours - please specify when ordering
please see inside back cover for full specification of our upholstery range

anti-bacterial vinyl upholstery colours

black	grey	beige	primrose	walnut	red wine	salmon	lilac	navy	mid blue	sky blue	cool blue	green	mint	white
-------	------	-------	----------	--------	----------	--------	-------	------	----------	----------	-----------	-------	------	-------

anti-bacterial inter/vene™ upholstery colours

black	grey	beige	yellow	plum	red	navy	sky blue	green	pastel green
-------	------	-------	--------	------	-----	------	----------	-------	--------------

Atlas Patient Arm Chair Range
Deluxe solid beech frame upholstered seating range with deep foam padded seat and back for optimum comfort.



Code: **Sun-CHA51VYL** - High back arm chair with anti-bacterial vinyl upholstery

Code: **Sun-CHA51IV** - High back arm chair with anti-bacterial inter/vene upholstery

Code: **Sun-CHA51DAVYL** - High back arm chair with drop arms in anti-bacterial vinyl upholstery

Code: **Sun-CHA51DAIV** - High back arm chair with drop arms in anti-bacterial inter/vene upholstery



Code: **Sun-CHA52VYL** - High back arm chair with wings in anti-bacterial vinyl upholstery

Code: **Sun-CHA52IV** - High back arm chair with wings in anti-bacterial inter/vene upholstery

Code: **Sun-CHA52DAVYL** - High back arm chair with wings and drop arms in anti-bacterial vinyl upholstery

Code: **Sun-CHA52DAIV** - High back arm chair with wings and drop arms in anti-bacterial inter/vene upholstery

dimensions

Sun-CHA51	Sun-CHA51DA	Sun-CHA52	Sun-CHA52DA
Seat Heights 48cm - Internal Seat Widths 49cm			

product information

254kg

Maximum user weight (40st.)

Comfortable deep foam upholstery

Solid beech frame

Solid beech arms

Removable seat pad for cleaning

Optional pressure relieving foam upgrade

Optional drop-arms to aid lateral transfer (mid & high back chairs)

Wings for head support (Sun-CHA56 only)

Optional adjustable height legs (10cm extension)

Optional house keeping wheels (fitted to rear legs only)

Vinyl or inter/vene Anti-bacterial upholstery

5 year upholstery wearability warranty (inter/vene only)

2 year warranty

Supplied fully assembled

Easy clean design and materials

BS EN 16139: 2016 Furniture Strength, Durability and Safety Requirements

Anti-bacterial Vinyl Strength Test BS EN ISO 13934-1: 2013 BS EN ISO 13973-3: 2000 BS EN ISO 12947-2: 1998

Fire (Flammability) Testing BS 5852 : 2006 Clause 11 (composites)

British made

accessories * these accessories can only be fitted when ordered with the chair

Sun-CHA55/PRF Pressure Relieving Foam Upgrade for Seat & Back

Sun-CHA56/PRF Pressure Relieving Foam Upgrade for Seat & Back

Sun-CHA50/AHL* Adjustable Height Legs (10cm extension)

Sun-CHA50/HKW* House Keeping Wheels (fitted to rear legs only)

upholstery colours - please specify when ordering
please see inside back cover for full specification of our upholstery range

anti-bacterial vinyl upholstery colours

black	grey	beige	primrose	walnut	red wine	salmon	lilac	navy	mid blue	sky blue	cool blue	green	mint	white
-------	------	-------	----------	--------	----------	--------	-------	------	----------	----------	-----------	-------	------	-------

anti-bacterial inter/vene™ upholstery colours

black	grey	beige	yellow	plum	red	navy	sky blue	green	pastel green
-------	------	-------	--------	------	-----	------	----------	-------	--------------

Atlas Bariatric High Back Arm Chair Range
Deluxe solid beech frame upholstered seating range with deep foam padded seat and back for optimum comfort.



Code: **Sun-CHA55VYL** - High back bariatric arm chair with anti-bacterial vinyl upholstery

Code: **Sun-CHA55DAVYL** - High back bariatric arm chair with drop arms in anti-bacterial vinyl upholstery

Code: **Sun-CHA55IV** - High back bariatric arm chair with anti-bacterial inter/vene upholstery

Code: **Sun-CHA55DAIV** - High back bariatric arm chair with drop arms in anti-bacterial inter/vene upholstery



Code: **Sun-CHA56VYL** - High back bariatric arm chair with wings in anti-bacterial vinyl upholstery

Code: **Sun-CHA56DAVYL** - High back bariatric arm chair with wings and drop arms in anti-bacterial vinyl upholstery

Code: **Sun-CHA56IV** - High back bariatric arm chair with wings in anti-bacterial inter/vene upholstery

Code: **Sun-CHA56DAIV** - High back bariatric arm chair with wings and drop arms in anti-bacterial inter/vene upholstery

dimensions

Sun-CHA55	Sun-CHA55DA	Sun-CHA56	Sun-CHA56DA
110cm(H) 93cm(W) 73cm(D)	110cm(H) 93cm(W) 73cm(D)	126cm(H) 93cm(W) 73cm(D)	126cm(H) 93cm(W) 73cm(D)

Seat Heights 51cm - Internal Seat Widths 81cm

save
when buying
quantities of
2+
see price list for details

Recommended Material Cleaning Guidelines:



We make every effort during our design process to ensure our products are easy to clean and utilise modern materials and processes where possible to assist in the prevention of cross infection. We recommend that all surfaces should be thoroughly cleaned and disinfected regularly, and have compiled the following guidelines to assist you with your cleaning process:

- Vinyl & Intervene Upholstery:**

 - Alcohol based wipes are acceptable as a cleaning method such as Clinell Universal Wipes. Vigorous wiping can be used on the vinyl upholstery but the use of abrasive materials is not recommended as it will lead to scratching and damage to the surface.
 - Hypochlorite based cleaning can be undertaken upto 10000ppm. Particular care must be taken to avoid pooling of cleaning substance as concentrations increase with evaporation and cause damage.
 - Hypochlorite cleaning can, over time, cause discolouration and eventually hardening and cracking. In real world use in hospitals Sunflower has not seen this to be a problem but it is a potential issue regardless of the manufacturer of the furniture.
 - All wetness must be dried after cleaning with care been taken to dry into joints.

KYDEX and Plastic Components:

 - Alcohol based wipes are acceptable as a cleaning method such as Clinell Universal Wipes. Vigorous wiping can be used on the Kydex and Plastic but the use of abrasive materials is not recommended as it will lead to scratching and damage to the surface.
 - Hypochlorite based cleaning can be undertaken up to 10000ppm. Particular care must be taken to avoid pooling of cleaning substance as concentrations increase with evaporation and cause damage.
 - All wetness must be dried after cleaning with care been taken to dry into joints.

Powder Coating (Paint):

 - Powder coated surfaces (paint) should be cleaned using a soft damp non-abrasive cloth & wiped dry. Increased or continual exposure to wet/damp conditions or environments of high humidity may result in rusting/or cause the coating to peel, and therefore is not recommended. Some dulling of the finish may occur over time if cleaning products are consistently used.
 - Hypochlorite based cleaning can be undertaken up to 10000ppm. Particular care must be taken to avoid pooling of cleaning substance as concentrations increase with evaporation and cause damage.
 - Bleach, bleach based sprays, solvents, abrasives, synthetic detergents, wax polishes or aerosols are not recommended.

Melamine, Laminates & Vinyl Wrapping:

 - Melamine & Laminate faced surfaces (used on MFC & MDF) should be cleaned using a soft damp non-abrasive cloth and wiped dry.
 - Vinyl wrapped surfaces (e.g. Doors & Drawers) should be cleaned using a soft damp non-abrasive cloth and wiped dry.
 - Hypochlorite based cleaning can be undertaken up to 10000ppm. Particular care must be taken to avoid pooling of cleaning substance as concentrations increase with evaporation and cause damage.
 - Bleach, bleach based sprays, solvents, abrasives, synthetic detergents, wax polishes or aerosols are not recommended.

Stainless Steel:

 - Stainless steel surfaces should be cleaned using a soft damp non-abrasive cloth and wiped dry. Discolouration may occur if allowed to dry naturally. Stainless steel can still rust and certain chemicals can still damage the material. Iodine is one particular chemical which is known to corrode/stain and should be avoided for cleaning. Increased or continual exposure to wet/damp conditions or environments of high humidity may result in some rusting and therefore is not recommended.
 - Hypochlorite based cleaning can be undertaken up to 10000ppm. Particular care must be taken to avoid pooling of cleaning substance as concentrations increase with evaporation and cause damage.
- Aluminium:**

 - Aluminium surfaces should be cleaned using a soft damp non-abrasive cloth and wiped dry. Discolouration may occur if allowed to dry naturally. Increased or continual exposure to wet/damp conditions or environments of high humidity may result in corrosion and therefore is not recommended.
 - Bleach, bleach based sprays, solvents, abrasives, synthetic detergents, wax polishes or aerosols are not recommended.

Plastics, PVC & Rubber Based Materials:

 - Plastics, PVC (e.g. screen panels) & rubber should be cleaned using a soft damp non-abrasive cloth and wiped dry.
 - Bleach, bleach based sprays, solvents, abrasives, synthetic detergents, wax polishes or aerosols are not recommended.

Electrical Components, Wiring & Actuators:

 - Any electrical components should not be exposed to water, liquids or areas of high humidity. Plastic casings, wires and handsets can be cleaned using a soft damp non-abrasive cloth and wiped dry.
 - Bleach, bleach based sprays, solvents, abrasives, synthetic detergents, wax polishes or aerosols are not recommended.

■ Always switch off the mains supply before cleaning.

Metal or MFC Components:

 - Alcohol based wipes are acceptable as a cleaning method such as Clinell Universal Wipes. Vigorous wiping can be used on the body but the use of abrasive materials is not recommended as it will lead to scratching and damage to the surface.
 - Hypochlorite based cleaning can be undertaken upto 10000ppm. Particular care must be taken to avoid pooling of cleaning substance as concentrations increase with evaporation and cause damage.
 - All wetness must be dried after cleaning with care been taken to dry into joints. Failure to dry full can lead to paint work damage and rusting over time.
 - Particular care must be taken no ensure that no liquids or cleaning materials enter the mechanism of locks. This will cause premature failure.

Castors:

 - Castors due to their nature pick up "fluff and bits". Initially before wet cleaning all "fluff and bits" should be removed by hand. Once these are removed then further cleaning can be undertaken as below.
 - Alcohol based wipes are acceptable as a cleaning method such as Clinell Universal Wipes. Vigorous wiping can be used on the castor but the use of abrasive materials is not recommended as it will lead to scratching and damage to the surface.
 - Hypochlorite based cleaning can be undertaken upto and 10000ppm. Particular care must be taken to avoid pooling of cleaning substance as concentrations increase with evaporation and cause damage.
 - All wetness must be dried after cleaning with care been taken to dry into joints.


Shelves and Trays:

 - The shelves and trays should be cleaned individually. They should be fully wiped on all surfaces and repositioned in the trolley as cleaning is undertaken.
 - Alcohol based wipes are acceptable as a cleaning method such as Clinell Universal Wipes. Vigorous wiping can be used on the shelves and trays but the use of abrasive materials is not recommended as it will lead to scratching and damage to the surface.
 - Hypochlorite based cleaning can be undertaken less upto and including 10000ppm. Particular care must be taken to avoid pooling of cleaning substance as concentrations increase with evaporation and cause damage.

We design medical furniture to help healthcare become more efficient, combined with high service levels, extensive warranties and consideration for the environment.

We believe passionately in innovative design and continuous development, utilise modern materials and technology to produce stylish and hygienic new products and solutions.

Ease of cleaning and hygienic design is carefully incorporated into everything we do.



Infection Control

The requirement for furniture that provides functional, ergonomic hygienic storage which assists in minimising infection control issues has always been a significant factor in our design and manufacturing processes and one of continuous improvement. Smart design such as sloping top wall units help facilitate cleaning and the build-up of dust as well as preventing the storage of items which should ideally be stored elsewhere. Materials which are durable and hard wearing, moisture resistant and anti-bacterial can all be utilised as part of our furniture and storage solutions.



Environment

We are committed to reducing our impact on the environment and utilise sustainable materials and sources wherever possible. Our cabinet and worktop manufacturing processes are optimised to be as efficient as possible, reducing energy costs and minimising waste. Reducing our carbon footprint through an on-going reduction programme is part of our company process and our ISO 14001 Environmental Quality Standard Certification and FSC Accreditation (the Mark of Responsible Forestry) shows our commitment to reducing our environmental impact even further.



Certified Number 5288
ISO 9001



Certified Number 5288
ISO 14001



FSC The mark of
responsible forestry



Membership No:
00334186



FIRA Standard 6250
2005 Specification:
Adhesion Performance
General/Severe



BS EN 16139: 2013
Furniture Strength,
Durability & Safety



Anti-bacterial Vinyl Strength Test
BS EN ISO 13934-1: 2013
BS EN ISO 13973-3: 2000
BS EN ISO 12947-2: 1998



Fire (Flammability) Testing
Flammability BS 7176 Medium Hazard
BS 5852 : 2006 Clause 11
(composites)

upholstery colours

Anti-bacterial Vinyl Upholstery

Composition:

- Top Coating 100% PVC
- Backing Cloth 100% Polyester

Flame Retardancy:

- BS5852 Ignition Source 5
- BS7176 Medium Hazard
- CRIB Test BS 585 : 2006

Effectiveness:

- Independent testing has verified effective reduction in the spread of MRSA E.coli and other similar bacteria
- Anti-microbial JIS Z 2801 : 2000

Abrasion:

- Martindale rubs 100,000+

Cleaning:

Alcohol based wipes are acceptable as a cleaning method such as Clinell Universal Wipes. Vigorous wiping can be used on the vinyl upholstery but the use of abrasive materials is not recommended as it will lead to scratching and damage to the surface.

Hypochlorite based cleaning can be undertaken upto 10000ppm. Particular care must be taken to avoid pooling of cleaning substance as concentrations increase with evaporation and cause damage.

Hypochlorite cleaning can, over time, cause discolouration and eventually hardening and cracking. In real world use in hospitals Sunflower has not seen this to be a problem but it is a potential issue regardless of the manufacturer of the furniture.

All wetness must be dried after cleaning with care been taken to dry into joints

colours - please specify when ordering



Inter/Vene™ Anti-bacterial Upholstery

Composition:

- 100% Recycled Polyester FR
- Non metallic dyestuffs

Flame Retardancy:

- Fire Retardant (Class 1) BS7176 (Crib 5)
- BS EN 1021 - 1:2006 (cigarette)
- BS EN 1021 - 2:2006 (match)
- BS 5852 : 2006 Section 4 Ignition Source 5
- BS 7176 : 2007 Medium Hazard

Abrasion:

- Martindale rubs 100,000
- 5 year wearability guarantee (8 hour usage)

Effectiveness:

- Anti-bacterial/ Anti-microbial, inhibits microbial growth of Staph Aureus, Altinata, salmonella, Aspergillus Niger, Penicillium, Staphylococcus Aureus, Cladosporium Herbarism,MRSA
- Completely Waterproof
- Urine Proof BS3424 part 19
- Topical stain repellent and stain resistant finish
- Breathable: Vapour permeable BS3546
- Bleach Cleanable: Colour fastness to chlorine bleach BS EN ISO E03

Cleaning:

Alcohol based wipes are acceptable as a cleaning method such as Clinell Universal Wipes. Vigorous wiping can be used on the vinyl upholstery but the use of abrasive materials is not recommended as it will lead to scratching and damage to the surface.

Hypochlorite based cleaning can be undertaken upto 10000ppm. Particular care must be taken to avoid pooling of cleaning substance as concentrations increase with evaporation and cause damage.

Hypochlorite cleaning can, over time, cause discolouration and eventually hardening and cracking. In real world use in hospitals Sunflower has not seen this to be a problem but it is a potential issue regardless of the manufacturer of the furniture.

All wetness must be dried after cleaning with care been taken to dry into joints

KYDEX® - Thermoplastic Covering

Antifungal & Antibacterial Properties:

KYDEX® sheet resists bacteriological and fungal development which makes it ideal for use in hospitals, food industry, etc. It does not readily provide a source of nutrients for bacteria and fungi. One source of these nutrients come from an additive known as a plasticizer, which in most cases will support microbial growth. Since it does not contain plasticizers, the main source of nutrients is removed and it performs very well when tested for bacterial and fungal growth.

KYDEX® sheet was subjected to fungus resistance testing in accordance with ASTM G-21 and bacteria resistance testing in accordance with ASTM G-22, Procedure B.

The KYDEX® sheet samples did not allow any fungus or bacteria growth and performs exceptionally when exposed to bacteria and fungus, without the addition of an anti-microbial additive.

Chemical Resistance:

From sulphuric acid to brake fluid to carbon tetrachloride, KYDEX® thermoplastic alloy is more resistant to a wider range of concentrated chemicals than any other thermoplastic. All grades of KYDEX® sheet meet the highest standards of chemical resistance for thermoplastic materials and pass the NHS's requirement for 1000ppm Hypochlorite solution cleaning.

FDA Compliance:

The Food and Drug Administration's (FDA's) regulations allows each KYDEX® 107 ingredient to be used in contact with food, making it ideal for overbed tables.


Cleaning:

Alcohol based wipes are acceptable as a cleaning method such as Clinell Universal Wipes. Vigorous wiping can be used on the Kydex but the use of abrasive materials is not recommended as it will lead to scratching and damage to the surface.

Hypochlorite based cleaning can be undertaken upto 10000ppm. Particular care must be taken to avoid pooling of cleaning substance as concentrations increase with evaporation and cause damage.

All wetness must be dried after cleaning with care been taken to dry into joints.

colours - please specify when ordering



colours - please specify when ordering



Sunflower Medical Limited

17 Roydsdale Way
Euroway Trading Estate
Bradford
BD4 6SE
United Kingdom

T: +44 (0)1274 684004

UK Enquiries

E: sales@sunflowermedical.co.uk

International Enquiries

E: export@sunflowermedical.co.uk

W: www.sunflowermedical.co.uk



Certificate Number 5288
ISO 9001



Certificate Number 5288
ISO 14001

