


product information

 Unique glass effect safety trays	 2 or 3 glass effect removable trays	 Lift-off top tray with integral handles	 Hygienically designed raised upstands	 Trays 'locate' around frame to prevent movement
 Tray protection buffer ring	 Durable non-marking 75mm castors	 Durable epoxy coated white frame	 Optional drawers available in a choice of 4 finishes	 4 Drawer finishes
 Maximum individual tray load	 Easy clean design and materials	 2 year warranty	 Supplied fully assembled	 British made

accessories * These accessories can only be fitted when ordered with the trolley

 Sun-DT4D* Single Drawer Unit (small trolleys)	 Sun-DTDT1 Drawer insert for small trolley drawer	 Sun-DT6D* Single Drawer Unit (medium trolleys)	 Sun-DTDT2 Drawer insert for medium trolley drawer	 Sun-DT8D* Double Drawer Unit (large trolleys)
 Sun-DTDT3 Drawer insert for large trolley drawer	 Sun-BUF1* 4 Protective buffers	 Sun-LC1* Two front braked castors	 Sun-LC1A* Antistatic castors (2 front braked)	

dimensions

small  86cm(H) 46cm(W) 52cm(D)	medium  86cm(H) 66cm(W) 52cm(D)	large  86cm(H) 86cm(W) 52cm(D)
---	--	---

drawer finishes - please specify when ordering

 White (high gloss)	 Titanium (high gloss)	 Beech effect	 Maple effect
---	--	---	--

Dressing Trolleys - Glass Effect Safety Trays

Innovative in design, and in use of materials, our range of dressing trolleys provide hygienic mobile support and storage, whilst being functional, durable and easy to clean.

small - 46cm wide dressing trolleys

 Code: Sun-DT4G2 2 glass effect safety trays	 Code: Sun-DT4G3 3 glass effect safety trays	 Sun-DT4G3 with optional Sun-DT4D single drawer unit
--	--	--

medium - 66cm wide dressing trolleys


 Code: Sun-DT6G2 2 glass effect safety trays	 Code: Sun-DT6G3 3 glass effect safety trays	 Sun-DT6G3 with optional Sun-DT6D single drawer unit
--	--	--

large - 86cm wide dressing trolleys

 Code: Sun-DT8G2 2 glass effect safety trays	 Code: Sun-DT8G3 3 glass effect safety trays	 Sun-DT8G3 with optional Sun-DT8D double drawer unit
--	--	--



product information

 Stainless steel trays (304 grade/DP2 finish)	 2 or 3 trays	 Lift-off top tray with integral handles	 Hygienically designed raised upstands	 Trays 'locate' around frame to prevent movement
 Tray protection buffer ring	 Durable non-marking 75mm castors	 Durable epoxy coated white frame	 Optional drawers available in a choice of 4 finishes	 4 Drawer finishes
 Maximum individual tray load	 Easy clean design and materials	 2 Year Warranty	 Supplied fully assembled	 British made

accessories * These accessories can only be fitted when ordered with the trolley

 Sun-DT4D* Single Drawer Unit (small trolleys)	 Sun-DTDT1 Drawer insert for small trolley drawer	 Sun-DT6D* Single Drawer Unit (medium trolleys)	 Sun-DTDT2 Drawer insert for medium trolley drawer	 Sun-DT8D* Double Drawer Unit (large trolleys)
 Sun-DTDT3 Drawer insert for large trolley drawer	 Sun-BUF1* 4 Protective buffers	 Sun-LC1* Two front braked castors	 Sun-LC1A* Antistatic castors (2 front braked)	

dimensions

small  86cm(H) 46cm(W) 52cm(D)	medium  86cm(H) 66cm(W) 52cm(D)	large  86cm(H) 86cm(W) 52cm(D)
---	--	---

drawer finishes - please specify when ordering

 White (high gloss)	 Titanium (high gloss)	 Beech effect	 Maple effect
---	--	---	--

Dressing Trolleys - Stainless Steel Trays

Innovative in design, and in use of materials, our range of dressing trolleys provide hygienic mobile support and storage, whilst being functional, durable and easy to clean.

small - 46cm wide dressing trolleys

 Code: Sun-DT4S2 2 stainless steel trays	 Code: Sun-DT4S3 3 stainless steel trays	 Sun-DT4S3 with optional Sun-DT4D single drawer unit
--	--	---

medium - 66cm wide dressing trolleys

 Code: Sun-DT6S2 2 stainless steel trays	 Code: Sun-DT6S3 3 stainless steel trays	 Sun-DT6S3 with optional Sun-DT6D single drawer unit
--	--	---

large - 86cm wide dressing trolleys

 Code: Sun-DT8S2 2 stainless steel trays	 Code: Sun-DT8S3 3 stainless steel trays	 Sun-DT8S3 with optional Sun-DT8D double drawer unit
--	--	---



product information



Stainless steel trays (304 grade/DP2 finish)



2 or 3 stainless steel removable trays



Lift-off top tray with integral handles



Hygienically designed raised upstands



Trays 'locate' around frame to prevent movement



Tray protection buffer ring



Durable non-marking 75mm castors



Polished stainless steel frame (304 grade)



Maximum individual tray load



Easy clean design and materials



2 year warranty



Supplied fully assembled



British made

accessories

* These accessories can only be fitted when ordered with the trolley



Sun-BUF1* 4 Protective buffers



Sun-LC1* Two front braked castors



Sun-LC1A* Antistatic castors (2 front braked)

dimensions

small



86cm(H) 46cm(W) 52cm(D)

medium



86cm(H) 66cm(W) 52cm(D)

large



86cm(H) 86cm(W) 52cm(D)

Surgical Trolleys - Stainless Steel Trays

Innovative in design, and in use of materials, our range of surgical trolleys provide hygienic mobile support for clinical environments, whilst being functional, durable and easy to clean.

small - 46cm wide surgical trolleys



Code: Sun-ST4S2
2 stainless steel trays



Code: Sun-ST4S3
3 stainless steel trays

medium - 66cm wide surgical trolleys



Code: Sun-ST6S2
2 stainless steel trays



Code: Sun-ST6S3
3 stainless steel trays

large - 86cm wide surgical trolleys



Code: Sun-ST8S2
2 stainless steel trays



Code: Sun-ST8S3
3 stainless steel trays

product information



Removable stainless steel shelf (304 grade/DP2 finish)



Lift-off stainless steel trays (304 grade/DP2 finish)



2 or 3 stainless steel shelf/tray combinations



Shelves & trays 'locate' around frame to prevent movement



Shelf/tray protection buffer ring



Durable non-marking 75mm castors



Polished stainless steel frame (304 grade)



Maximum individual tray load



Easy clean design and materials



Supplied fully assembled



Supplied fully assembled



British made

accessories

* These accessories can only be fitted when ordered with the trolley



Sun-BUF1*
4 Protective buffers




Sun-LC1*
Two front braked castors



Sun-LC1A*
Antistatic castors (2 front braked)

dimensions

small




86cm(H) 46cm(W) 52cm(D)

medium



86cm(H) 66cm(W) 52cm(D)

large



86cm(H) 86cm(W) 52cm(D)

Surgical Trolleys - Stainless Steel Flat Top Shelf & Trays

Innovative in design, and in use of materials, our range of surgical trolleys provide hygienic mobile support for clinical environments, whilst being functional, durable and easy to clean.

small - 46cm wide surgical trolleys



Code: Sun-ST4S2SF
1 stainless steel shelf and 1 tray



Code: Sun-ST4S3SF
1 stainless steel shelf and 2 trays

medium - 66cm wide surgical trolleys



Code: Sun-ST6S2SF
1 stainless steel shelf and 1 tray



Code: Sun-ST6S3SF
1 stainless steel shelf and 2 trays

large - 86cm wide surgical trolleys



Code: Sun-ST8S2SF
1 stainless steel shelf and 1 tray



Code: Sun-ST8S3SF
1 stainless steel shelf and 2 trays

product information

 premium quality Stainless steel shelves, trays and frame (304 grade)	 Mirror polished stainless steel shelves, trays and frame	 Fully welded seamless frame	 BS 4068 : 1977 Compliant *	 3 shelf options Flat, fixed or removable/reversible shelves
 Mirror polished lift off flat shelves (Sun-STFW/RFS2 only)	 Removable, reversible folded shelves (Sun-STFW/RRFS2 only)	 Fully welded fixed shelves (Sun-STFW/FFS2 only)	 Anti-static shelf/tray protection buffers	 75mm antistatic castors
 4 Protective buffers	 Maximum individual tray load 10Kg	 Easy clean design and materials	 2 year warranty	 Supplied fully assembled



* Designed to comply with the requirements of the hospital market (BS4068: 1977)

accessories * These accessories can only be fitted when ordered with the trolley



Sun-LC1A *
Antistatic Castors
(2 front braked)

dimensions

small	medium	large
		
86cm(H) 45cm(W) 45cm(D) 53cm(W) 53cm(D) with buffers	86cm(H) 61cm(W) 45cm(D) 69cm(W) 53cm(D) with buffers	86cm(H) 75cm(W) 45cm(D) 83cm(W) 53cm(D) with buffers

Premium use of materials, and finished to the highest standard, our range of 'mirror polished' surgical trolleys provide hygienic mobile support for clinical environments and theatres, whilst being functional, durable and easy to clean.

small - 45cm wide surgical trolleys (mirror polished)



Code: Sun-STFW4/RFS2
2 removable flat shelves



Code: Sun-STFW4/RRFS2
2 removable reversible folded shelves (flange up or down)



Code: Sun-STFW4/FFS2
2 fully welded fixed shelves (flange down)

medium - 61cm wide surgical trolleys (mirror polished)



Code: Sun-STFW6/RFS2
2 removable flat shelves



Code: Sun-STFW6/RRFS2
2 removable reversible folded shelves (flange up or down)



Code: Sun-STFW6/FFS2
2 fully welded fixed shelves (flange down)

large - 75cm wide surgical trolleys (mirror polished)



Code: Sun-STFW7/RFS2
2 removable flat shelves



Code: Sun-STFW7/RRFS2
2 removable reversible folded shelves (flange up or down)



Code: Sun-STFW7/FFS2
2 fully welded fixed shelves (flange down)

product information

 premium quality Stainless steel shelves, trays and frame (304 grade)	 Mirror polished stainless steel shelves, trays and frame	 Fully welded seamless frame	 BS 4068 : 1977 Compliant *	 3 shelf options Flat, fixed or removable/reversible shelves
 Mirror polished lift off flat shelves (Sun-STFW/RFS2 only)	 Removable, reversible folded shelves (Sun-STFW/RRFS2 only)	 Fully welded fixed shelves (Sun-STFW/FFS2 only)	 Anti-static shelf/tray protection buffers	 75mm antistatic castors
 4 Protective buffers	 Maximum individual tray load	 Easy clean design and materials	 2 year warranty	 Supplied fully assembled

* Designed to comply with the requirements of the hospital market (BS4068: 1977)




accessories

* These accessories can only be fitted when ordered with the trolley



Sun-LC1A *
Antistatic Castors
(2 front braked)

dimension

 86cm(H) 91cm(W) 45cm(D) 99cm(W) 53cm(D) with buffers	 89cm(H) 91cm(W) 61cm(D) 99cm(W) 69cm(D) with buffers	 89cm(H) 122cm(W) 61cm(D) 130cm(W) 69cm(D) with buffers
--	--	---

Premium use of materials, and finished to the highest standard, our range of 'mirror polished' surgical trolleys provide hygienic mobile support for clinical environments and theatres, whilst being functional, durable and easy to clean.

extra large - 91cm wide surgical trolleys (mirror polished)



Code: **Sun-STFW9/RFS2**
2 removable flat shelves

Code: **Sun-STFW9/RRFS2**
2 removable reversible folded shelves (flange up or down)

Code: **Sun-STFW9/FFS2**
2 fully welded fixed shelves (flange down)

extra large/extra deep - 91cm wide surgical trolley (mirror polished)



Code: **Sun-STFW69/FFS2**
2 fully welded fixed shelves (flange down)

extra extra large/extra deep - 122cm wide surgical trolley (mirror polished)








Code: **Sun-STFW612/FFS2**
2 fully welded fixed shelves (flange down)

product information

 Premium quality (304 grade) stainless steel	 Fully welded seamless frame	 Polished (304 grade) stainless steel bowls, trays and buckets	 HEIGHT RANGE 72-108cm Adjustable height bowl stands (Sun-TBSV01/02 only)	 HEIGHT RANGE 74-112cm Adjustable height mayo tables (Sun-TMAY01/02 only)
 Five castor base for stability and manoeuvrability	 Removable bowls	 Removable spiral bowl and bucket buffers	 Removable bucket with handle (18 litre capacity)	 Removable stainless steel instrument tray
 Height adjustment locking handle	 50mm antistatic castors	 eze Easy clean design and materials	 2 Year Warranty 2 Year Warranty	 Supplied Assembled Supplied Assembled

dimensions

 77cm - 108cm(H) 74cm(W) 58cm(D) Bowls 31cm ∅ - 9cm(H)	 77cm - 108cm(H) 58cm(W) 58cm(D) Bowl 31cm ∅ - 9cm(H)	 74cm - 112cm(H) 61cm(W) 45cm(D) Tray 2cm(H) 47cm(W) 36cm(D)	 74cm - 112cm(H) 67cm(W) 47cm(D) Tray 2cm(H) 47cm(W) 36cm(D)
 85cm(H) 75cm(W) 39cm(D) Bowls 31cm ∅ - 9cm(H)	 85cm(H) 45cm(W) 45cm(D) Bowls 31cm ∅ - 9cm(H)	 52cm(H) 38cm(W) 38cm(D) Bucket 29cm ∅ - 28cm(H) 18 Litre capacity	 19cm(H) 31cm(W) 31cm(D) Bowl 31cm ∅ - 9cm(H)

Premium use of materials, and finished to the highest standard, our range of polished stainless steel theatre equipment provide hygienic mobile support and solutions for clinical environments and theatres, whilst being functional, durable and easy to clean.



Code: **Sun-TBSV01**
Single variable height bowl stand with bowl



Code: **Sun-TBSV02**
Double variable height bowl stand with 2 bowls



Code: **Sun-TMAY01**
Mayo table, adjustable height with 2 castors and skid feet



Code: **Sun-TMAY02**
Mayo table, adjustable height with 4 castors



Code: **Sun-TBKS01**
Bucket stand with bucket 18 litre capacity



Code: **Sun-TBSF01**
Two tier double bowl stand with 2 bowls



Code: **Sun-TBSF02**
Side by side double bowl stand with 2 bowls



Code: **Sun-TKAB01**
Kickabout stand with bowl