

product information

Fully powered height adjustment 45-96cm	Fully powered multi-position head section +70° / -35°	Fully powered multi-position foot section +70° / -35°	Maximum user weight (35st.)	Easy to use handset
Retractable positioning castors	Adjustable height levelling foot	Fully upholstered head & foot sections	Double layer density foam for support & comfort	Luxury 65mm deep contoured upholstery
Anti-bacterial/anti-microbial vinyl	Paper roll holder (supplied as standard)	Easy clean design and materials	3 year warranty on all moving components	Lifetime warranty on frame

accessories \* these accessories can only be fitted when ordered with the couch

Sun-FVHC/LEG1 Leg Supports	Sun-FVHC/PHB1 Phlebotomy Arm	Sun-FVHC/FSW1* Foot Switch	Sun-FVHC/FSW2* Foot Switch 2 Section	Sun-FVHC/FSW3* Foot Switch 3 Section
Sun-FVHC/SUP1 Gynae Stirrups	Sun-CHP1 Head Pillow	Sun-FVHC/CBH1* Breathing Hole	Sun-FVHC/CSR Folding Side Rails (pair)	Sun-FVHC/IV2 IV Pole Attachment
Sun-FVHC/RAIL Medi-Rail Attachment	Sun-INST1* Installation Service & Packaging Recycling	Sun-DISPOSAL1* Removal & Disposal of Old Couch		

dimensions

185cm(L) 65cm(W)

**HEAVY ITEM**

anti-bacterial vinyl upholstery colours - please specify when ordering

please see inside back cover for full specification of our vinyl range

black	grey	beige	primrose	walnut	red wine	salmon	lavender	navy	mid blue	sky blue	cool blue	green	mint	white

Fusion Classic Plus+  
2 & 3 Section Fully Powered Electric Couches

The new Fusion Classic Plus+ extends our already popular Fusion range, providing enhanced user functionality and easier patient positioning through its fully powered height, head and foot sections. Fusion Classic Plus+ is simply, strong, stable, aesthetically pleasing and easy to clean.



**Fusion Classic Plus+**

Fully powered height adjustment

Fully powered head section

Fully powered foot section

Code: **Sun-FVHC+3E**  
Fusion Classic Plus+ 3 Section Electric Couch with powered height, head & foot section (Shown with optional accessories)

Normally delivered within 7-10 working days (UK Only)



Code: **Sun-FVHC+2E**  
Fusion Classic Plus+ 2 Section Electric Couch with powered height & head section